

daM3System

Generated by Doxygen 1.8.7

Fri May 29 2015 15:51:52

Contents

1	daM3System code documentation	1
1.1	Introduction	1
2	Hierarchical Index	3
2.1	Class Hierarchy	3
3	Class Index	5
3.1	Class List	5
4	File Index	7
4.1	File List	7
5	Class Documentation	9
5.1	AnalogIn_diligent Class Reference	9
5.1.1	Detailed Description	10
5.1.2	Constructor & Destructor Documentation	10
5.1.2.1	AnalogIn_diligent	10
5.1.3	Member Function Documentation	10
5.1.3.1	_delete	10
5.1.3.2	acquire_adc	10
5.1.3.3	cleanUp	10
5.1.3.4	getHWinfoADC	10
5.1.3.5	getInputRanges	10
5.1.3.6	setPickupVrange	11
5.1.4	Member Data Documentation	11
5.1.4.1	adc_initialized	11
5.1.4.2	adc_meas_duration	11
5.1.4.3	adc_samples_per_second	11
5.1.4.4	ADCCard_device_id	11
5.1.4.5	ADCinputType	11
5.1.4.6	ai	11
5.1.4.7	filter	11
5.1.4.8	hallprobe_volt	11

5.1.4.9	hallprobe_vrange	11
5.1.4.10	ini_hallprobe_vrange	11
5.1.4.11	ini_pickup_vrange	11
5.1.4.12	pickup	11
5.1.4.13	pickup_channel	11
5.1.4.14	pickup_q	11
5.1.4.15	pickup_q_channel	11
5.1.4.16	pickup_vrange	11
5.2	AnalogIn_mcc Class Reference	11
5.2.1	Detailed Description	14
5.2.2	Constructor & Destructor Documentation	14
5.2.2.1	AnalogIn_mcc	14
5.2.3	Member Function Documentation	14
5.2.3.1	_delete	14
5.2.3.2	acquire_adc	14
5.2.3.3	cleanUp	14
5.2.3.4	fetchBuffer_adc	14
5.2.3.5	getHWinfoADC	15
5.2.3.6	getInputRanges	15
5.2.3.7	setHallprobeVrange	15
5.2.3.8	setPickupVrange	15
5.2.3.9	startBackground_adc	15
5.2.3.10	stopBackground_adc	15
5.2.4	Member Data Documentation	15
5.2.4.1	adc_initialized	15
5.2.4.2	adc_meas_duration	15
5.2.4.3	adc_meas_duration_snapshot	15
5.2.4.4	adc_samples_per_second	15
5.2.4.5	ADCCard_ID	16
5.2.4.6	ADCinputType	16
5.2.4.7	ai	16
5.2.4.8	background_acq	16
5.2.4.9	customval_adc	16
5.2.4.10	customval_channel	16
5.2.4.11	feedback_channel	16
5.2.4.12	filter	16
5.2.4.13	hallprobe_channel	16
5.2.4.14	hallprobe_volt	16
5.2.4.15	hallprobe_vrange	16
5.2.4.16	initial_magnet_vout	16

5.2.4.17	pickup	16
5.2.4.18	pickup_channel	17
5.2.4.19	pickup_q	17
5.2.4.20	pickup_q_channel	17
5.2.4.21	pickup_vrange	17
5.2.4.22	stddev_hallprobe_volt	17
5.2.4.23	stddev_pickup	17
5.2.4.24	stddev_pickup_q	17
5.2.4.25	VOut_feedback_vrange	17
5.3	AnalogIO Class Reference	17
5.3.1	Detailed Description	21
5.3.2	Constructor & Destructor Documentation	21
5.3.2.1	AnalogIO	21
5.3.3	Member Function Documentation	21
5.3.3.1	_delete	21
5.3.3.2	_delete	21
5.3.3.3	acquire	21
5.3.3.4	acquire_adc	22
5.3.3.5	acquire_snapshot	22
5.3.3.6	checkFieldOffset	22
5.3.3.7	constrain	22
5.3.3.8	failSafe	22
5.3.3.9	fetchBuffer	22
5.3.3.10	fetchBuffer_adc	23
5.3.3.11	getDACout	23
5.3.3.12	getHWinfo	23
5.3.3.13	getHWinfoADC	23
5.3.3.14	getHWinfoDAC	23
5.3.3.15	getInputRanges	23
5.3.3.16	getOutputRanges	23
5.3.3.17	getReachTime	23
5.3.3.18	getReachTimeField	24
5.3.3.19	interrupt	24
5.3.3.20	reachMagnet	24
5.3.3.21	resetFieldOffset	24
5.3.3.22	setHallprobeVrange	24
5.3.3.23	setMField	24
5.3.3.24	setPickupVrange	25
5.3.3.25	startBackground_adc	25
5.3.3.26	stopBackground_adc	25

5.3.4	Member Data Documentation	25
5.3.4.1	adc_initialized	25
5.3.4.2	adc_meas_duration	25
5.3.4.3	adc_meas_duration_snapshot	25
5.3.4.4	adc_samples_per_second	25
5.3.4.5	ADCCard_ID	25
5.3.4.6	ADCinputType	25
5.3.4.7	ai	26
5.3.4.8	ao	26
5.3.4.9	background_acq	26
5.3.4.10	customval_adc	26
5.3.4.11	customval_channel	26
5.3.4.12	dac_initialized	26
5.3.4.13	dac_samples_per_second	26
5.3.4.14	feedback_channel	26
5.3.4.15	filter	26
5.3.4.16	hallprobe_channel	26
5.3.4.17	hallprobe_volt	26
5.3.4.18	hallprobe_vrange	26
5.3.4.19	initial_magnet_vout	26
5.3.4.20	initialized	26
5.3.4.21	magnet_powersupply_vrange	26
5.3.4.22	magnet_vout	27
5.3.4.23	magnetic_field	27
5.3.4.24	magnetic_field_setpoint	27
5.3.4.25	MField	27
5.3.4.26	pickup	27
5.3.4.27	pickup_channel	27
5.3.4.28	pickup_q	27
5.3.4.29	pickup_q_channel	27
5.3.4.30	pickup_vrange	27
5.3.4.31	putdata_minsamples	27
5.3.4.32	stddev_hallprobe_volt	27
5.3.4.33	stddev_pickup	27
5.3.4.34	stddev_pickup_q	27
5.3.4.35	VOut_feedback_vrange	28
5.3.4.36	yoke_offsets	28
5.4	AnalogOut_mcc Class Reference	28
5.4.1	Detailed Description	30
5.4.2	Constructor & Destructor Documentation	30

5.4.2.1	AnalogOut_mcc	30
5.4.3	Member Function Documentation	30
5.4.3.1	_delete	30
5.4.3.2	cleanUp	30
5.4.3.3	constrain	30
5.4.3.4	getDACout	30
5.4.3.5	getHWinfoDAC	30
5.4.3.6	getOutputRanges	31
5.4.3.7	getReachTime	31
5.4.3.8	reachMagnet	32
5.4.4	Member Data Documentation	32
5.4.4.1	ao	32
5.4.4.2	dac_initialized	32
5.4.4.3	dac_ramp	32
5.4.4.4	dac_samplerate	32
5.4.4.5	dac_samples_per_second	32
5.4.4.6	magnet_powersupply_channel	32
5.4.4.7	magnet_powersupply_vrange	32
5.4.4.8	magnet_ramp	32
5.4.4.9	magnet_volt_setpoint	32
5.4.4.10	magnet_vout	33
5.4.4.11	magnetic_field_setpoint	33
5.4.4.12	putdata_minsamples	33
5.5	ArduinoIO Class Reference	33
5.5.1	Constructor & Destructor Documentation	34
5.5.1.1	ArduinoIO	34
5.5.2	Member Function Documentation	34
5.5.2.1	_delete	34
5.5.2.2	processCommands	34
5.5.2.3	processReadbuffer	34
5.5.2.4	sendData	34
5.5.2.5	startConnection	34
5.5.2.6	stopConnection	34
5.5.2.7	timedReceive	34
5.5.2.8	waitforAnswerTo	34
5.5.3	Member Data Documentation	34
5.5.3.1	await_respond	34
5.5.3.2	baud_rate	34
5.5.3.3	Connection	34
5.5.3.4	port	35

5.5.3.5	read_buffer	35
5.5.3.6	ReceiveTimer	35
5.5.3.7	updateInterval	35
5.6	CfgBoolean Class Reference	35
5.6.1	Detailed Description	36
5.6.2	Constructor & Destructor Documentation	36
5.6.2.1	CfgBoolean	36
5.6.3	Member Function Documentation	36
5.6.3.1	chkbox_Callback	36
5.6.3.2	chkboxChanged	36
5.6.3.3	getValue	36
5.6.3.4	setValue	36
5.6.4	Member Data Documentation	36
5.6.4.1	cfg_property	36
5.6.4.2	CfgObj	36
5.6.4.3	chkbox	36
5.6.4.4	ori_chkbox_callback	36
5.6.4.5	value	37
5.7	CfgRange Class Reference	37
5.7.1	Detailed Description	38
5.7.2	Constructor & Destructor Documentation	38
5.7.2.1	CfgRange	38
5.7.3	Member Function Documentation	38
5.7.3.1	_delete	39
5.7.3.2	checkValue	39
5.7.3.3	getValue	40
5.7.3.4	internal_Callback	40
5.7.3.5	setValue	40
5.7.4	Member Data Documentation	40
5.7.4.1	cfg_obj	40
5.7.4.2	cfg_varname	40
5.7.4.3	default_val	40
5.7.4.4	hObject	40
5.7.4.5	ori_Callback	40
5.7.4.6	value	40
5.8	CfgStr Class Reference	41
5.8.1	Detailed Description	42
5.8.2	Constructor & Destructor Documentation	42
5.8.2.1	CfgStr	42
5.8.3	Member Function Documentation	42

5.8.3.1	_delete	42
5.8.3.2	checkValue	42
5.8.3.3	getValue	43
5.8.3.4	internal_Callback	43
5.8.3.5	setValue	43
5.8.4	Member Data Documentation	43
5.8.4.1	cfg_obj	43
5.8.4.2	cfg_varname	43
5.8.4.3	default_val	43
5.8.4.4	hObject	43
5.8.4.5	ori_Callback	43
5.8.4.6	value	43
5.9	CfgStrOrNum Class Reference	43
5.9.1	Detailed Description	45
5.9.2	Constructor & Destructor Documentation	45
5.9.2.1	CfgStrOrNum	45
5.9.3	Member Function Documentation	45
5.9.3.1	_delete	45
5.9.3.2	checkValue	45
5.9.3.3	getValue	46
5.9.3.4	internal_Callback	46
5.9.3.5	setValue	46
5.9.4	Member Data Documentation	46
5.9.4.1	cfg_obj	46
5.9.4.2	cfg_varname	46
5.9.4.3	default_val	46
5.9.4.4	hObject	46
5.9.4.5	ori_Callback	46
5.9.4.6	value	46
5.10	Config Class Reference	46
5.10.1	Detailed Description	51
5.10.2	Constructor & Destructor Documentation	51
5.10.2.1	Config	51
5.10.3	Member Function Documentation	51
5.10.3.1	getSingletonData	51
5.10.3.2	getTest	51
5.10.3.3	hasChanged	51
5.10.3.4	instance	51
5.10.3.5	instance	51
5.10.3.6	loadFrom	52

5.10.3.7	saveDialog	52
5.10.3.8	saveMe	52
5.10.3.9	saveMeAs	52
5.10.3.10	setSingletonData	52
5.10.4	Member Data Documentation	52
5.10.4.1	adc_customval_channel	52
5.10.4.2	adc_customval_factor	52
5.10.4.3	adc_duration	52
5.10.4.4	adc_feedback_channel	52
5.10.4.5	adc_hallprobe_channel	52
5.10.4.6	adc_hallprobe_vrange	52
5.10.4.7	adc_lockin_sensitivity	52
5.10.4.8	adc_pickup_channel	53
5.10.4.9	adc_pickup_q_channel	53
5.10.4.10	adc_pickup_vrange	53
5.10.4.11	adc_samplerate	53
5.10.4.12	adc_samplerate_cont	53
5.10.4.13	adc_snapshot_duration	53
5.10.4.14	adc_use_customval	53
5.10.4.15	adc_use_feedback_channel	53
5.10.4.16	ADCCard_ID	53
5.10.4.17	ADCDriver	53
5.10.4.18	autoload_on_first_use	53
5.10.4.19	check_lockin_sensitivity	53
5.10.4.20	check_saveconfig_before_meas	53
5.10.4.21	comment	53
5.10.4.22	config_filename	53
5.10.4.23	dac_magnet_powersupply_channel	54
5.10.4.24	dac_ramp	54
5.10.4.25	dac_samplerate	54
5.10.4.26	DACCard_ID	54
5.10.4.27	DACDriver	54
5.10.4.28	daq_filter_idx	54
5.10.4.29	datfile_path	54
5.10.4.30	err_cancel_measure	54
5.10.4.31	err_IO_init	54
5.10.4.32	file_ext	54
5.10.4.33	gui_bar_magnet_field	54
5.10.4.34	gui_bar_magnet_vout	54
5.10.4.35	gui_dataPlot	54

5.10.4.36 gui_figCfgMain	54
5.10.4.37 gui_figFieldControl	54
5.10.4.38 gui_fmPlot	55
5.10.4.39 gui_loop_pause	55
5.10.4.40 gui_SweepPanel	55
5.10.4.41 gui_txtCustomADC	55
5.10.4.42 gui_txtFieldSetpoint	55
5.10.4.43 gui_txtPickup	55
5.10.4.44 gui_txtPickupQ	55
5.10.4.45 gui_yoke_offsets	55
5.10.4.46 guiInterlock	55
5.10.4.47 hallprobe_factor	55
5.10.4.48 hallprobe_offset	55
5.10.4.49 intermediate_datasaving	55
5.10.4.50 large_field_step	56
5.10.4.51 magnet_calibration_file_nr	56
5.10.4.52 magnet_calibration_max_volt	56
5.10.4.53 magnet_calibration_min_volt	56
5.10.4.54 magnet_calibration_step_volt	56
5.10.4.55 magnet_dynamic_offset	56
5.10.4.56 magnet_max_field	56
5.10.4.57 magnet_min_field	56
5.10.4.58 manual_temperature_val	56
5.10.4.59 mplotview_magnetcalibration	56
5.10.4.60 mplotview_measure	56
5.10.4.61 mplotview_std	56
5.10.4.62 mplotview_viewcalibration	57
5.10.4.63 mplotview_viewdatafile	57
5.10.4.64 pause_between_datapoints	57
5.10.4.65 pause_large_fieldstep	57
5.10.4.66 plot_show_average	57
5.10.4.67 plot_show_points	57
5.10.4.68 plot_show_quadrature	57
5.10.4.69 plot_window	57
5.10.4.70 RangeTable	57
5.10.4.71 runstate	57
5.10.4.72 runstate_pleasestop	57
5.10.4.73 runstate_pleasestop_user	57
5.10.4.74 runstate_running	58
5.10.4.75 runstate_stopped	58

5.10.4.76	sw_version	58
5.10.4.77	SweepTable	58
5.10.4.78	test	58
5.10.4.79	timestamp_end	58
5.10.4.80	timestamp_start	58
5.10.4.81	use_customcontroller	58
5.10.4.82	use_tempcontroller	58
5.10.4.83	winpos_FieldCtl	58
5.10.4.84	winpos_MainCfg	58
5.10.4.85	winpos_MPlot	58
5.11	CustomController Class Reference	59
5.11.1	Constructor & Destructor Documentation	60
5.11.1.1	CustomController	60
5.11.2	Member Function Documentation	60
5.11.2.1	_delete	60
5.11.2.2	abortMeasurement	60
5.11.2.3	afterAllLoops	60
5.11.2.4	afterLoop	60
5.11.2.5	beforeAllLoops	60
5.11.2.6	beforeLoop	60
5.11.2.7	cyclicTimer	60
5.11.2.8	getValue	60
5.11.2.9	setupTimer	60
5.11.2.10	setValue	60
5.11.2.11	stopTimer	60
5.11.2.12	timedFunction	60
5.11.3	Member Data Documentation	60
5.11.3.1	controller_setpoint	60
5.11.3.2	controller_value	61
5.11.3.3	CyclicTimer	61
5.11.3.4	filename_suffix	61
5.12	DataFilter Class Reference	61
5.12.1	Detailed Description	62
5.12.2	Member Enumeration Documentation	62
5.12.2.1	DataFilter	62
5.12.3	Member Function Documentation	62
5.12.3.1	filter	62
5.12.4	Member Data Documentation	62
5.12.4.1	trimmean_percentage	62
5.13	DataPlot Class Reference	62

5.13.1 Detailed Description	64
5.13.2 Constructor & Destructor Documentation	65
5.13.2.1 DataPlot	65
5.13.3 Member Function Documentation	65
5.13.3.1 _delete	65
5.13.3.2 configMagnetView	65
5.13.3.3 configPickupView	65
5.13.3.4 updateAverage	65
5.13.3.5 updateData	65
5.13.4 Member Data Documentation	66
5.13.4.1 dataLoopAvg	66
5.13.4.2 dataObject	66
5.13.4.3 grid	66
5.13.4.4 label_Magnet_X	66
5.13.4.5 label_Magnet_Y	66
5.13.4.6 label_Pickup_X	66
5.13.4.7 label_Pickup_Y	66
5.13.4.8 legend_magnet	66
5.13.4.9 legend_pickup	66
5.13.4.10 lineWidth	66
5.13.4.11 marker	66
5.13.4.12 pickup_xlim	66
5.13.4.13 pickup_ylim	66
5.13.4.14 plotAxes	66
5.13.4.15 plotIdx	66
5.13.4.16 plotObj1	67
5.13.4.17 plotObj2	67
5.13.4.18 plotObj3	67
5.13.4.19 plotType	67
5.13.4.20 show_average	67
5.13.4.21 show_points	67
5.13.4.22 show_quadrature	67
5.13.4.23 typeMagnet	67
5.13.4.24 typeNone	67
5.13.4.25 typePickup	67
5.13.4.26 xlabel_nonexp	67
5.13.4.27 ylim_factor	67
5.14 ErrorHandler Class Reference	68
5.14.1 Detailed Description	68
5.14.2 Member Function Documentation	68

5.14.2.1	getList	68
5.14.2.2	getLog	69
5.14.2.3	logError	69
5.14.2.4	report	69
5.14.3	Member Data Documentation	69
5.14.3.1	logfile	69
5.14.3.2	maxlogs	69
5.15	Fake_AnalogInAndOut Class Reference	69
5.15.1	Detailed Description	72
5.15.2	Constructor & Destructor Documentation	72
5.15.2.1	Fake_AnalogInAndOut	72
5.15.3	Member Function Documentation	72
5.15.3.1	_delete	72
5.15.3.2	acquire_adc	72
5.15.3.3	cleanUp	72
5.15.3.4	fakeHallSig	72
5.15.3.5	fetchBuffer_adc	72
5.15.3.6	getDACout	72
5.15.3.7	getHWinfoADC	72
5.15.3.8	getHWinfoDAC	72
5.15.3.9	getInputRanges	72
5.15.3.10	getOutputRanges	72
5.15.3.11	getReachTime	72
5.15.3.12	reachMagnet	72
5.15.3.13	setHallprobeVrange	72
5.15.3.14	setPickupVrange	72
5.15.3.15	startBackground_adc	72
5.15.3.16	stopBackground_adc	72
5.15.4	Member Data Documentation	72
5.15.4.1	adc_initialized	72
5.15.4.2	adc_meas_duration	72
5.15.4.3	adc_meas_duration_snapshot	72
5.15.4.4	adc_samples_per_second	72
5.15.4.5	ADCCard_ID	72
5.15.4.6	ADCinputType	73
5.15.4.7	ai	73
5.15.4.8	ao	73
5.15.4.9	background_acq	73
5.15.4.10	customval_adc	73
5.15.4.11	dac_initialized	73

5.15.4.12 DAC_ramp	73
5.15.4.13 DAC_samplerate	73
5.15.4.14 dac_samples_per_second	73
5.15.4.15 DACCard_ID	73
5.15.4.16 fake_mag_curve_dn	73
5.15.4.17 fake_mag_curve_up	73
5.15.4.18 fake_yoke	73
5.15.4.19 filter	73
5.15.4.20 hallprobe_channel	73
5.15.4.21 hallprobe_volt	73
5.15.4.22 hallprobe_vrange	73
5.15.4.23 ini_hallprobe_vrange	73
5.15.4.24 ini_pickup_vrange	73
5.15.4.25 initial_magnet_vout	73
5.15.4.26 magnet_powersupply_channel	73
5.15.4.27 magnet_powersupply_vrange	73
5.15.4.28 magnet_ramp	73
5.15.4.29 magnet_volt_setpoint	74
5.15.4.30 magnet_vout	74
5.15.4.31 magnetic_field_setpoint	74
5.15.4.32 pickup	74
5.15.4.33 pickup_channel	74
5.15.4.34 pickup_q	74
5.15.4.35 pickup_q_channel	74
5.15.4.36 pickup_vrange	74
5.15.4.37 SamplesOutput	74
5.15.4.38 stddev_hallprobe_volt	74
5.15.4.39 stddev_pickup	74
5.15.4.40 stddev_pickup_q	74
5.15.4.41 use_feedback_channel	74
5.16 GuiCfg Class Reference	75
5.16.1 Member Data Documentation	76
5.16.1.1 cfgmain_show_ask_lockin	76
5.16.1.2 cfgmain_show_ask_saveconfig_before_meas	76
5.16.1.3 cfgmain_signal1ch	76
5.16.1.4 cfgmain_signal2ch	76
5.16.1.5 fctl_customADC	76
5.16.1.6 fctl_customADC_unit	76
5.16.1.7 fctl_field_dn	76
5.16.1.8 fctl_field_dndn	76

5.16.1.9	fctl_field_up	76
5.16.1.10	fctl_field_upup	76
5.16.1.11	fctl_fielddtxt_dn	76
5.16.1.12	fctl_fielddtxt_dndn	76
5.16.1.13	fctl_fielddtxt_up	76
5.16.1.14	fctl_fielddtxt_upup	76
5.16.1.15	fctl_pickup	76
5.16.1.16	fctl_quadr	77
5.16.1.17	file_dataformat_header	77
5.16.1.18	plot_label_Pickup_Y	77
5.16.1.19	plot_legend_magnet	77
5.16.1.20	plot_legend_pickup	77
5.16.1.21	plotw_aux	77
5.16.1.22	plotw_lockin	77
5.16.1.23	plotw_pickup	77
5.16.1.24	plotw_quadr	77
5.16.1.25	plotw_showquadr	77
5.17	hgsetget Class Reference	77
5.18	HWController Class Reference	78
5.18.1	Constructor & Destructor Documentation	79
5.18.1.1	HWController	79
5.18.2	Member Function Documentation	79
5.18.2.1	_delete	79
5.18.2.2	abortMeasurement	79
5.18.2.3	afterAllLoops	79
5.18.2.4	afterLoop	79
5.18.2.5	beforeAllLoops	79
5.18.2.6	beforeLoop	79
5.18.2.7	cyclicTimer	79
5.18.2.8	getValue	79
5.18.2.9	setupTimer	79
5.18.2.10	setValue	79
5.18.2.11	stopTimer	79
5.18.2.12	timedFunction	79
5.18.3	Member Data Documentation	79
5.18.3.1	controller_setpoint	79
5.18.3.2	controller_value	79
5.18.3.3	CyclicTimer	80
5.18.3.4	filename_suffix	80
5.19	InterlockGUIElements Class Reference	80

5.19.1 Detailed Description	81
5.19.2 Member Function Documentation	81
5.19.2.1 _delete	81
5.19.2.2 addControlEl	81
5.19.2.3 addPanel	81
5.19.2.4 disableItems	81
5.19.2.5 enableItems	81
5.19.2.6 garbageCollect	81
5.19.3 Member Data Documentation	81
5.19.3.1 Items	81
5.19.3.2 items_enabled	81
5.19.3.3 states	81
5.20 LevelBar Class Reference	81
5.20.1 Detailed Description	83
5.20.2 Constructor & Destructor Documentation	83
5.20.2.1 LevelBar	83
5.20.3 Member Function Documentation	83
5.20.3.1 _delete	83
5.20.3.2 setMinMaxDisplay	83
5.20.3.3 setTextDisplay	84
5.20.3.4 setValue	84
5.20.4 Member Data Documentation	84
5.20.4.1 laxes	84
5.20.4.2 lbar	84
5.20.4.3 MaxDisplay	84
5.20.4.4 maxval	84
5.20.4.5 MinDisplay	84
5.20.4.6 minval	84
5.20.4.7 TextDisplay	85
5.21 MagField Class Reference	85
5.21.1 Detailed Description	86
5.21.2 Constructor & Destructor Documentation	86
5.21.2.1 MagField	86
5.21.3 Member Function Documentation	87
5.21.3.1 _delete	87
5.21.3.2 calibrationFiles	87
5.21.3.3 FieldToVolt	88
5.21.3.4 getDate	88
5.21.3.5 getDescription	88
5.21.3.6 HallprobeToField	88

5.21.3.7	polyFitMatrix	89
5.21.3.8	save_calib	89
5.21.3.9	setDescription	89
5.21.3.10	VoltToField	89
5.21.4	Member Data Documentation	89
5.21.4.1	calibration_file_nr	89
5.21.4.2	description	89
5.21.4.3	field_matrix	89
5.21.4.4	field_poly	89
5.21.4.5	max_yoke_offset	90
5.21.4.6	not_calibrated	90
5.21.4.7	num_calibration_files	90
5.21.4.8	poly_order	90
5.21.4.9	yoke_offset	90
5.22	MagnetizationCurve Class Reference	90
5.22.1	Detailed Description	92
5.22.2	Constructor & Destructor Documentation	92
5.22.2.1	MagnetizationCurve	92
5.22.3	Member Function Documentation	92
5.22.3.1	addDataPoint	92
5.22.3.2	addElementStruct	92
5.22.3.3	averageLoop	93
5.22.3.4	datapoints	94
5.22.3.5	datapoints	94
5.22.4	Member Data Documentation	94
5.22.4.1	customval	94
5.22.4.2	customval_adc	94
5.22.4.3	dataformat_descriptor	94
5.22.4.4	dataformat_descriptor_read	94
5.22.4.5	dataformat_header	94
5.22.4.6	datapoints	94
5.22.4.7	field	94
5.22.4.8	field_setpoint	94
5.22.4.9	field_stdev	95
5.22.4.10	signal	95
5.22.4.11	signal_quadr	95
5.22.4.12	signal_stdev	95
5.22.4.13	temperature	95
5.23	Measurement Class Reference	95
5.23.1	Detailed Description	95

5.23.2	Member Function Documentation	95
5.23.2.1	degaussYoke	95
5.23.2.2	ipause	96
5.23.2.3	updateMagnetVout	96
5.23.2.4	viewFile	96
5.24	MeasurementFile Class Reference	96
5.24.1	Detailed Description	97
5.24.2	Constructor & Destructor Documentation	97
5.24.2.1	MeasurementFile	97
5.24.3	Member Function Documentation	98
5.24.3.1	_delete	98
5.24.3.2	getHeaderValue	98
5.24.3.3	loadMagnetizationCurve	99
5.24.3.4	readConfigFromHeader	99
5.24.3.5	saveMagnetizationCurve	99
5.24.3.6	setFilenameSuffix	99
5.24.4	Member Data Documentation	99
5.24.4.1	fileId_csv	99
5.24.4.2	fileId_vsmdat	100
5.24.4.3	filename	100
5.24.4.4	filename_csv	100
5.24.4.5	filename_vsmdat	100
5.25	NetIO Class Reference	100
5.25.1	Constructor & Destructor Documentation	102
5.25.1.1	NetIO	102
5.25.2	Member Function Documentation	102
5.25.2.1	_delete	102
5.25.2.2	processCommands	102
5.25.2.3	processReadbuffer	102
5.25.2.4	sendData	102
5.25.2.5	startConnection	102
5.25.2.6	stopConnection	102
5.25.2.7	timedReceive	102
5.25.2.8	waitForAnswerTo	102
5.25.3	Member Data Documentation	102
5.25.3.1	await_respond	102
5.25.3.2	Connection	102
5.25.3.3	ip_addr	102
5.25.3.4	port	102
5.25.3.5	read_buffer	102

5.25.3.6	ReceiveTimer	102
5.25.3.7	updateInterval	102
5.26	NoController Class Reference	103
5.26.1	Constructor & Destructor Documentation	104
5.26.1.1	NoController	104
5.26.2	Member Function Documentation	104
5.26.2.1	_delete	104
5.26.2.2	abortMeasurement	104
5.26.2.3	afterAllLoops	104
5.26.2.4	afterLoop	104
5.26.2.5	beforeAllLoops	104
5.26.2.6	beforeLoop	104
5.26.2.7	cyclicTimer	104
5.26.2.8	getValue	104
5.26.2.9	setupTimer	104
5.26.2.10	setValue	104
5.26.2.11	stopTimer	104
5.26.2.12	timedFunction	104
5.26.3	Member Data Documentation	104
5.26.3.1	controller_setpoint	104
5.26.3.2	controller_value	105
5.26.3.3	CyclicTimer	105
5.26.3.4	filename_suffix	105
5.27	PhyEditBox Class Reference	105
5.27.1	Detailed Description	106
5.27.2	Constructor & Destructor Documentation	106
5.27.2.1	PhyEditBox	106
5.27.3	Member Function Documentation	107
5.27.3.1	_delete	107
5.27.3.2	checkVal	107
5.27.3.3	checkValue	107
5.27.3.4	getValue	107
5.27.3.5	internal_Callback	107
5.27.3.6	setValue	107
5.27.4	Member Data Documentation	107
5.27.4.1	cfg_obj	107
5.27.4.2	cfg_varname	107
5.27.4.3	default_val	107
5.27.4.4	hObject	107
5.27.4.5	max_val	107

5.27.4.6	min_val	108
5.27.4.7	ori_Callback	108
5.27.4.8	unit_str	108
5.27.4.9	value	108
5.28	RangeTable Class Reference	108
5.28.1	Detailed Description	111
5.28.2	Constructor & Destructor Documentation	111
5.28.2.1	RangeTable	111
5.28.3	Member Function Documentation	111
5.28.3.1	checkFields	111
5.28.3.2	checkFrameSize	111
5.28.3.3	checkValue	112
5.28.3.4	createSweepTable	112
5.28.3.5	fillTable	112
5.28.3.6	firstSweep	112
5.28.3.7	generateFieldArray	112
5.28.3.8	getData	112
5.28.3.9	getFieldRange	113
5.28.3.10	getJavaHandle	113
5.28.3.11	internal_Callback	113
5.28.3.12	internal_KeyPress_Callback	113
5.28.3.13	nextSweep	113
5.28.3.14	setData	113
5.28.3.15	setTable	113
5.28.3.16	stopSweep	114
5.28.4	Member Data Documentation	114
5.28.4.1	active_table	114
5.28.4.2	autoextend	114
5.28.4.3	back_to_zero_field	114
5.28.4.4	cycles	114
5.28.4.5	data	114
5.28.4.6	dataformat_descriptor	114
5.28.4.7	dataformat_descriptor_read	114
5.28.4.8	dataformat_header	114
5.28.4.9	fieldarray	114
5.28.4.10	frame_autosize	114
5.28.4.11	from	115
5.28.4.12	got_resized	115
5.28.4.13	ignore_bounds	115
5.28.4.14	lastedit_row	115

5.28.4.15	maxval	115
5.28.4.16	minval	115
5.28.4.17	ori_Callback	115
5.28.4.18	ori_keypress_Callback	115
5.28.4.19	scrollbar_rows	115
5.28.4.20	scrollbar_size	115
5.28.4.21	step	115
5.28.4.22	sweep_idx	115
5.28.4.23	table	115
5.28.4.24	table_initialized	116
5.28.4.25	tableJavaHandle	116
5.28.4.26	to	116
5.28.4.27	unit_str	116
5.28.4.28	use_rangetable	116
5.28.4.29	use_unofficial_java	116
5.29	Singleton Class Reference	116
5.29.1	Member Function Documentation	117
5.29.1.1	getSingletonData	117
5.29.1.2	instance	117
5.29.1.3	setSingletonData	117
5.29.2	Member Data Documentation	117
5.29.2.1	singletonData	117
5.30	SweepPanel Class Reference	117
5.30.1	Detailed Description	119
5.30.2	Constructor & Destructor Documentation	119
5.30.2.1	SweepPanel	119
5.30.3	Member Function Documentation	120
5.30.3.1	chkbox_Callback	120
5.30.3.2	chkboxChanged	120
5.30.3.3	chkboxZeroField_Callback	120
5.30.3.4	chkboxZeroFieldChanged	120
5.30.3.5	setCurrentCycle	120
5.30.3.6	showSweepSequence	120
5.30.4	Member Data Documentation	120
5.30.4.1	chkbox_advanced	120
5.30.4.2	chkbox_zerofield	120
5.30.4.3	ori_chkbox_callback	120
5.30.4.4	ori_chkboxZeroField_Callback	120
5.30.4.5	pan_advanced	120
5.30.4.6	pan_easy	120

5.30.4.7	RangeTable	121
5.30.4.8	str_current_cycle	121
5.30.4.9	SweepTable	121
5.30.4.10	txtCurrentCycle	121
5.30.4.11	txtCycles	121
5.31	SweepTable Class Reference	121
5.31.1	Detailed Description	124
5.31.2	Constructor & Destructor Documentation	124
5.31.2.1	SweepTable	124
5.31.3	Member Function Documentation	124
5.31.3.1	checkFields	124
5.31.3.2	checkFrameSize	125
5.31.3.3	checkValue	125
5.31.3.4	fillTable	125
5.31.3.5	firstSweep	125
5.31.3.6	generateFieldArray	125
5.31.3.7	getData	125
5.31.3.8	getFieldRange	126
5.31.3.9	getJavaHandle	126
5.31.3.10	internal_Callback	126
5.31.3.11	internal_KeyPress_Callback	126
5.31.3.12	nextSweep	126
5.31.3.13	setData	126
5.31.3.14	setTable	126
5.31.3.15	stopSweep	127
5.31.4	Member Data Documentation	127
5.31.4.1	active_table	127
5.31.4.2	autoextend	127
5.31.4.3	back_to_zero_field	127
5.31.4.4	cycles	127
5.31.4.5	data	127
5.31.4.6	dataformat_descriptor	127
5.31.4.7	dataformat_descriptor_read	127
5.31.4.8	dataformat_header	127
5.31.4.9	fieldarray	127
5.31.4.10	frame_autosize	127
5.31.4.11	from	128
5.31.4.12	got_resized	128
5.31.4.13	ignore_bounds	128
5.31.4.14	lastedit_row	128

5.31.4.15 maxval	128
5.31.4.16 minval	128
5.31.4.17 ori_Callback	128
5.31.4.18 ori_keypress_Callback	128
5.31.4.19 scrollbar_rows	128
5.31.4.20 scrollbar_size	128
5.31.4.21 step	128
5.31.4.22 sweep_idx	128
5.31.4.23 table	128
5.31.4.24 table_initialized	129
5.31.4.25 tableJavaHandle	129
5.31.4.26 to	129
5.31.4.27 unit_str	129
5.31.4.28 use_rangetable	129
5.31.4.29 use_unofficial_java	129
5.32 TempController Class Reference	129
5.32.1 Constructor & Destructor Documentation	131
5.32.1.1 TempController	131
5.32.2 Member Function Documentation	131
5.32.2.1 _delete	131
5.32.2.2 abortMeasurement	131
5.32.2.3 afterAllLoops	131
5.32.2.4 afterLoop	131
5.32.2.5 beforeAllLoops	131
5.32.2.6 beforeLoop	131
5.32.2.7 cyclicTimer	131
5.32.2.8 getValue	131
5.32.2.9 processCommands	131
5.32.2.10 processReadbuffer	132
5.32.2.11 sendData	132
5.32.2.12 setupTimer	132
5.32.2.13 setValue	132
5.32.2.14 startConnection	132
5.32.2.15 stopConnection	132
5.32.2.16 stopTimer	132
5.32.2.17 timedFunction	132
5.32.2.18 timedReceive	132
5.32.2.19 waitForAnswerTo	132
5.32.3 Member Data Documentation	132
5.32.3.1 await_respond	132

5.32.3.2	Connection	132
5.32.3.3	controller_setpoint	132
5.32.3.4	controller_value	132
5.32.3.5	CyclicTimer	132
5.32.3.6	filename_suffix	132
5.32.3.7	ip_addr	132
5.32.3.8	meas_continue	132
5.32.3.9	port	132
5.32.3.10	read_buffer	132
5.32.3.11	ReceiveTimer	132
5.32.3.12	updateInterval	133
6	File Documentation	135
6.1	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIn_digilent.cpp File Reference	135
6.2	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIn_mcc.cpp File Reference	135
6.3	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIO.cpp File Reference	135
6.4	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogOut_mcc.cpp File Reference	135
6.5	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/ArduinoIO.cpp File Reference	136
6.6	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/calibrateMagnet.cpp File Reference	136
6.6.1	Function Documentation	136
6.6.1.1	calibrateMagnet	136
6.7	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgBoolean.cpp File Reference	136
6.8	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgRange.cpp File Reference	136
6.9	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStr.cpp File Reference	136
6.10	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStrOrNum.cpp File Reference	137
6.11	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Config.cpp File Reference	137
6.12	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CustomController.cpp File Reference	137
6.13	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/dam3System.cpp File Reference	137
6.14	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/DataFilter.cpp File Reference	137
6.15	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/DataPlot.cpp File Reference	137
6.16	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/ErrorHandler.cpp File Reference	138
6.17	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/export2karamoke.cpp File Reference	138
6.18	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Fake_AnalogInAndOut.cpp File Reference	138
6.19	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/figCfgMain.cpp File Reference	138
6.19.1	Function Documentation	141
6.19.1.1	chkADCUseFeedback_Callback	141
6.19.1.2	chkCustomCtl_Callback	141
6.19.1.3	chkCustomValADC_Callback	141
6.19.1.4	chkDynamicYokeOffset_Callback	141

6.19.1.5	chkLockinSensitivity_Callback	141
6.19.1.6	chkSaveconfigMeas_Callback	141
6.19.1.7	chkTempCtl_Callback	141
6.19.1.8	cmdCalibrateMagnet_Callback	141
6.19.1.9	cmdCfgAdvanced_Callback	141
6.19.1.10	cmdCfgHardware_Callback	141
6.19.1.11	cmdManualField_Callback	141
6.19.1.12	cmdMeasureWindow_Callback	141
6.19.1.13	cmdSaveDescr_Callback	141
6.19.1.14	cmdSaveSettings_Callback	141
6.19.1.15	cmdTestADC_Callback	141
6.19.1.16	cmdViewCalib_Callback	141
6.19.1.17	cmdViewDataFile_Callback	142
6.19.1.18	cmdViewFieldoffset_Callback	142
6.19.1.19	disableGUI	142
6.19.1.20	enableGUI	142
6.19.1.21	figCfgMain	142
6.19.1.22	figCfgMain_CloseRequestFcn	142
6.19.1.23	figCfgMain_OpeningFcn	142
6.19.1.24	figCfgMain_OutputFcn	142
6.19.1.25	IstADCCardID_Callback	142
6.19.1.26	IstADCCardID_CreateFcn	142
6.19.1.27	IstADCDriver_Callback	142
6.19.1.28	IstADCDriver_CreateFcn	142
6.19.1.29	IstDACCardID_Callback	142
6.19.1.30	IstDACCardID_CreateFcn	142
6.19.1.31	IstDACDriver_Callback	142
6.19.1.32	IstDACDriver_CreateFcn	142
6.19.1.33	IstFilter_Callback	142
6.19.1.34	IstFilter_CreateFcn	142
6.19.1.35	IstMFieldCalib_Callback	142
6.19.1.36	IstMFieldCalib_CreateFcn	142
6.19.1.37	mnuAdvanced_Callback	143
6.19.1.38	mnuConfigView_Callback	143
6.19.1.39	mnuHardware_Callback	143
6.19.1.40	mnuLoadFrom_Callback	143
6.19.1.41	mnuLoadSave_Callback	143
6.19.1.42	mnuNormal_Callback	143
6.19.1.43	mnuSave_Callback	143
6.19.1.44	mnuSaveAs_Callback	143

6.19.1.45 onoff	143
6.19.1.46 txtADCCardID_Callback	143
6.19.1.47 txtADCCardID_CreateFcn	143
6.19.1.48 txtADCCustomValCalib_Callback	143
6.19.1.49 txtADCCustomValCalib_CreateFcn	143
6.19.1.50 txtADCCustomValCh_Callback	143
6.19.1.51 txtADCCustomValCh_CreateFcn	143
6.19.1.52 txtADCDriver_Callback	144
6.19.1.53 txtADCDriver_CreateFcn	144
6.19.1.54 txtADCduration_Callback	144
6.19.1.55 txtADCduration_CreateFcn	144
6.19.1.56 txtADCFeedbackCh_Callback	144
6.19.1.57 txtADCFeedbackCh_CreateFcn	144
6.19.1.58 txtADCHallCh_Callback	144
6.19.1.59 txtADCHallCh_CreateFcn	144
6.19.1.60 txtADCHallVRange_Callback	144
6.19.1.61 txtADCHallVRange_CreateFcn	144
6.19.1.62 txtADCPickupCh_Callback	144
6.19.1.63 txtADCPickupCh_CreateFcn	144
6.19.1.64 txtADCPickupQCh_Callback	144
6.19.1.65 txtADCPickupQCh_CreateFcn	144
6.19.1.66 txtADCPickupVRange_Callback	144
6.19.1.67 txtADCPickupVRange_CreateFcn	144
6.19.1.68 txtADCsamples_Callback	144
6.19.1.69 txtADCsamples_CreateFcn	144
6.19.1.70 txtADCsamplescont_Callback	144
6.19.1.71 txtADCsamplescont_CreateFcn	144
6.19.1.72 txtDACCardID_Callback	144
6.19.1.73 txtDACCardID_CreateFcn	144
6.19.1.74 txtDACDriver_Callback	144
6.19.1.75 txtDACDriver_CreateFcn	144
6.19.1.76 txtDACPowerSupplCh_Callback	144
6.19.1.77 txtDACPowerSupplCh_CreateFcn	144
6.19.1.78 txtDACramp_Callback	144
6.19.1.79 txtDACramp_CreateFcn	144
6.19.1.80 txtDACsamplerate_Callback	145
6.19.1.81 txtDACsamplerate_CreateFcn	145
6.19.1.82 txtFileExt_Callback	145
6.19.1.83 txtFileExt_CreateFcn	145
6.19.1.84 txtHallFactor_Callback	145

6.19.1.85 txtHallFactor_CreateFcn	145
6.19.1.86 txtLargeStep_Callback	145
6.19.1.87 txtLargeStep_CreateFcn	145
6.19.1.88 txtLockin_Callback	145
6.19.1.89 txtLockin_CreateFcn	145
6.19.1.90 txtmcal_maxfield_Callback	145
6.19.1.91 txtmcal_maxfield_CreateFcn	145
6.19.1.92 txtmcal_minfield_Callback	145
6.19.1.93 txtmcal_minfield_CreateFcn	145
6.19.1.94 txtmcal_Vmax_Callback	145
6.19.1.95 txtmcal_Vmax_CreateFcn	145
6.19.1.96 txtmcal_Vmin_Callback	145
6.19.1.97 txtmcal_Vmin_CreateFcn	145
6.19.1.98 txtmcal_Vstep_Callback	145
6.19.1.99 txtmcal_Vstep_CreateFcn	145
6.19.1.100 txtMFieldDescr_Callback	145
6.19.1.101 txtMFieldDescr_CreateFcn	145
6.19.1.102 txtPauseBetween_Callback	145
6.19.1.103 txtPauseBetween_CreateFcn	145
6.19.1.104 txtPauseLargeStep_Callback	145
6.19.1.105 txtPauseLargeStep_CreateFcn	145
6.19.1.106 txtTimePerDatapoint_Callback	145
6.19.1.107 txtTimePerDatapoint_CreateFcn	145
6.19.1.108 updateTimePerDataPoint	146
6.20 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/figFieldControl.cpp File Reference	147
6.20.1 Function Documentation	148
6.20.1.1 chkWriteProtect_Callback	148
6.20.1.2 cmdAdjustHall_Callback	148
6.20.1.3 cmdAdjustHallOffset_Callback	148
6.20.1.4 cmdAdjustVOut_Callback	148
6.20.1.5 cmdDegauss_Callback	148
6.20.1.6 cmdDn_Callback	149
6.20.1.7 cmdDnDn_Callback	149
6.20.1.8 cmdFDn_Callback	149
6.20.1.9 cmdFDnDn_Callback	149
6.20.1.10 cmdFUp_Callback	149
6.20.1.11 cmdFUpUp_Callback	149
6.20.1.12 cmdFZero_Callback	149
6.20.1.13 cmdInterrupt_Callback	149
6.20.1.14 cmdSetField_Callback	149

6.20.1.15 cmdSetVOut_Callback	149
6.20.1.16 cmdUp_Callback	149
6.20.1.17 cmdUpUp_Callback	149
6.20.1.18 cmdZero_Callback	150
6.20.1.19 disableButtons	150
6.20.1.20 driveMagnet	150
6.20.1.21 driveMagnet_relative	150
6.20.1.22 enableButtons	150
6.20.1.23 figFieldControl	150
6.20.1.24 figFieldControl_OpeningFcn	150
6.20.1.25 figFieldControl_OutputFcn	151
6.20.1.26 figFieldCtl_CloseRequestFcn	151
6.20.1.27 setField	151
6.20.1.28 setField_relative	151
6.20.1.29 txtCustomADC_Callback	151
6.20.1.30 txtCustomADC_CreateFcn	151
6.20.1.31 txtField_Callback	151
6.20.1.32 txtField_CreateFcn	151
6.20.1.33 txtFieldSetpoint_Callback	151
6.20.1.34 txtFieldSetpoint_CreateFcn	151
6.20.1.35 txtHallFct_Callback	151
6.20.1.36 txtHallFct_CreateFcn	151
6.20.1.37 txtHallOff_Callback	151
6.20.1.38 txtHallOff_CreateFcn	151
6.20.1.39 txtHallVolt_Callback	151
6.20.1.40 txtHallVolt_CreateFcn	151
6.20.1.41 txtPickup_Callback	151
6.20.1.42 txtPickup_CreateFcn	151
6.20.1.43 txtPickupQ_Callback	151
6.20.1.44 txtPickupQ_CreateFcn	151
6.20.1.45 txtVOut_Callback	152
6.20.1.46 txtVOut_CreateFcn	152
6.20.1.47 txtVOutNow_Callback	152
6.20.1.48 txtVOutNow_CreateFcn	152
6.20.1.49 updateFields	152
6.21 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/fileReadSample.cpp File Reference	152
6.22 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/findjobj.cpp File Reference	152
6.22.1 Function Documentation	153
6.22.1.1 btCheckFex_Callback	153
6.22.1.2 btExport_Callback	153

6.22.1.3	btFocus_Callback	153
6.22.1.4	btInspect_Callback	153
6.22.1.5	btRefresh_Callback	153
6.22.1.6	btWebsite_Callback	153
6.22.1.7	cbHideStdCbs_Callback	153
6.22.1.8	charizeData	153
6.22.1.9	checkVersion	153
6.22.1.10	createJButton	153
6.22.1.11	disableDbstopError	153
6.22.1.12	dispError	153
6.22.1.13	expandNode	154
6.22.1.14	findjobj	154
6.22.1.15	flashComponent	154
6.22.1.16	getCbsData	154
6.22.1.17	getChildrenNodes	154
6.22.1.18	getCurrentFigure	154
6.22.1.19	getLabelsJavaPos	154
6.22.1.20	getNodeName	154
6.22.1.21	getNodeTitleStr	154
6.22.1.22	getNumMenuComponents	154
6.22.1.23	getPropsHtml	154
6.22.1.24	getRalativeDivlocation	154
6.22.1.25	getRootPanel	154
6.22.1.26	getTopSelectedObject	154
6.22.1.27	getTreeData	154
6.22.1.28	getXY	154
6.22.1.29	iff	154
6.22.1.30	ischar	154
6.22.1.31	menuRemoveItem	154
6.22.1.32	nodeExpanded	154
6.22.1.33	nodeSelected	154
6.22.1.34	paramSupplied	154
6.22.1.35	presentObjectTree	154
6.22.1.36	processArgs	154
6.22.1.37	processClassArgs	154
6.22.1.38	processDepthArgs	154
6.22.1.39	processPositionArgs	154
6.22.1.40	processPrintArgs	154
6.22.1.41	processPropertyArgs	155
6.22.1.42	processSizeArgs	155

6.22.1.43	removeDuplicateNode	155
6.22.1.44	requestFocus	155
6.22.1.45	resizeImg	155
6.22.1.46	restoreDbstopError	155
6.22.1.47	revertCbTableModification	155
6.22.1.48	setIconSize	155
6.22.1.49	setProp	155
6.22.1.50	setTreeContextMenu	155
6.22.1.51	setTreeNodeIcon	155
6.22.1.52	stripStdCbs	155
6.22.1.53	tbCallbacksChanged	155
6.22.1.54	traverseContainer	155
6.22.1.55	traverseHGContainer	155
6.22.1.56	treeMouseMovedCallback	155
6.22.1.57	treeMousePressedCallback	155
6.22.1.58	updateNodeTooltip	155
6.22.1.59	warnInvisible	155
6.23	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/fMPlot.cpp File Reference	155
6.23.1	Function Documentation	157
6.23.1.1	chkAdvanced_Callback	157
6.23.1.2	chkCustomCtl_Callback	157
6.23.1.3	chkLoopPause_Callback	157
6.23.1.4	chkShowAverage_Callback	157
6.23.1.5	chkShowPoints_Callback	157
6.23.1.6	chkShowQuadrature_Callback	157
6.23.1.7	chkTemperatureCtl_Callback	157
6.23.1.8	chkZeroField_Callback	157
6.23.1.9	cmdCreateSequenceAndMeasure_Callback	157
6.23.1.10	cmdMeasure_Callback	157
6.23.1.11	cmdStopMeasurement_Callback	157
6.23.1.12	figMPlot_ButtonDownFcn	157
6.23.1.13	figMPlot_CloseRequestFcn	157
6.23.1.14	figMPlot_KeyPressFcn	157
6.23.1.15	figMPlot_SizeChangedFcn	157
6.23.1.16	figMPlot_WindowButtonDownFcn	157
6.23.1.17	fMPlot	157
6.23.1.18	fMPlot_OpeningFcn	158
6.23.1.19	fMPlot_OutputFcn	158
6.23.1.20	grpAdvancedSweeps_CreateFcn	158
6.23.1.21	lstLoopPlots_Callback	158

6.23.1.22	lstLoopPlots_CreateFcn	158
6.23.1.23	resizeGUI	158
6.23.1.24	resizeOriginal	158
6.23.1.25	setUnitNormalized	158
6.23.1.26	setUnitPixels	158
6.23.1.27	togglebutton1_Callback	158
6.23.1.28	txtComment_Callback	158
6.23.1.29	txtComment_CreateFcn	158
6.23.1.30	txtCustomADC_Callback	158
6.23.1.31	txtCustomADC_CreateFcn	158
6.23.1.32	txtCycles_Callback	158
6.23.1.33	txtCycles_CreateFcn	158
6.23.1.34	txtEasyStep_Callback	158
6.23.1.35	txtEasyStep_CreateFcn	158
6.23.1.36	txtEasyTo_Callback	158
6.23.1.37	txtEasyTo_CreateFcn	158
6.23.1.38	txtFieldOut_Callback	159
6.23.1.39	txtFieldOut_CreateFcn	159
6.23.1.40	txtFieldSetpoint_Callback	159
6.23.1.41	txtFieldSetpoint_CreateFcn	159
6.23.1.42	txtLockin_Callback	159
6.23.1.43	txtLockin_CreateFcn	159
6.23.1.44	txtPickup_Callback	159
6.23.1.45	txtPickup_CreateFcn	159
6.23.1.46	txtPickupQ_Callback	159
6.23.1.47	txtPickupQ_CreateFcn	159
6.23.1.48	txtTempManual_Callback	159
6.23.1.49	txtTempManual_CreateFcn	159
6.23.1.50	txtVOut_Callback	159
6.23.1.51	txtVOut_CreateFcn	159
6.24	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/GuiCfg.cpp File Reference	159
6.25	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/HWController.cpp File Reference	159
6.26	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/InterlockGUIElements.cpp File Reference	159
6.27	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/LevelBar.cpp File Reference	160
6.28	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/M3System.cpp File Reference	160
6.29	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MagField.cpp File Reference	160
6.30	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MagnetizationCurve.cpp File Reference	160
6.31	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/magnetizationLoop.cpp File Reference	160

6.31.1	Function Documentation	160
6.31.1.1	magnetizationLoop	160
6.32	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Measurement.cpp File Reference	160
6.33	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MeasurementFile.cpp File Reference	161
6.34	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/NetIO.cpp File Reference	161
6.35	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/NoController.cpp File Reference	161
6.36	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/PhyEditBox.cpp File Reference	161
6.37	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/prepareDocu.cpp File Reference	161
6.37.1	Function Documentation	162
6.37.1.1	postprocessCPP	162
6.37.1.2	prepareDocu	162
6.37.1.3	prepareMFiles	162
6.37.1.4	recursdir	162
6.38	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/RangeTable.cpp File Reference	162
6.39	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/show_stddev.cpp File Reference	162
6.40	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Singleton.cpp File Reference	162
6.41	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/sinum.cpp File Reference	162
6.41.1	Function Documentation	162
6.41.1.1	sinum	162
6.42	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/sipre.cpp File Reference	162
6.42.1	Function Documentation	163
6.42.1.1	sipre	163
6.43	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/SweepPanel.cpp File Reference	163
6.44	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/SweepTable.cpp File Reference	163
6.45	C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/TempController.cpp File Reference	163

Chapter 1

daM3System code documentation

1.1 Introduction

daM3System stands for data acquisition magnetic moments and MOKE.

This software is designed for controlling VSM and MOKE setups. Please find more abstract documentation here: (</doc/ppt>). You need MATLAB® and the Data Acquisition Toolbox™ and optional the Statistics Toolbox™ in order to compile this software.

Files

```
\VSM-Prog

+PkgAdvGUI\           : Class package
+PkgAnalogIO\         : Class package
+PkgHardwExt\         : Class package

@Config\             : Class in separate directory
@DataFilter\          : Class in separate directory
@ErrorHandler\        : Class in separate directory
@MagField\            : Class in separate directory
@MagnetizationCurve\  : Class in separate directory
@Measurement\         : Class in separate directory
@MeasurementFile\     : Class in separate directory
daM3System\           : Class in separate directory
daM3System_resources\ : Class in separate directory
doc\                  : Code documentation
doxygen\              : Temporary path for doxygen generation
.git\                 : GitHub distributed code revision

.gitattributes        : GitHub
.gitignore            : GitHub
Doxyfile              : Doxygen configuration file
daM3System.m          : Application main file
daM3System.prj        : Compiler project
daqtest.txt           : Hardware DAQ test report
error_log.mat         : ErrorHandler log file (track some exceptions using ErrorHandler class)
fake_pickup.dn.dat    : Data file for Fake_AnalogInAndOut hardware emulation
fake_pickup_dn.dat    : Data file for Fake_AnalogInAndOut hardware emulation
fake_pickup_up.dat    : Data file for Fake_AnalogInAndOut hardware emulation
figCfgMain.fig        : Main configuration window
figCfgMain.m          : Main configuration window
figFieldControl.fig   : Manual field control and hall probe calibration window
figFieldControl.m     : Manual field control and hall probe calibration window
fileReadSample.m      : Example for loading measurement data files for further processing throu
Matlab
findjob.m             : Access Java handles
findjob_license.txt   : Access Java handles licence file from findjobj-author
fMPlot.fig            : Measurement setup window
fMPlot.m              : Measurement setup window
gen_cpp.bat           : Doxygen: generate cpp header files
GuiCfg.m              : Central GUI configuration for VSM and MOKE (string table and some
default values)
Logo.png              : Splash
Logo3_k.png           : Icon
mag_curce_sample.dat  : Data file for Fake_AnalogInAndOut hardware emulation
MFieldCalibration_1.mat : Magnet calibration file
MFieldCalibration_2.mat : Magnet calibration file
MFieldCalibration_3.mat : Magnet calibration file
MFieldCalibration_4.mat : Magnet calibration file
```

prepareDocu.m	: Prepare code for Doxygen
Singleton.m	: Singleton design pattern
VSM_Config.mat	: Default configuration file

All product and company names are trademarks TM or registered ® trademarks of their respective holders. Use of them does not imply any affiliation with or endorsement by them.

Chapter 2

Hierarchical Index

2.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

AnalogIn_digilent	9
AnalogIn_mcc	11
AnalogIO	17
AnalogOut_mcc	28
AnalogIO	17
ArduinoIO	33
CfgBoolean	35
CfgStrOrNum	43
CfgRange	37
CfgStr	41
DataFilter	61
DataPlot	62
ErrorHandler	68
Fake_AnalogInAndOut	69
GuiCfg	75
hgsetget	77
MagnetizationCurve	90
HWController	78
CustomController	59
NoController	103
TempController	129
InterlockGUIElements	80
LevelBar	81
MagField	85
Measurement	95
MeasurementFile	96
NetIO	100
TempController	129
PhyEditBox	105
Singleton	116
Config	46
SweepPanel	117
SweepTable	121
RangeTable	108

Chapter 3

Class Index

3.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

AnalogIn_diligent	
AnalogIO class	9
AnalogIn_mcc	
Handles analog input for cards using DAQ-toolbox legacy interface	11
AnalogIO	
Anaalog In/Out implementation using DAQ-Toolbox legacy interface (tested with Measurement	
Computing and National Instruments Hardware)	17
AnalogOut_mcc	
AnalogIO class	28
ArduinoIO	33
CfgBoolean	
Extends a check box in order to automatically update the corresponding Config parameter	35
CfgRange	
Extends a text box in order to accept numeric ranges (2 dim arrays like "-10 10") and automati-	
cally update the corresponding Config parameter	37
CfgStr	
Extends a text box in order to accept string values and automatically update the corresponding	
Config parameter	41
CfgStrOrNum	
Extends a text box in order to accept numeric an string values and automatically update the	
corresponding Config parameter	43
Config	
Configuration class	46
CustomController	59
DataFilter	
DataFilter class (filters noise of acquired data)	61
DataPlot	
DataPlot class (graphical representation of data)	62
ErrorHandler	
Logs Exceptions and stores them in error_log.mat Most IO and other Exceptions are saved in	
error_log.mat and hence can be reviewed using static functions of the ErrorHandler class	68
Fake_AnalogInAndOut	
AnalogIO class	69
GuiCfg	75
hgsetget	77
HWController	78

InterlockGUIElements	
You can add GUI elements (eg. at window initialization) and enable/disable them all together anywhre else	80
LevelBar	
Extends an Axes object to act as a level bar	81
MagField	
MagField class; magnetic field calculation, calibration	85
MagnetizationCurve	
MagnetizationCurve class (represents data of entire magnetization curve)	90
Measurement	
Main functions for Measurement and Calibration; Contains main parts of measurement logic .	95
MeasurementFile	
MeasurementFile class (saves/loads measured data); See m for an example	96
NetIO	100
NoController	103
PhyEditBox	
Extends a text box in order to show SI values with units (kg, mA, ...)	105
RangeTable	
Extends a Uitable with PhyEditBox capabilities	108
Singleton	116
SweepPanel	
Holds and organizes all tables options and buttons for the sweep configuration inside the fMPlot window	117
SweepTable	
Extends a Uitable with PhyEditBox capabilities and and provides the sweep sequence	121
TempController	129

Chapter 4

File Index

4.1 File List

Here is a list of all files with brief descriptions:

C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIn_diligent.cpp	135
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIn_mcc.cpp	135
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIO.cpp	135
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogOut_mcc.cpp	135
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/ArduinoIO.cpp	136
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/calibrateMagnet.cpp	136
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgBoolean.cpp	136
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgRange.cpp	136
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStr.cpp	136
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStrOrNum.cpp	137
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Config.cpp	137
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CustomController.cpp	137
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/dam3System.cpp	137
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/DataFilter.cpp	137
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/DataPlot.cpp	137
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/ErrorHandler.cpp	138
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/export2karamoke.cpp	138
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Fake_AnalogInAndOut.cpp	138
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/figCfgMain.cpp	138
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/figFieldControl.cpp	147
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/fileReadSample.cpp	152
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/findjobj.cpp	152
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/fmPlot.cpp	155
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/GuiCfg.cpp	159
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/HWController.cpp	159
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/InterlockGUIElements.cpp	159
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/LevelBar.cpp	160
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/M3System.cpp	160
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MagField.cpp	160
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MagnetizationCurve.cpp	160
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/magnetizationLoop.cpp	160
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Measurement.cpp	160
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MeasurementFile.cpp	161
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/NetIO.cpp	161
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/NoController.cpp	161
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/PhyEditBox.cpp	161
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/prepareDocu.cpp	161
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/RangeTable.cpp	162

C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/show_stddev.cpp	162
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Singleton.cpp	162
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/sinum.cpp	162
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/sipre.cpp	162
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/SweepPanel.cpp	163
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/SweepTable.cpp	163
C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/TempController.cpp	163


Chapter 5

Class Documentation

5.1 AnalogIn_digilent Class Reference

[AnalogIO](#) class.

Collaboration diagram for AnalogIn_digilent:



Public Member Functions

- function [AnalogIn_digilent](#) ()
- function [_delete](#) (in obj)
- function [getHWinfoADC](#) (in obj)

- function [getInputRanges](#) (in obj)
- function [acquire_adc](#) (in obj)
- function [setPickupVrange](#) (in obj, in range)

Public Attributes

- Constant Property [ADCCard_device_id](#)
- Constant Property [ADCinputType](#)
- Constant Property [pickup_channel](#)
- Constant Property [pickup_q_channel](#)
- Constant Property [ini_pickup_vrange](#)
- Constant Property [ini_hallprobe_vrange](#)
- Property [adc_initialized](#)
- Property [pickup](#)
- Property [pickup_q](#)
- Property [hallprobe_volt](#)
- [DataFilter](#) filter

Protected Attributes

- Property [ai](#)
- Property [pickup_vrange](#)
- Property [hallprobe_vrange](#)
- Property [adc_samples_per_second](#)
- Property [adc_meas_duration](#)

Private Member Functions

- function [cleanUp](#) (in obj)

5.1.1 Detailed Description

[AnalogIO](#) class.

reads and filters from ADC-hardware + writes to DAC

5.1.2 Constructor & Destructor Documentation

5.1.2.1 function [AnalogIn_digilent](#) ()

5.1.3 Member Function Documentation

5.1.3.1 function [_delete](#) (in *obj*)

5.1.3.2 function [acquire_adc](#) (in *obj*)

5.1.3.3 function [cleanUp](#) (in *obj*) [private]

5.1.3.4 function [getHWinfoADC](#) (in *obj*)

5.1.3.5 function [getInputRanges](#) (in *obj*)

5.1.3.6 function setPickupVrange (in *obj*, in *range*)

5.1.4 Member Data Documentation

5.1.4.1 Property adc_initialized

5.1.4.2 Property adc_meas_duration [protected]

5.1.4.3 Property adc_samples_per_second [protected]

5.1.4.4 Constant Property ADCCard_device_id

5.1.4.5 Constant Property ADCinputType

5.1.4.6 Property ai [protected]

5.1.4.7 DataFilter filter

5.1.4.8 Property hallprobe_volt

5.1.4.9 Property hallprobe_vrange [protected]

5.1.4.10 Constant Property ini_hallprobe_vrange

5.1.4.11 Constant Property ini_pickup_vrange

5.1.4.12 Property pickup

5.1.4.13 Constant Property pickup_channel

5.1.4.14 Property pickup_q

5.1.4.15 Constant Property pickup_q_channel

5.1.4.16 Property pickup_vrange [protected]


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[AnalogIn_diligent.cpp](#)


5.2 AnalogIn_mcc Class Reference

Handles analog input for cards using DAQ-toolbox legacy interface.

Inheritance diagram for AnalogIn_mcc:



Collaboration diagram for AnalogIn_mcc:



Public Member Functions

- function [AnalogIn_mcc](#) ()
- function [_delete](#) (in obj)
- function [getHWinfoADC](#) (in obj)
Returns short hardware configuration report.
- function [getInputRanges](#) (in obj)
- function [acquire_adc](#) (in obj, in meas_duration)

- Captures all configured analog input channels.*

 - function [fetchBuffer_adc](#) (in obj)

Used with background capture mode; Processes analog data captured so far and clears data buffer.
- function [startBackground_adc](#) (in obj)

Starts continuous data acquisition in the background.
- function [stopBackground_adc](#) (in obj)

Stops continuous data acquisition in the background.
- function [setPickupVrange](#) (in obj, in range)
- function [setHallprobeVrange](#) (in obj, in range)

Public Attributes

- Constant Property [ADCinputType](#)

Input type to verify; use InstaCal to change.
- Constant Property [VOut_feedback_vrange](#)

has to be \geq magnet_powersupply_vrange !!
- Property [ai](#)
- Property [adc_initialized](#)

successfully initialized analog input
- Property [pickup](#)

Pickup/Signal ADC input returned by .acquire()
- Property [pickup_q](#)

Quadrature/Signal2 ADC input returned by .acquire()
- Property [hallprobe_volt](#)

Hall probe amplifier ADC input returned by .acquire()
- Property [customval_adc](#)

Custom ADC input returned by .acquire()
- Property [stddev_pickup](#)

Standard deviation of filtered signal input.
- Property [stddev_pickup_q](#)

Standard deviation of filtered signal2 input.
- Property [stddev_hallprobe_volt](#)

Standard deviation of magnetic field readings.
- [DataFilter](#) filter

Raw measurement data Filter. Filters data for a single data point.
- Property [background_acq](#)

continuously background acquisition is running

Protected Attributes

- Property [ADCCard_ID](#)

Unique Card identification used by Matlab ('DEV1',0,...)
- Property [pickup_channel](#)

channel of pickup/signal
- Property [pickup_q_channel](#)

channel of quadrature/another signal to capture
- Property [hallprobe_channel](#)

channel of hall probe amplifier
- Property [feedback_channel](#)

connected to V-Out for correct startup of DAC output; optional

- Property [customval_channel](#)
optional channel for additional analog data (eg. angle)
- Property [pickup_vrange](#)
Input range [Volt] for pickup/signal channels.
- Property [hallprobe_vrange](#)
Input range [Volt] for hallprobe channels.
- Property [adc_samples_per_second](#)
- Property [adc_meas_duration](#)
meas. duration of .acquire()
- Property [adc_meas_duration_snapshot](#)
duration of .acquire_snapshot()
- Property [initial_magnet_vout](#)
initial output of DAC

Private Member Functions

- function [cleanUp](#) (in *obj*)

5.2.1 Detailed Description

Handles analog input for cards using DAQ-toolbox legacy interface.
reads and filters from ADC-hardware

5.2.2 Constructor & Destructor Documentation

5.2.2.1 function [AnalogIn_mcc](#) ()

5.2.3 Member Function Documentation

5.2.3.1 function [_delete](#) (in *obj*)

5.2.3.2 function [acquire_adc](#) (in *obj*, in *meas_duration*)

Captures all configured analog input channels.

Parameters

<i>meas_duration</i>	[s] Measurement time for capturing analog data
----------------------	--

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.2.3.3 function [cleanUp](#) (in *obj*) [private]

5.2.3.4 function [fetchBuffer_adc](#) (in *obj*)

Used with background capture mode; Processes analog data captured so far and clears data buffer.

Data is stored in properties pickup, pickup_q, ...

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.2.3.5 function getHWinfoADC (in *obj*)

Returns short hardware configuration report.

Return values

<i>hw_info</i>	short hardware info as string
----------------	-------------------------------

5.2.3.6 function getInputRanges (in *obj*)

5.2.3.7 function setHallprobeVrange (in *obj*, in *range*)

5.2.3.8 function setPickupVrange (in *obj*, in *range*)

5.2.3.9 function startBackground_adc (in *obj*)

Starts continuous data acquisition in the background.

Sample rate is defined in [Config.adc_samplerate_cont](#) Use .fetchBuffer_adc to receive data.

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.2.3.10 function stopBackground_adc (in *obj*)

Stops continuous data acquisition in the background.

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.2.4 Member Data Documentation

5.2.4.1 Property adc_initialized

successfully initialized analog input

5.2.4.2 Property adc_meas_duration [protected]

meas. duration of .acquire()

5.2.4.3 Property adc_meas_duration_snapshot [protected]

duration of .acquire_snapshot()

5.2.4.4 Property adc_samples_per_second [protected]

5.2.4.5 Property ADCCard_ID [protected]

Unique Card identification used by Matlab ('DEV1',0,...)

5.2.4.6 Constant Property ADCinputType

Input type to verify; use InstaCal to change.

5.2.4.7 Property ai

5.2.4.8 Property background_acq

continuously background acquisition is running

5.2.4.9 Property customval_adc

Custom ADC input returned by .acquire()

5.2.4.10 Property customval_channel [protected]

optional channel for additional analog data (eg. angle)

5.2.4.11 Property feedback_channel [protected]

connected to V-Out for correct startup of DAC output; optional

5.2.4.12 DataFilter filter

Raw measurement data Filter. Filters data for a single data point.

5.2.4.13 Property hallprobe_channel [protected]

channel of hall probe amplifier

5.2.4.14 Property hallprobe_volt

Hall probe amplifier ADC input returned by .acquire()

5.2.4.15 Property hallprobe_vrange [protected]

Input range [Volt] for hallprobe channels.

5.2.4.16 Property initial_magnet_vout [protected]

initial output of DAC

5.2.4.17 Property pickup

Pickup/Signal ADC input returned by .acquire()

5.2.4.18 Property pickup_channel [protected]

channel of pickup/signal

5.2.4.19 Property pickup_q

Quadrature/Signal2 ADC input returned by .acquire()

5.2.4.20 Property pickup_q_channel [protected]

channel of quadrature/another signal to capture

5.2.4.21 Property pickup_vrange [protected]

Input range [Volt] for pickup/signal channels.

5.2.4.22 Property stddev_hallprobe_volt

Standard deviation of magnetic field readings.

5.2.4.23 Property stddev_pickup

Standard deviation of filtered signal input.

5.2.4.24 Property stddev_pickup_q

Standard deviation of filtered signal2 input.

5.2.4.25 Constant Property VOut_feedback_vrange

has to be \geq magnet_powersupply_vrange !!


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[AnalogIn_mcc.cpp](#)


5.3 AnalogIO Class Reference

Anaolog In/Out implementation using DAQ-Toolbox legacy interface (tested with [Measurement](#) Computing and National Instruments Hardware)

Inheritance diagram for AnalogIO:



Collaboration diagram for AnalogIO:



Public Member Functions

- function [AnalogIO](#) ()
Class constructor Establishes connection to Analog In/Out Hardware.
- function [acquire](#) (in obj)
Captures all configured analog channels Data is stored in properties pickup, pickup_q, ...
- function [acquire_snapshot](#) (in obj)
Captures analog data just for a short time interval ADC duration: [Config.adc_snapshot_duration](#) Data is stored in properties pickup, pickup_q, ...
- function [fetchBuffer](#) (in obj)
Used with background capture mode; Processes analog data captured so far and clears data buffer.

- function [getHWinfo](#) (in obj)
Returns short hardware configuration report.
- function [setMField](#) (in obj, in field, in monotone)
Drives magnet power supply using yoke calibration in order to reach a specified field strength.
- function [getReachTimeField](#) (in obj, in field)
Get time amount it takes to reach a specified field strength.
- function [resetFieldOffset](#) (in obj)
Clears recorded offset list (field setpoint vs.
- function [checkFieldOffset](#) (in obj)
Get time amount it takes to reach a specified field strength.
- function [constrain](#) (in obj, in inval, in range)
- function [failSafe](#) (in obj)
Try to sweep down field in case of an error; if connection to hardware is broken reestablishes connection and tries again.
- function [_delete](#) (in obj)
- function [getHWinfoADC](#) (in obj)
Returns short hardware configuration report.
- function [getInputRanges](#) (in obj)
- function [acquire_adc](#) (in obj, in meas_duration)
Captures all configured analog input channels.
- function [fetchBuffer_adc](#) (in obj)
Used with background capture mode; Processes analog data captured so far and clears data buffer.
- function [startBackground_adc](#) (in obj)
Starts continuous data acquisition in the background.
- function [stopBackground_adc](#) (in obj)
Stops continuous data acquisition in the background.
- function [setPickupVrange](#) (in obj, in range)
- function [setHallprobeVrange](#) (in obj, in range)
- function [_delete](#) (in obj)
- function [getHWinfoDAC](#) (in obj)
Returns short hardware configuration report.
- function [getOutputRanges](#) (in obj)
- function [getReachTime](#) (in obj, in volt_setpoint)
Get time amount it takes to reach a specified output voltage by DAC.
- function [reachMagnet](#) (in obj, in volt_setpoint)
Drives DAC output for magnet power supply.
- function [getDACout](#) (in obj)
Currently expected output voltage of DAC.

Static Public Member Functions

- static function [interrupt](#) (in cmd, in AnaIO)
Interrupts data acquisition; used for measurement abortion
`interrupt('set', AnaIO)` : Set AnaIO as current Analog in/out connection that should be able to be interrupted
`interrupt('ai')` : Interrupt current analog output
`interrupt('ao')` : Interrupt current analog input acquisition
`interrupt()` : Interrupt current analog input and output
`interrupt('reset')` : Remove current AnaIO object

Public Attributes

- Property [initialized](#)
- Property [magnetic_field](#)
[Oe] determined by hallprobe
- [MagField](#) [MField](#)
Functions for magnet and field.
- Property [yoke_offsets](#)
:field_sp :offset
- Constant Property [ADCinputType](#)
Input type to verify; use InstaCal to change.
- Constant Property [VOut_feedback_vrange](#)
has to be >= magnet_powersupply_vrange !!
- Property [ai](#)
- Property [adc_initialized](#)
sucessfully initialized analog input
- Property [pickup](#)
Pickup/Signal ADC input returned by .acquire()
- Property [pickup_q](#)
Quadrature/Signal2 ADC input returned by .acquire()
- Property [hallprobe_volt](#)
Hall probe amplifier ADC input returned by .acquire()
- Property [customval_adc](#)
Custom ADC input returned by .acquire()
- Property [stddev_pickup](#)
Standard deviation of filtered signal input.
- Property [stddev_pickup_q](#)
Standard deviation of filtered signal2 input.
- Property [stddev_hallprobe_volt](#)
Standard deviation of magnetic field readings.
- [DataFilter](#) [filter](#)
Raw measurement data Filter. Filters data for a single data point.
- Property [background_acq](#)
continuously background acquisition is running
- Constant Property [magnet_powersupply_vrange](#)
DACCard_ID = 0 DAC_samplerate = 1000 magnet_powersupply_channel = 0.
- Constant Property [putdata_minsamples](#)
Moke IO card requires at least 6 samples.
- Property [ao](#)
- Property [dac_samples_per_second](#)
current rate
- Property [dac_initialized](#)
- Property [magnet_vout](#)
- Property [magnetic_field_setpoint](#)

Protected Attributes

- Property [ADCCard_ID](#)
Unique Card identification used by Matlab ('DEV1',0,...)
- Property [pickup_channel](#)
channel of pickup/signal
- Property [pickup_q_channel](#)
channel of quadrature/another signal to capture
- Property [hallprobe_channel](#)
channel of hall probe amplifier
- Property [feedback_channel](#)
connected to V-Out for correct startup of DAC output; optional
- Property [customval_channel](#)
optional channel for additional analog data (eg. angle)
- Property [pickup_vrange](#)
Input range [Volt] for pickup/signal channels.
- Property [hallprobe_vrange](#)
Input range [Volt] for hallprobe channels.
- Property [adc_samples_per_second](#)
- Property [adc_meas_duration](#)
meas. duration of .acquire()
- Property [adc_meas_duration_snapshot](#)
duration of .acquire_snapshot()
- Property [initial_magnet_vout](#)
initial output of DAC

5.3.1 Detailed Description

Analog In/Out implementation using DAQ-Toolbox legacy interface (tested with [Measurement Computing](#) and National Instruments Hardware)

reads and filters from ADC-hardware + writes to DAC

5.3.2 Constructor & Destructor Documentation

5.3.2.1 function `AnalogIO ()`

Class constructor Establishes connection to Analog In/Out Hardware.

Returns

instance of [AnalogIO](#).

5.3.3 Member Function Documentation

5.3.3.1 function `_delete (in obj)` [inherited]

5.3.3.2 function `_delete (in obj)` [inherited]

5.3.3.3 function `acquire (in obj)`

Captures all configured analog channels Data is stored in properties pickup, pickup_q, ...

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.3.4 function `acquire_adc (in obj, in meas_duration)` [inherited]

Captures all configured analog input channels.

Parameters

<i>meas_duration</i>	[s] Measurement time for capturing analog data
----------------------	--

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.3.5 function `acquire_snapshot (in obj)`

Captures analog data just for a short time interval ADC duration: [Config.adc_snapshot_duration](#) Data is stored in properties pickup, pickup_q, ...

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.3.6 function `checkFieldOffset (in obj)`

Get time amount it takes to reach a specified field strength.

Parameters

<i>field</i>	Magnetic field setpoint
--------------	-------------------------

Return values

<i>rtime</i>	Time in seconds
--------------	-----------------

5.3.3.7 function `constrain (in obj, in inval, in range)`5.3.3.8 function `failSafe (in obj)`

Try to sweep down field in case of an error; if connection to hardware is broken reestablishes connection and tries again.

5.3.3.9 function `fetchBuffer (in obj)`

Used with background capture mode; Processes analog data captured so far and clears data buffer.

Data is stored in properties pickup, pickup_q, ...

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.3.10 function `fetchBuffer_adc (in obj)` [inherited]

Used with background capture mode; Processes analog data captured so far and clears data buffer.

Data is stored in properties `pickup`, `pickup_q`, ...

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.3.11 function `getDACout (in obj)` [inherited]

Currently expected output voltage of DAC.

Return values

<i>volt</i>	Currently expected output voltage of DAC
-------------	--

5.3.3.12 function `getHWinfo (in obj)`

Returns short hardware configuration report.

Return values

<i>hw_info</i>	short hardware info as string
----------------	-------------------------------

5.3.3.13 function `getHWinfoADC (in obj)` [inherited]

Returns short hardware configuration report.

Return values

<i>hw_info</i>	short hardware info as string
----------------	-------------------------------

5.3.3.14 function `getHWinfoDAC (in obj)` [inherited]

Returns short hardware configuration report.

Return values

<i>hw_info</i>	short hardware info as string
----------------	-------------------------------

5.3.3.15 function `getInputRanges (in obj)` [inherited]**5.3.3.16** function `getOutputRanges (in obj)` [inherited]**5.3.3.17** function `getReachTime (in obj, in volt_setpoint)` [inherited]

Get time amount it takes to reach a specified output voltage by DAC.

Parameters

<i>volt_setpoint</i>	DAC voltage setpoint
----------------------	----------------------

Return values

<i>rtime</i>	Time in seconds
--------------	-----------------

5.3.3.18 function getReachTimeField (in *obj*, in *field*)

Get time amount it takes to reach a specified field strength.

Parameters

<i>field</i>	Magnetic field setpoint
--------------	-------------------------

Return values

<i>rtime</i>	Time in seconds
--------------	-----------------

5.3.3.19 static function interrupt (in *cmd*, in *AnaIO*) [static]

Interrupts data acquisition; used for measurement abortion

`interrupt('set', AnaIO)` : Set AnaIO as current Analog in/out connection that should be able to be interrupted

`interrupt('ai')` : Interrupt current analog output

`interrupt('ao')` : Interrupt current analog input acquisition

`interrupt()` : Interrupt current analog input and output

`interrupt('reset')` : Remove current AnaIO object

Parameters

<i>cmd</i>	What to do (see description)
<i>AnaIO</i>	AnalogIO object

5.3.3.20 function reachMagnet (in *obj*, in *volt_setpoint*) [inherited]

Drives DAC output for magnet power supply.

Speed is defined in [Config.dac_ramp](#)

Parameters

<i>volt_setpoint</i>	DAC voltage setpoint
----------------------	----------------------

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.3.21 function resetFieldOffset (in *obj*)

Clears recorded offset list (field setpoint vs. measure field value)

5.3.3.22 function setHallprobeVrange (in *obj*, in *range*) [inherited]5.3.3.23 function setMField (in *obj*, in *field*, in *monotone*)

Drives magnet power supply using yoke calibration in order to reach a specified field strength.

Parameters

<i>field</i>	Magnetic field setpoint
<i>monotone</i>	Either 'monotone' or left open; Ensures monotone raise/fall of DAC_out for magnet power supply (regardless of any yoke calibration)

5.3.3.24 function `setPickupVrange (in obj, in range)` [inherited]

5.3.3.25 function `startBackground_adc (in obj)` [inherited]

Starts continuous data acquisition in the background.

Sample rate is defined in [Config.adc_samplerate_cont](#) Use `.fetchBuffer_adc` to receive data.

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.3.26 function `stopBackground_adc (in obj)` [inherited]

Stops continuous data acquisition in the background.

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.3.4 Member Data Documentation

5.3.4.1 Property `adc_initialized` [inherited]

successfully initialized analog input

5.3.4.2 Property `adc_meas_duration` [protected],[inherited]

meas. duration of `.acquire()`

5.3.4.3 Property `adc_meas_duration_snapshot` [protected],[inherited]

duration of `.acquire_snapshot()`

5.3.4.4 Property `adc_samples_per_second` [protected],[inherited]

5.3.4.5 Property `ADCCard_ID` [protected],[inherited]

Unique Card identification used by Matlab ('DEV1',0,...)

5.3.4.6 Constant Property `ADCInputType` [inherited]

Input type to verify; use `InstaCal` to change.

5.3.4.7 Property ai [inherited]

5.3.4.8 Property ao [inherited]

5.3.4.9 Property background_acq [inherited]

continuously background acquisition is running

5.3.4.10 Property customval_adc [inherited]

Custom ADC input returned by .acquire()

5.3.4.11 Property customval_channel [protected],[inherited]

optional channel for additional analog data (eg. angle)

5.3.4.12 Property dac_initialized [inherited]

5.3.4.13 Property dac_samples_per_second [inherited]

current rate

5.3.4.14 Property feedback_channel [protected],[inherited]

connected to V-Out for correct startup of DAC output; optional

5.3.4.15 DataFilter filter [inherited]

Raw measurement data Filter. Filters data for a single data point.

5.3.4.16 Property hallprobe_channel [protected],[inherited]

channel of hall probe amplifier

5.3.4.17 Property hallprobe_volt [inherited]

Hall probe amplifier ADC input returned by .acquire()

5.3.4.18 Property hallprobe_vrange [protected],[inherited]

Input range [Volt] for hallprobe channels.

5.3.4.19 Property initial_magnet_vout [protected],[inherited]

initial output of DAC

5.3.4.20 Property initialized

5.3.4.21 Constant Property magnet_powersupply_vrange [inherited]

DACCard_ID = 0 DAC_samplerate = 1000 magnet_powersupply_channel = 0.

5.3.4.22 Property `magnet_vout` [inherited]

5.3.4.23 Property `magnetic_field`

[Oe] determined by hallprobe

5.3.4.24 Property `magnetic_field_setpoint` [inherited]

5.3.4.25 **MagField** `MField`

Functions for magnet and field.

5.3.4.26 Property `pickup` [inherited]

Pickup/Signal ADC input returned by `.acquire()`

5.3.4.27 Property `pickup_channel` [protected],[inherited]

channel of pickup/signal

5.3.4.28 Property `pickup_q` [inherited]

Quadrature/Signal2 ADC input returned by `.acquire()`

5.3.4.29 Property `pickup_q_channel` [protected],[inherited]

channel of quadrature/another signal to capture

5.3.4.30 Property `pickup_vrange` [protected],[inherited]

Input range [Volt] for pickup/signal channels.

5.3.4.31 Constant Property `putdata_minsamples` [inherited]

Moke IO card requires at least 6 samples.

5.3.4.32 Property `stddev_hallprobe_volt` [inherited]

Standard deviation of magnetic field readings.

5.3.4.33 Property `stddev_pickup` [inherited]

Standard deviation of filtered signal input.

5.3.4.34 Property `stddev_pickup_q` [inherited]

Standard deviation of filtered signal2 input.

5.3.4.35 Constant Property VOut_feedback_vrange [inherited]

has to be \geq magnet_powersupply_vrange !!

5.3.4.36 Property yoke_offsets

:field_sp :offset


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[AnalogIO.cpp](#)


5.4 AnalogOut_mcc Class Reference

[AnalogIO](#) class.

Inheritance diagram for AnalogOut_mcc:



Collaboration diagram for AnalogOut_mcc:



Public Member Functions

- function [AnalogOut_mcc](#) ()
- function [_delete](#) (in obj)
- function [getHWinfoDAC](#) (in obj)
Returns short hardware configuration report.
- function [getOutputRanges](#) (in obj)
- function [getReachTime](#) (in obj, in volt_setpoint)
Get time amount it takes to reach a specified output voltage by DAC.
- function [reachMagnet](#) (in obj, in volt_setpoint)
Drives DAC output for magnet power supply.
- function [getDACout](#) (in obj)
Currently expected output voltage of DAC.

Public Attributes

- Constant Property [magnet_powersupply_vrange](#)
DACCard_ID = 0 DAC_samplerate = 1000 magnet_powersupply_channel = 0.
- Constant Property [putdata_minsamples](#)
MoKE IO card requires at least 6 samples.
- Property [ao](#)
- Property [dac_samples_per_second](#)
current rate
- Property [dac_initialized](#)
- Property [magnet_vout](#)
- Property [magnetic_field_setpoint](#)

Private Member Functions

- function [cleanUp](#) (in *obj*)
- function [constrain](#) (in *obj*, in *inval*, in *range*)

Private Attributes

- Property [magnet_volt_setpoint](#)
used by reachMagnet & getDACout
- Property [magnet_ramp](#)
"
- Property [magnet_powersupply_channel](#)
- Property [dac_samplerate](#)
Hz.
- Property [dac_ramp](#)
V/s.

5.4.1 Detailed Description

[AnalogIO](#) class.

reads and filters from ADC-hardware + writes to DAC

5.4.2 Constructor & Destructor Documentation

5.4.2.1 function [AnalogOut_mcc](#) ()

5.4.3 Member Function Documentation

5.4.3.1 function [_delete](#) (in *obj*)

5.4.3.2 function [cleanUp](#) (in *obj*) [private]

5.4.3.3 function [constrain](#) (in *obj*, in *inval*, in *range*) [private]

5.4.3.4 function [getDACout](#) (in *obj*)

Currently expected output voltage of DAC.

Return values

<i>volt</i>	Currently expected output voltage of DAC
-------------	--

5.4.3.5 function [getHWinfoDAC](#) (in *obj*)

Returns short hardware configuration report.

Return values

<i>hw_info</i>	short hardware info as string
----------------	-------------------------------

5.4.3.6 function getOutputRanges (in *obj*)

5.4.3.7 function getReachTime (in *obj*, in *volt_setpoint*)

Get time amount it takes to reach a specified output voltage by DAC.

Parameters

<i>volt_setpoint</i>	DAC voltage setpoint
----------------------	----------------------

Return values

<i>rtime</i>	Time in seconds
--------------	-----------------

5.4.3.8 function reachMagnet (in *obj*, in *volt_setpoint*)

Drives DAC output for magnet power supply.

Speed is defined in [Config.dac_ramp](#)

Parameters

<i>volt_setpoint</i>	DAC voltage setpoint
----------------------	----------------------

Return values

<i>ok</i>	1=success, 0=error
-----------	--------------------

5.4.4 Member Data Documentation

5.4.4.1 Property *ao*5.4.4.2 Property *dac_initialized*5.4.4.3 Property *dac_ramp* [private]

V/s.

5.4.4.4 Property *dac_samplerate* [private]

Hz.

5.4.4.5 Property *dac_samples_per_second*

current rate

5.4.4.6 Property *magnet_powersupply_channel* [private]5.4.4.7 Constant Property *magnet_powersupply_vrange*

DACCard_ID = 0 DAC_samplerate = 1000 magnet_powersupply_channel = 0.

5.4.4.8 Property *magnet_ramp* [private]

"

5.4.4.9 Property *magnet_volt_setpoint* [private]

used by reachMagnet & getDACout

5.4.4.10 Property magnet_vout

5.4.4.11 Property magnetic_field_setpoint

5.4.4.12 Constant Property putdata_minsamples


Moke IO card requires at least 6 samples.

The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogOut_mcc.cpp

5.5 ArduinoIO Class Reference

Collaboration diagram for ArduinoIO:



Public Member Functions

- function [ArduinoIO \(\)](#)
- function [_delete](#) (in obj)
- virtual function [processCommands](#) (in obj, in commands)

Public Attributes

- Constant Property [updateInterval](#)
seconds

Protected Member Functions

- function [startConnection](#) (in obj)
- function [stopConnection](#) (in obj)
- function [processReadbuffer](#) (in obj)

- function [sendData](#) (in *obj*, in *txt*)
- function [waitForAnswerTo](#) (in *obj*, in *txt*)
- function [timedReceive](#) (in *obj*, in *event*)

Protected Attributes

- Property [Connection](#)
serial object
- Property [port](#)
- Property [baud_rate](#)
- Property [read_buffer](#)
- Property [ReceiveTimer](#)
- Property [await_respond](#)
wait for receiving answer

5.5.1 Constructor & Destructor Documentation

5.5.1.1 function [ArduinoIO](#) ()

5.5.2 Member Function Documentation

5.5.2.1 function [_delete](#) (in *obj*)

5.5.2.2 virtual function [processCommands](#) (in *obj*, in *commands*) [virtual]

5.5.2.3 function [processReadbuffer](#) (in *obj*) [protected]

5.5.2.4 function [sendData](#) (in *obj*, in *txt*) [protected]

5.5.2.5 function [startConnection](#) (in *obj*) [protected]

5.5.2.6 function [stopConnection](#) (in *obj*) [protected]

5.5.2.7 function [timedReceive](#) (in *obj*, in *event*) [protected]

5.5.2.8 function [waitForAnswerTo](#) (in *obj*, in *txt*) [protected]

5.5.3 Member Data Documentation

5.5.3.1 Property [await_respond](#) [protected]

wait for receiving answer

5.5.3.2 Property [baud_rate](#) [protected]

5.5.3.3 Property [Connection](#) [protected]

serial object

5.5.3.4 Property port [protected]

5.5.3.5 Property read_buffer [protected]

5.5.3.6 Property ReceiveTimer [protected]

5.5.3.7 Constant Property updateInterval

seconds


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/ArduinoIO.cpp

5.6 CfgBoolean Class Reference

Extends a check box in order to automatically update the corresponding [Config](#) parameter.

Collaboration diagram for CfgBoolean:



Public Member Functions

- function [CfgBoolean](#) (in [checkbox](#), in [CfgObj](#), in [cfg_property](#))
- function [checkbox_Callback](#) (in [obj](#), in [hObject](#), in [eventdata](#), in [handles](#))
Internal callback function.
- function [checkboxChanged](#) (in [obj](#))
- function [getValue](#) (in [obj](#))
- function [setValue](#) (in [obj](#), in [val](#))

Protected Attributes

- Property [checkbox](#)
check box object used
- Property [ori_checkbox_callback](#)
original callback of check box

- Property [CfgObj](#)
config object
- Property [cfg_property](#)
parameter name of config
- Property [value](#)
value of check box

5.6.1 Detailed Description

Extends a check box in order to automatically update the corresponding [Config](#) parameter.

The Axes and Text objects can be placed using Mathworks's GUIDE(TM)

usage: `guiobj = PkgAdvGUI.CfgBoolean(hObject, cfg_obj, cfg_varname)` `hObject` Handle of text box object `cfg_obj` Optional; An object or structure that should get automatically updated on value changes of the [PhyEditBox](#) `cfg_varname` Optional; The property of `cfg_obj` which will be updated

The value can be changed with: `guiobj.SetValue(value);`

5.6.2 Constructor & Destructor Documentation

5.6.2.1 `function CfgBoolean (in chkbox, in CfgObj, in cfg_property)`

5.6.3 Member Function Documentation

5.6.3.1 `function chkbox_Callback (in obj, in hObject, in eventdata, in handles)`

Internal callback function.

5.6.3.2 `function chkboxChanged (in obj)`

5.6.3.3 `function getValue (in obj)`

5.6.3.4 `function setValue (in obj, in val)`

5.6.4 Member Data Documentation

5.6.4.1 **Property `cfg_property`** `[protected]`

parameter name of config

5.6.4.2 **Property `CfgObj`** `[protected]`

config object

5.6.4.3 **Property `chkbox`** `[protected]`

check box object used

5.6.4.4 **Property `ori_chkbox_callback`** `[protected]`

original callback of check box

5.6.4.5 Property value [protected]

value of check box


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgBoolean.cpp


5.7 CfgRange Class Reference

Extends a text box in order to accept numeric ranges (2 dim arrays like "-10 10") and automatically update the corresponding [Config](#) parameter.

Inheritance diagram for CfgRange:



Collaboration diagram for CfgRange:



Public Member Functions

- function [CfgRange](#) (in [hObject](#), in [default_val](#), in [cfg_obj](#), in [cfg_varname](#))
- function [_delete](#) (in [obj](#))
- function [internal_Callback](#) (in [obj](#), in [hObject](#), in [eventdata](#), in [handles](#))
Internal callback function for the text change event; Regular callback function will be called after this method.
- function [getValue](#) (in [obj](#))
- function [setValue](#) (in [obj](#), in [sval](#))
Interpret sval as string and as number.

Public Attributes

- Property [default_val](#)

Protected Member Functions

- function [checkValue](#) (in [obj](#), in [str](#), in [val](#))
Interpret and accept or reject entered text.

Protected Attributes

- Property [hObject](#)
text box object used
- Property [ori_Callback](#)
original callback of text box
- Property [value](#)
construed value of entered text
- Property [cfg_obj](#)
config object
- Property [cfg_varname](#)
parameter name of config

5.7.1 Detailed Description

Extends a text box in order to accept numeric ranges (2 dim arrays like "-10 10") and automatically update the corresponding [Config](#) parameter.

The Axes and Text objects can be placed using Mathworks's GUIDE(TM)

usage: `guiobj = PkgAdvGUI.CfgRange(hObject, default_val, cfg_obj, cfg_varname)`
[hObject](#) Handle of text box object [default_val](#) Default value; ignored whe [cfg_obj](#) is used [cfg_obj](#) Optional; An object or structure that should get automatically updated on value changes of the [PhyEditBox](#)
[cfg_varname](#) Optional; The property of [cfg_obj](#) which will be updated

The value can be changed with: `guiobj.setValue(value);`

5.7.2 Constructor & Destructor Documentation

5.7.2.1 function [CfgRange](#) (in [hObject](#), in [default_val](#), in [cfg_obj](#), in [cfg_varname](#))

5.7.3 Member Function Documentation

5.7.3.1 function `_delete` (in *obj*) [inherited]

5.7.3.2 function `checkValue` (in *obj*, in *str*, in *val*) [protected]

Interpret and accept or reject entered text.

Parameters

<i>str</i>	Entered text interpreted as string
<i>val</i>	Entered text interpreted as numeric value

Return values

<i>rval</i>	accepted value of accepted datatype or [] if entered text is rejected
-------------	---

5.7.3.3 function `getValue (in obj)` [inherited]

5.7.3.4 function `internal_Callback (in obj, in hObject, in eventdata, in handles)` [inherited]

Internal callback function for the text change event; Regular callback function will be called after this method.

5.7.3.5 function `setValue (in obj, in sval)` [inherited]

Interpret sval as string and as number.

Calls `checkValue` to verify entered value

Parameters

<i>sval</i>	entered text as string
-------------	------------------------

5.7.4 Member Data Documentation

5.7.4.1 Property `cfg_obj` [protected],[inherited]

config object

5.7.4.2 Property `cfg_varname` [protected],[inherited]

parameter name of config

5.7.4.3 Property `default_val` [inherited]

5.7.4.4 Property `hObject` [protected],[inherited]

text box object used

5.7.4.5 Property `ori_Callback` [protected],[inherited]

original callback of text box

5.7.4.6 Property `value` [protected],[inherited]

construed value of entered text


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgRange.cpp


5.8 CfgStr Class Reference

Extends a text box in order to accept string values and automatically update the corresponding [Config](#) parameter.

Inheritance diagram for CfgStr:



Collaboration diagram for CfgStr:



Public Member Functions

- function [CfgStr](#) (in [hObject](#), in [default_val](#), in [cfg_obj](#), in [cfg_varname](#))
- function [_delete](#) (in [obj](#))
- function [internal_Callback](#) (in [obj](#), in [hObject](#), in [eventdata](#), in [handles](#))
Internal callback function for the text change event; Regular callback function will be called after this method.
- function [getValue](#) (in [obj](#))
- function [setValue](#) (in [obj](#), in [sval](#))
Interpret sval as string and as number.

Public Attributes

- Property [default_val](#)

Protected Member Functions

- function [checkValue](#) (in *obj*, in *str*, in *val*)
Interpret and accept or reject entered text.

Protected Attributes

- Property [hObject](#)
text box object used
- Property [ori_Callback](#)
original callback of text box
- Property [value](#)
construed value of entered text
- Property [cfg_obj](#)
config object
- Property [cfg_varname](#)
parameter name of config

5.8.1 Detailed Description

Extends a text box in order to accept string values and automatically update the corresponding [Config](#) parameter.

The Axes and Text objects can be placed using Mathworks's GUIDE(TM)

usage: `guiobj = PkgAdvGUI.CfgStr(hObject, default_val, cfg_obj, cfg_varname)`
hObject Handle of text box object *default_val* Default value; ignored whe *cfg_obj* is used *cfg_obj* Optional; An object or structure that should get automatically updated on value changes of the [PhyEditBox](#)
cfg_varname Optional; The property of *cfg_obj* which will be updated

The value can be changed with: `guiobj.SetValue(value);`

5.8.2 Constructor & Destructor Documentation

5.8.2.1 function `CfgStr (in hObject, in default_val, in cfg_obj, in cfg_varname)`

5.8.3 Member Function Documentation

5.8.3.1 function `_delete (in obj)` [inherited]

5.8.3.2 function `checkValue (in obj, in str, in val)` [protected]

Interpret and accept or reject entered text.

Parameters

<i>str</i>	Entered text interpreted as string
<i>val</i>	Entered text interpreted as numeric value

Return values

<i>rval</i>	accepted value of accepted datatype or [] if entered text is rejected
-------------	---

5.8.3.3 function `getValue (in obj)` [inherited]

5.8.3.4 function `internal_Callback (in obj, in hObject, in eventdata, in handles)` [inherited]

Internal callback function for the text change event; Regular callback function will be called after this method.

5.8.3.5 function `setValue (in obj, in sval)` [inherited]

Interpret sval as string and as number.

Calls `checkValue` to verify entered value

Parameters

<i>sval</i>	entered text as string
-------------	------------------------

5.8.4 Member Data Documentation

5.8.4.1 Property `cfg_obj` [protected], [inherited]

config object

5.8.4.2 Property `cfg_varname` [protected], [inherited]

parameter name of config

5.8.4.3 Property `default_val` [inherited]

5.8.4.4 Property `hObject` [protected], [inherited]

text box object used

5.8.4.5 Property `ori_Callback` [protected], [inherited]

original callback of text box

5.8.4.6 Property `value` [protected], [inherited]

construed value of entered text


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStr.cpp


5.9 CfgStrOrNum Class Reference

Extends a text box in order to accept numeric an string values and automatically update the corresponding [Config](#) parameter.

Inheritance diagram for CfgStrOrNum:



Collaboration diagram for CfgStrOrNum:



Public Member Functions

- function [CfgStrOrNum](#) (in [hObject](#), in [default_val](#), in [cfg_obj](#), in [cfg_varname](#))
- function [_delete](#) (in [obj](#))
- function [internal_Callback](#) (in [obj](#), in [hObject](#), in [eventdata](#), in [handles](#))

Internal callback function for the text change event; Regular callback function will be called after this method.

- function [getValue](#) (in [obj](#))
- function [setValue](#) (in [obj](#), in [sval](#))

Interpret sval as string and as number.

Public Attributes

- Property [default_val](#)

Protected Member Functions

- function [checkValue](#) (in obj, in str, in val)
Interpret and accept or reject entered text.

Protected Attributes

- Property [hObject](#)
text box object used
- Property [ori_Callback](#)
original callback of text box
- Property [value](#)
construed value of entered text
- Property [cfg_obj](#)
config object
- Property [cfg_varname](#)
parameter name of config

5.9.1 Detailed Description

Extends a text box in order to accept numeric an string values and automatically update the corresponding [Config](#) parameter.

The Axes and Text objects can be placed using Mathworks's GUIDE(TM)

usage: `guiobj = PkgAdvGUI.CfgStrOrNum(hObject, default_val, cfg_obj, cfg_varname)`
 hObject Handle of text box object default_val Default value; ignored whe cfg_obj is used
 cfg_obj Optional; An object or structure that should get automatically updated on value changes of the [Phy↔](#)
[EditBox](#) cfg_varname Optional; The property of cfg_obj which will be updated

cfg_varname will be an integer if entered value (by user) is numeric or string otherwise

The value can be changed with: `guiobj.setValue(value);`

5.9.2 Constructor & Destructor Documentation

5.9.2.1 function `CfgStrOrNum (in hObject, in default_val, in cfg_obj, in cfg_varname)`

5.9.3 Member Function Documentation

5.9.3.1 function `_delete (in obj)`

5.9.3.2 function `checkValue (in obj, in str, in val)` [protected]

Interpret and accept or reject entered text.

Parameters

<i>str</i>	Entered text interpreted as string
<i>val</i>	Entered text interpreted as numeric value

Return values

<i>rval</i>	accepted value of accepted datatype or [] if entered text is rejected
-------------	---

5.9.3.3 function `getValue (in obj)`

5.9.3.4 function `internal_Callback (in obj, in hObject, in eventdata, in handles)`

Internal callback function for the text change event; Regular callback function will be called after this method.

5.9.3.5 function `setValue (in obj, in sval)`

Interpret *sval* as string and as number.

Calls `checkValue` to verify entered value

Parameters

<i>sval</i>	entered text as string
-------------	------------------------

5.9.4 Member Data Documentation

5.9.4.1 Property `cfg_obj` [`protected`]

config object

5.9.4.2 Property `cfg_varname` [`protected`]

parameter name of config

5.9.4.3 Property `default_val`

5.9.4.4 Property `hObject` [`protected`]

text box object used

5.9.4.5 Property `ori_Callback` [`protected`]

original callback of text box

5.9.4.6 Property `value` [`protected`]



construed value of entered text

The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStrOrNum.cpp

5.10 Config Class Reference

Configuration class.



- function **saveMe** (in cfg)

- function `saveMe` (in `cfg`)
Save configuration to default file.
- function `saveMeAs` (in `cfg`)
Save configuration to different file (with file dialog)
- function `hasChanged` (in `obj`)
Check if there are any unsaved parameters (changed window position will also apply)
- function `saveDialog` (in `cfg`, in `dialog_text`)
Message box asking if configuration shall be saved.
- function `getSingletonData` (in `obj`)
- function `setSingletonData` (in `obj`, in `singletonData`)

Static Public Member Functions

- static function [instance](#) (in loaded_obj)
Returns singleton instance of [Config](#) and loads parameters from file if needed.
- static function [loadFrom](#) ()
Loads configuration from different file (closes all windows)
- static function [getTest](#) ()
- virtual static [instance](#) ()

Public Attributes

- Constant Property [autoload_on_first_use](#)
automatically load [Config](#) on first access of [Config.instance\(\)](#)
- Constant Property [config_filename](#)
name of default config file
- Constant Property [sw_version](#)
- Constant Property [runstate_stopped](#)
state definition
- Constant Property [runstate_running](#)
state definition
- Constant Property [runstate_pleasestop](#)
state definition
- Constant Property [runstate_pleasestop_user](#)
state definition
- Constant Property [mplotview_std](#)
appearance of fmpplot window
- Constant Property [mplotview_viewcalibration](#)
view magnet calibration
- Constant Property [mplotview_magnetcalibration](#)
calibrate magnet
- Constant Property [mplotview_measure](#)
measure magnetization curve
- Constant Property [mplotview_viewdatafile](#)
view magnetization curve
- Constant Property [adc_snapshot_duration](#)
seconds; fixed to fit in Measurements timer interval
- Constant Property [err_cancel_measure](#)
Exceptions.
- Constant Property [err_IO_init](#)
- Property [daq_filter_idx](#)
nr of DAQ filter
- Property [datfile_path](#)
last used filepath to measurement files
- Property [use_tempcontroller](#)
whether to use the temperature controller HW extension class
- Property [use_customcontroller](#)
whether to use the custom controller HW extension class
- Property [file_ext](#)
file extension for data files
- Property [intermediate_datasaving](#)
save Loops after every completed Loop

- Property [pause_large_fieldstep](#)
seconds
- Property [large_field_step](#)
Oe.
- Property [pause_between_datapoints](#)
seconds
- Property [ADCDriver](#)
Hardware properties :
- Property [ADCCard_ID](#)
- Property [adc_samplerate](#)
samples/sec
- Property [adc_duration](#)
seconds
- Property [adc_samplerate_cont](#)
continuous background sample rate
- Property [adc_pickup_channel](#)
- Property [adc_pickup_q_channel](#)
- Property [adc_hallprobe_channel](#)
- Property [adc_feedback_channel](#)
- Property [adc_customval_channel](#)
- Property [adc_pickup_vrange](#)
- Property [adc_hallprobe_vrange](#)
- Property [adc_use_customval](#)
- Property [adc_customval_factor](#)
- Property [adc_use_feedback_channel](#)
feedback channel to get last VOut value
- Property [adc_lockin_sensitivity](#)
0: always start at zero field (VOut)
- Property [DACDriver](#)
- Property [DACCard_ID](#)
- Property [dac_magnet_powersupply_channel](#)
- Property [dac_samplerate](#)
Hz.
- Property [dac_ramp](#)
V/s.
- Property [check_lockin_sensitivity](#)
ask if lock-in sensitivity is set correctly
- Property [check_saveconfig_before_meas](#)
ask if config shall be saved before measurement starts
- Property [hallprobe_factor](#)
Oe/V.
- Property [hallprobe_offset](#)
V.
- Property [magnet_min_field](#)
Oe.
- Property [magnet_max_field](#)
Oe.
- Property [magnet_calibration_min_volt](#)
volt
- Property [magnet_calibration_max_volt](#)
volt
- Property [magnet_calibration_step_volt](#)

- volt*
- Property [magnet_calibration_file_nr](#)
1-9
- Property [magnet_dynamic_offset](#)
determine current offset due to yoke remanence
- Property [plot_show_quadrature](#)
show in signal2 data in plot window
- Property [plot_show_average](#)
show in average data in plot window
- Property [plot_show_points](#)
show data points in plot window instead of line
- [SweepTable](#) [SweepTable](#)
initalize sweep table, %
- [RangeTable](#) [RangeTable](#)
initalize range table, %
- Property [winpos_MPlot](#)
mFieldSweeps = struct('from', [], Window positions: Window position
- Property [winpos_MainCfg](#)
Window position.
- Property [winpos_FieldCtl](#)
Window position.
- Property [comment](#)
a comment for measurement files
- Property [test](#)
delme
- Property [runstate](#)
set to 0 to abort measurements
- Property [plot_window](#)
store plot window object
- Property [manual_temperature_val](#)
temperature set by user when no [HWController](#) is used
- Property [timestamp_start](#)
start of measurement
- Property [timestamp_end](#)
end of measurement
- Property [gui_loop_pause](#)
*gui helper variables (used as global variables)
pause button [CfgBoolean](#) object*
- Property [gui_yoke_offsets](#)
missmatch between setpoints and real field value
- Property [gui_fMPlot](#)
handles to opened window
- Property [gui_figFieldControl](#)
handles to opened window
- Property [gui_figCfgMain](#)
handles to opened window
- [LevelBar](#) [gui_bar_magnet_vout](#)
vout level bar (fMPlot window)
- [LevelBar](#) [gui_bar_magnet_field](#)
magn. field level bar (fMPlot window)

- [DataPlot](#) [gui_dataPlot](#)
data plot object (fMPlot window)
- [SweepPanel](#) [gui_SweepPanel](#)
sweep sequence (points to advanced sweep table in fMPlot)
- [PhyEditBox](#) [gui_txtFieldSetpoint](#)
numerical display (fMPlot window)
- [PhyEditBox](#) [gui_txtPickup](#)
numerical display (fMPlot window)
- [PhyEditBox](#) [gui_txtPickupQ](#)
numerical display (fMPlot window)
- [PhyEditBox](#) [gui_txtCustomADC](#)
numerical display (fMPlot window)
- Property [guiInterlock](#)
Disable (grey) GUI elements during DAQ.

Private Member Functions

- function [Config](#) ()

5.10.1 Detailed Description

Configuration class.

loads / saves parameters for measurement cycles

5.10.2 Constructor & Destructor Documentation

5.10.2.1 function [Config](#) () [private]

5.10.3 Member Function Documentation

5.10.3.1 function [getSingletonData](#) (in *obj*) [inherited]

5.10.3.2 static function [getTest](#) () [static]

5.10.3.3 function [hasChanged](#) (in *obj*)

Check if there are any unsaved parameters (changed window position will also apply)

Return values

<i>changed</i>	1 if parameters have changed; 0 otherwise
----------------	---

5.10.3.4 virtual static [instance](#) () [static],[virtual],[inherited]

5.10.3.5 static function [instance](#) (in *loaded_obj*) [static]

Returns singleton instance of [Config](#) and loads parameters from file if needed.

Return values

<i>obj</i>	Singleton instance of Config
------------	--

5.10.3.6 static function loadFrom () [static]

Loads configuration from different file (closes all windows)

Return values

<i>did_it</i>	1=new config has been loaded; 0 otherwise
---------------	---

5.10.3.7 function saveDialog (in *cfg*, in *dialog_text*)

Message box asking if configuration shall be saved.

Parameters

<i>dialog_text</i>	Optional. Text shown in question dialog
--------------------	---

5.10.3.8 function saveMe (in *cfg*)

Save configuration to default file.

5.10.3.9 function saveMeAs (in *cfg*)

Save configuration to different file (with file dialog)

5.10.3.10 function setSingletonData (in *obj*, in *singletonData*) [inherited]

5.10.4 Member Data Documentation

5.10.4.1 Property *adc_customval_channel*5.10.4.2 Property *adc_customval_factor*5.10.4.3 Property *adc_duration*

seconds

5.10.4.4 Property *adc_feedback_channel*5.10.4.5 Property *adc_hallprobe_channel*5.10.4.6 Property *adc_hallprobe_vrange*5.10.4.7 Property *adc_lockin_sensitivity*

0: always start at zero field (VOut)

5.10.4.8 Property `adc_pickup_channel`

5.10.4.9 Property `adc_pickup_q_channel`

5.10.4.10 Property `adc_pickup_vrange`

5.10.4.11 Property `adc_samplerate`

samples/sec

5.10.4.12 Property `adc_samplerate_cont`

continuous background sample rate

5.10.4.13 Constant Property `adc_snapshot_duration`

seconds; fixed to fit in Measurements timer interval

5.10.4.14 Property `adc_use_customval`

5.10.4.15 Property `adc_use_feedback_channel`

feedback channel to get last VOut value

5.10.4.16 Property `ADCCard_ID`

5.10.4.17 Property `ADCDriver`

Hardware properties :

5.10.4.18 Constant Property `autoload_on_first_use`

automatically load [Config](#) on first access of [Config.instance\(\)](#)

5.10.4.19 Property `check_lockin_sensitivity`

ask if lock-in sensitivity is set correctly

5.10.4.20 Property `check_saveconfig_before_meas`

ask if config shall be saved before measurement starts

5.10.4.21 Property `comment`

a comment for measurement files

5.10.4.22 Constant Property `config_filename`

name of default config file

5.10.4.23 **Property** `dac_magnet_powersupply_channel`

5.10.4.24 **Property** `dac_ramp`

V/s.

5.10.4.25 **Property** `dac_samplerate`

Hz.

5.10.4.26 **Property** `DACCard_ID`

5.10.4.27 **Property** `DACDriver`

5.10.4.28 **Property** `daq_filter_idx`

nr of DAQ filter

5.10.4.29 **Property** `datfile_path`

last used filepath to measurement files

5.10.4.30 **Constant Property** `err_cancel_measure`

Exceptions.

5.10.4.31 **Constant Property** `err_IO_init`

5.10.4.32 **Property** `file_ext`

file extension for data files

5.10.4.33 **LevelBar** `gui_bar_magnet_field`

magn. field level bar (fMPlot window)

5.10.4.34 **LevelBar** `gui_bar_magnet_vout`

vout level bar (fMPlot window)

5.10.4.35 **DataPlot** `gui_dataPlot`

data plot object (fMPlot window)

5.10.4.36 **Property** `gui_figCfgMain`

handles to opened window

5.10.4.37 **Property** `gui_figFieldControl`

handles to opened window

5.10.4.38 Property gui_fmPlot

handles to opened window

5.10.4.39 Property gui_loop_pause

gui helper variables (used as global variables)

pause button [CfgBoolean](#) object

5.10.4.40 SweepPanel gui_SweepPanel

sweep sequence (points to advanced sweep table in fmPlot)

5.10.4.41 PhyEditBox gui_txtCustomADC

numerical display (fmPlot window)

5.10.4.42 PhyEditBox gui_txtFieldSetpoint

numerical display (fmPlot window)

5.10.4.43 PhyEditBox gui_txtPickup

numerical display (fmPlot window)

5.10.4.44 PhyEditBox gui_txtPickupQ

numerical display (fmPlot window)

5.10.4.45 Property gui_yoke_offsets

missmatch between setpoints and real field value

5.10.4.46 Property guiInterlock

Disable (grey) GUI elements during DAQ.

5.10.4.47 Property hallprobe_factor

Oe/V.

5.10.4.48 Property hallprobe_offset

V.

5.10.4.49 Property intermediate_datasaving

save Loops after every completed Loop

5.10.4.50 Property large_field_step

Oe.

5.10.4.51 Property magnet_calibration_file_nr

1-9

5.10.4.52 Property magnet_calibration_max_volt

volt

5.10.4.53 Property magnet_calibration_min_volt

volt

5.10.4.54 Property magnet_calibration_step_volt

volt

5.10.4.55 Property magnet_dynamic_offset

determine current offset due to yoke remanence

5.10.4.56 Property magnet_max_field

Oe.

5.10.4.57 Property magnet_min_field

Oe.

5.10.4.58 Property manual_temperature_val

temperature set by user when no [HWController](#) is used

5.10.4.59 Constant Property mplotview_magnetcalibration

calibrate magnet

5.10.4.60 Constant Property mplotview_measure

measure magnetization curve

5.10.4.61 Constant Property mplotview_std

appearance of fmpplot window

5.10.4.62 Constant Property mplotview_viewcalibration

view magnet calibration

5.10.4.63 Constant Property mplotview_viewdatafile

view magnetization curve

5.10.4.64 Property pause_between_datapoints

seconds

5.10.4.65 Property pause_large_fieldstep

seconds

5.10.4.66 Property plot_show_average

show in average data in plot window

5.10.4.67 Property plot_show_points

show data points in plot window instead of line

5.10.4.68 Property plot_show_quadrature

show in signal2 data in plot window

5.10.4.69 Property plot_window

store plot window object

5.10.4.70 RangeTable RangeTable

italize range table, %

5.10.4.71 Property runstate

set to 0 to abort measurements

5.10.4.72 Constant Property runstate_pleasestop

state definition

5.10.4.73 Constant Property runstate_pleasestop_user

state definition

5.10.4.74 Constant Property runstate_running

state definition

5.10.4.75 Constant Property runstate_stopped

state definition

5.10.4.76 Constant Property sw_version

5.10.4.77 SweepTable SweepTable

inititalize sweep table, %

5.10.4.78 Property test

delme

5.10.4.79 Property timestamp_end

end of measurement

5.10.4.80 Property timestamp_start

start of measurement

5.10.4.81 Property use_customcontroller

whether to use the custom controller HW extension class

5.10.4.82 Property use_tempcontroller

whether to use the temperature controller HW extension class

5.10.4.83 Property winpos_FieldCtl

Window position.

5.10.4.84 Property winpos_MainCfg

Window position.

5.10.4.85 Property winpos_MPlot


mFieldSweeps = struct('from', [] Window positions: Window position

The documentation for this class was generated from the following file:


- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[Config.cpp](#)

5.11 CustomController Class Reference

Inheritance diagram for CustomController:



Collaboration diagram for CustomController:



Public Member Functions

- function [CustomController](#) ()
- function [_delete](#) (in obj)
- function [beforeAllLoops](#) (in obj, in cycles)
- function [afterAllLoops](#) (in obj, in restart_measure)
- function [beforeLoop](#) (in obj)
- function [afterLoop](#) (in obj)
- function [getValue](#) (in obj)
- function [setValue](#) (in obj, in val)

Public Attributes

- Property [filename_suffix](#)

Protected Member Functions

- function [timedFunction](#) (in obj)
- function [abortMeasurement](#) (in obj)
- function [cyclicTimer](#) (in obj, in event)
- function [setupTimer](#) (in obj, in time_interval_s)
- function [stopTimer](#) (in obj)

Protected Attributes

- Property [controller_value](#)
internal storage
- Property [controller_setpoint](#)
- Property [CyclicTimer](#)
Timer for timedFunction.

5.11.1 Constructor & Destructor Documentation

5.11.1.1 function **CustomController** ()

5.11.2 Member Function Documentation

5.11.2.1 function **_delete** (in *obj*)

5.11.2.2 function **abortMeasurement** (in *obj*) [protected],[inherited]

5.11.2.3 function **afterAllLoops** (in *obj*, in *restart_measure*)

5.11.2.4 function **afterLoop** (in *obj*)

5.11.2.5 function **beforeAllLoops** (in *obj*, in *cycles*)

5.11.2.6 function **beforeLoop** (in *obj*)

5.11.2.7 function **cyclicTimer** (in *obj*, in *event*) [protected],[inherited]

5.11.2.8 function **getValue** (in *obj*)

5.11.2.9 function **setupTimer** (in *obj*, in *time_interval_s*) [protected],[inherited]

5.11.2.10 function **setValue** (in *obj*, in *val*) [inherited]

5.11.2.11 function **stopTimer** (in *obj*) [protected],[inherited]

5.11.2.12 function **timedFunction** (in *obj*) [protected],[inherited]

5.11.3 Member Data Documentation

5.11.3.1 Property **controller_setpoint** [protected],[inherited]

5.11.3.2 Property controller_value [protected],[inherited]

internal storage

5.11.3.3 Property CyclicTimer [protected],[inherited]

Timer for timedFunction.

5.11.3.4 Property filename_suffix [inherited]


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CustomController.cpp

5.12 DataFilter Class Reference

[DataFilter](#) class (filters noise of acquired data)

Collaboration diagram for DataFilter:



Public Types

- enum [DataFilter](#) {
[median_filter](#), [mean_filter](#), [trimmed_mean_filter](#), [cryo_moke_1](#),
[cryo_moke_2](#), [histogram](#) }

Public Member Functions

- function [filter](#) (in obj, in data_array)
Filter function.

Public Attributes

- Constant Property [trimmean_percentage](#)
'mean value', ...

5.12.1 Detailed Description

[DataFilter](#) class (filters noise of acquired data)

holds differend types of filters to improve measured data

5.12.2 Member Enumeration Documentation

5.12.2.1 enum DataFilter

Enumerator

median_filter

mean_filter

trimmed_mean_filter

cryo_moke_1

cryo_moke_2

histogram

5.12.3 Member Function Documentation

5.12.3.1 function filter (in *obj*, in *data_array*)

Filter function.

Parameters

<i>data_array</i>	acquired data matrix
-------------------	----------------------

5.12.4 Member Data Documentation

5.12.4.1 Constant Property trimmean_percentage

'mean value', ...

'trimmed mean'} percentage of expected outliers


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[DataFilter.cpp](#)

5.13 DataPlot Class Reference

[DataPlot](#) class (graphical representation of data)

Collaboration diagram for DataPlot:



Public Member Functions

- function [DataPlot](#) (in `parent_axes`)
Class constructor.
- function [_delete](#) (in `obj`)
- function [updateData](#) (in `obj`)
Updates plot data.
- function [configMagnetView](#) (in `obj`, in `dataObject`)
Set up axes for viewing magnet calibration data.
- function [configPickupView](#) (in `obj`, in `dataObject`, in `pxlim`, in `pylim`)
Set up axes for viewing measurement data.
- function [updateAverage](#) (in `obj`)
Calculate average of loops for further plotting.

Public Attributes

- Constant Property [lineWidth](#)
Line width of line plot.
- Constant Property [marker](#)
Marker type for `show_points = 1`
- Constant Property [xylim_factor](#)
10%> larger x- & ylim
- Property [show_quadrature](#)
at `typePickup-View` plot quadrature signal
- Property [show_average](#)
at `typePickup-View` plot average of all loops

- Property [show_points](#)
plot datapoints instead of line
- Property [grid](#)
show grid
- Property [xlabel_nonexp](#)
use non exponential labels on X-Axis
- Property [dataObject](#)
the data object containing data for plotting (of type [MagField](#) or [MagnetizationCurve](#))
- Property [dataLoopAvg](#)
stores averaged loop of dataObject if `show_average == 1`
- Property [plotIdx](#)
if dataObject is an array of [MagnetizationCurve](#) this indicates which loop to plot index of dataObject if dataObject is an array

- Constant Property [typeNone](#)
data type of plot
- Constant Property [typeMagnet](#)
- Constant Property [typePickup](#)

- Constant Property [label_Magnet_X](#)
axes labels and legends for each data type
- Constant Property [label_Magnet_Y](#)
- Constant Property [label_Pickup_X](#)
- Constant Property [label_Pickup_Y](#)
- Constant Property [legend_pickup](#)
- Constant Property [legend_magnet](#)

- Property [pickup_xlim](#)
view area of plot in pickup view
- Property [pickup_ylim](#)

Protected Attributes

- Property [plotAxes](#)
axes object
- Property [plotType](#)
plot type eg. `typePickup`
- Property [plotObj1](#)
handle to 1st plot (layer) eg. pickup
- Property [plotObj2](#)
handle to 2nd plot (layer)
- Property [plotObj3](#)
handle to 3rd plot (layer)

5.13.1 Detailed Description

[DataPlot](#) class (graphical representation of data)

Controls an Axes object and handles the graphical representation of data. PlotType specifies the data type:

- `typeMagnet`: Calibration data where dataObject is of type [MagField](#)

- Raw data is grey-colored while polyomial is red
- `typePickup`: [Measurement](#) data of type [MagnetizationCurve](#)
 - Pickup signal shows up in grey, Quadrature in green and the average of recently measured Pickup loops is shown in red

The colors grey, green and red are defined in the axes' `ColorOrder` index 1,2 and 3. It is therefore possible to redefine them there

5.13.2 Constructor & Destructor Documentation

5.13.2.1 function `DataPlot` (in *parent_axes*)

Class constructor.

Parameters

<i>parent_axes</i>	axes object which should be used for plotting
--------------------	---

Returns

instance of [DataPlot](#) class.

5.13.3 Member Function Documentation

5.13.3.1 function `_delete` (in *obj*)

5.13.3.2 function `configMagnetView` (in *obj*, in *dataObject*)

Set up axes for viewing magnet calibration data.

Parameters

<i>dataObject</i>	the data object of type MagField
-------------------	--

5.13.3.3 function `configPickupView` (in *obj*, in *dataObject*, in *pxlim*, in *pylim*)

Set up axes for viewing measurement data.

Parameters

<i>dataObject</i>	The data object of type MagnetizationCurve
<i>pxlim</i>	Optional; Range of x-axis eg. [-100 100] or 'auto'; Set to 'auto' if omitted
<i>pylim</i>	Optional; Range of y-axis eg. [-100 100] or 'auto'; Set to lock-in sensitivity if omitted

5.13.3.4 function `updateAverage` (in *obj*)

Calculate average of loops for further plotting.

5.13.3.5 function `updateData` (in *obj*)

Updates plot data.

5.13.4 Member Data Documentation

5.13.4.1 Property dataLoopAvg

stores averaged loop of dataObject if `show_average == 1`

5.13.4.2 Property dataObject

the data object containing data for plotting (of type [MagField](#) or [MagnetizationCurve](#))

5.13.4.3 Property grid

show grid

5.13.4.4 Constant Property label_Magnet_X

axes labels and legends for each data type

5.13.4.5 Constant Property label_Magnet_Y

5.13.4.6 Constant Property label_Pickup_X

5.13.4.7 Constant Property label_Pickup_Y

5.13.4.8 Constant Property legend_magnet

5.13.4.9 Constant Property legend_pickup

5.13.4.10 Constant Property lineWidth

Line width of line plot.

5.13.4.11 Constant Property marker

Marker type for `show_points = 1`

5.13.4.12 Property pickup_xlim

view area of plot in pickup view

5.13.4.13 Property pickup_ylim

5.13.4.14 Property plotAxes [protected]

axes object

5.13.4.15 Property plotIdx

if `dataObject` is an array of [MagnetizationCurve](#) this indicates which loop to plot index of dataObject if dataObject is an array

5.13.4.16 Property plotObj1 [protected]

handle to 1st plot (layer) eg. pickup

5.13.4.17 Property plotObj2 [protected]

handle to 2nd plot (layer)

5.13.4.18 Property plotObj3 [protected]

handle to 3rd plot (layer)

5.13.4.19 Property plotType [protected]

plot type eg. typePickup

5.13.4.20 Property show_average

at typePickup-View plot average of all loops

5.13.4.21 Property show_points

plot datapoints instead of line

5.13.4.22 Property show_quadrature

at typePickup-View plot quadrature signal

5.13.4.23 Constant Property typeMagnet**5.13.4.24 Constant Property typeNone**

data type of plot

5.13.4.25 Constant Property typePickup**5.13.4.26 Property xlabel_nonexp**

use non exponential labels on X-Axis

5.13.4.27 Constant Property ylim_factor

10%> larger x- & ylim


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[DataPlot.cpp](#)

5.14 ErrorHandler Class Reference

Logs Exceptions and stores them in error_log.mat Most IO and other Exceptions are saved in error_log.mat and hence can be reviewed using static functions of the [ErrorHandler](#) class.

Collaboration diagram for ErrorHandler:



Static Public Member Functions

- static function [logError](#) (in Exception, in error_info, in comment)
Logs an Exception.
- static function [getLog](#) ()
- static function [getList](#) ()
Briefly lists all saved exceptions.
- static function [report](#) (in nr)
Get detailed information about the exception.

Public Attributes

- Constant Property [logfile](#)
- Constant Property [maxlogs](#)
maximum number of errors in log file

5.14.1 Detailed Description

Logs Exceptions and stores them in error_log.mat Most IO and other Exceptions are saved in error_log.mat and hence can be reviewed using static functions of the [ErrorHandler](#) class.

5.14.2 Member Function Documentation

5.14.2.1 static function `getList ()` [static]

Briefly lists all saved exceptions.

Return values

<i>A</i>	String containing a list of exceptions
----------	--

5.14.2.2 static function `getLog ()` [static]

5.14.2.3 static function `logError (in Exception, in error_info, in comment)` [static]

Logs an Exception.

Parameters

<i>Exception</i>	The error as an MException class
<i>error_info</i>	Short info e.g. what were we doing, where did it happen
<i>comment</i>	Additional comment

5.14.2.4 static function `report (in nr)` [static]

Get detailed information about the exception.

Parameters

<i>nr</i>	The number in the exception list
-----------	----------------------------------

Return values

<i>description</i>	of the exception
--------------------	------------------

5.14.3 Member Data Documentation

5.14.3.1 Constant Property logfile

5.14.3.2 Constant Property maxlogs

maximum number of errors in log file


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[ErrorHandler.cpp](#)

5.15 Fake_AnalogInAndOut Class Reference

[AnalogIO](#) class.

Collaboration diagram for Fake_AnalogInAndOut:



Public Member Functions

- function [Fake_AnalogInAndOut](#) ()
- function [_delete](#) (in obj)
- function [getReachTime](#) (in obj, in volt)
- function [fakeHallSig](#) (in obj, in volt)
- function [getHWinfoADC](#) (in obj)
- function [getInputRanges](#) (in obj)
- function [acquire_adc](#) (in obj, in meas_duration)
- function [fetchBuffer_adc](#) (in obj)
- function [startBackground_adc](#) (in obj)
- function [stopBackground_adc](#) (in obj)
- function [setPickupVrange](#) (in obj, in range)
- function [setHallprobeVrange](#) (in obj, in range)
- function [getHWinfoDAC](#) (in obj)
- function [getOutputRanges](#) (in obj)
- function [reachMagnet](#) (in obj, in volt_setpoint)
- function [getDACout](#) (in obj)

Public Attributes

- Constant Property [ADCCard_ID](#)

- Constant Property [ADCinputType](#)
 - Constant Property [pickup_channel](#)
 - Constant Property [pickup_q_channel](#)
 - Constant Property [hallprobe_channel](#)
 - Constant Property [use_feedback_channel](#)
 - Constant Property [ini_pickup_vrange](#)
 - Constant Property [ini_hallprobe_vrange](#)
 - Constant Property [DACCard_ID](#)
 - Constant Property [DAC_samplerate](#)
 - Constant Property [DAC_ramp](#)
- V/s.
- Constant Property [magnet_powersupply_channel](#)
 - Constant Property [magnet_powersupply_vrange](#)
 - Property [adc_initialized](#)
 - Property [dac_initialized](#)
 - Property [magnet_vout](#)
 - Property [pickup](#)
 - Property [pickup_q](#)
 - Property [hallprobe_volt](#)
 - Property [customval_adc](#)
 - Property [stddev_pickup](#)
 - Property [stddev_pickup_q](#)
 - Property [stddev_hallprobe_volt](#)
 - [DataFilter](#) filter
 - Property [initial_magnet_vout](#)
 - Property [background_acq](#)
 - Property [magnetic_field_setpoint](#)

Protected Attributes

- Property [ai](#)
 \wedge *vrange used in MPlot Levelbars; has to be const (or redefine in [Config](#))*
- Property [pickup_vrange](#)
- Property [hallprobe_vrange](#)
- Property [adc_samples_per_second](#)
- Property [adc_meas_duration](#)
- Property [adc_meas_duration_snapshot](#)
- Property [ao](#)
- Property [dac_samples_per_second](#)
current rate
- Property [fake_yoke](#)
- Property [fake_mag_curve_up](#)
- Property [fake_mag_curve_dn](#)

Private Member Functions

- function [cleanUp](#) (in obj)

Private Attributes

- Property [magnet_volt_setpoint](#)
used by reachMagnet & getDACout
- Property [magnet_ramp](#)
"
- Property [SamplesOutput](#)

5.15.1 Detailed Description

[AnalogIO](#) class.

reads and filters from ADC-hardware + writes to DAC

5.15.2 Constructor & Destructor Documentation

5.15.2.1 function Fake_AnalogInAndOut ()

5.15.3 Member Function Documentation

5.15.3.1 function _delete (in *obj*)

5.15.3.2 function acquire_adc (in *obj*, in *meas_duration*)

5.15.3.3 function cleanUp (in *obj*) [private]

5.15.3.4 function fakeHallSig (in *obj*, in *volt*)

5.15.3.5 function fetchBuffer_adc (in *obj*)

5.15.3.6 function getDACout (in *obj*)

5.15.3.7 function getHWinfoADC (in *obj*)

5.15.3.8 function getHWinfoDAC (in *obj*)

5.15.3.9 function getInputRanges (in *obj*)

5.15.3.10 function getOutputRanges (in *obj*)

5.15.3.11 function getReachTime (in *obj*, in *volt*)

5.15.3.12 function reachMagnet (in *obj*, in *volt_setpoint*)

5.15.3.13 function setHallprobeVrange (in *obj*, in *range*)

5.15.3.14 function setPickupVrange (in *obj*, in *range*)

5.15.3.15 function startBackground_adc (in *obj*)

5.15.3.16 function stopBackground_adc (in *obj*)

5.15.4 Member Data Documentation

5.15.4.1 Property adc_initialized

5.15.4.2 Property adc_meas_duration [protected]

5.15.4.3 Property adc_meas_duration_snapshot [protected]

5.15.4.4 Property adc_samples_per_second [protected]

5.15.4.5 Constant Property ADCCard_ID

5.15.4.6 Constant Property ADCInputType

5.15.4.7 Property ai [protected]

^ vrange used in MPlot Levelbars; has to be const (or redefine in [Config](#))

5.15.4.8 Property ao [protected]

5.15.4.9 Property background_acq

5.15.4.10 Property customval_adc

5.15.4.11 Property dac_initialized

5.15.4.12 Constant Property DAC_ramp

V/s.

5.15.4.13 Constant Property DAC_samplerate

5.15.4.14 Property dac_samples_per_second [protected]

current rate

5.15.4.15 Constant Property DACCard_ID

5.15.4.16 Property fake_mag_curve_dn [protected]

5.15.4.17 Property fake_mag_curve_up [protected]

5.15.4.18 Property fake_yoke [protected]

5.15.4.19 DataFilter filter

5.15.4.20 Constant Property hallprobe_channel

5.15.4.21 Property hallprobe_volt

5.15.4.22 Property hallprobe_vrange [protected]

5.15.4.23 Constant Property ini_hallprobe_vrange

5.15.4.24 Constant Property ini_pickup_vrange

5.15.4.25 Property initial_magnet_vout

5.15.4.26 Constant Property magnet_powersupply_channel

5.15.4.27 Constant Property magnet_powersupply_vrange

5.15.4.28 Property magnet_ramp [private]

"

5.15.4.29 Property magnet_volt_setpoint [private]

used by reachMagnet & getDACout

5.15.4.30 Property magnet_vout

5.15.4.31 Property magnetic_field_setpoint

5.15.4.32 Property pickup

5.15.4.33 Constant Property pickup_channel

5.15.4.34 Property pickup_q

5.15.4.35 Constant Property pickup_q_channel

5.15.4.36 Property pickup_vrange [protected]

5.15.4.37 Property SamplesOutput [private]

5.15.4.38 Property stddev_hallprobe_volt

5.15.4.39 Property stddev_pickup

5.15.4.40 Property stddev_pickup_q


5.15.4.41 Constant Property use_feedback_channel

The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[Fake_AnalogInAndOut.cpp](#)

5.16 GuiCfg Class Reference

Collaboration diagram for GuiCfg:



Public Attributes

- Constant Property [cfgmain_signal1ch](#)
Config window.
- Constant Property [cfgmain_signal2ch](#)
- Constant Property [cfgmain_show_ask_lockin](#)
whether option 'ask if lockin sensitivity is set up correctly' is shown
- Constant Property [cfgmain_show_ask_saveconfig_before_meas](#)
whether option 'ask if config shall be saved before measurement starts' is shown
- Constant Property [plotw_showquadr](#)
Plot window.
- Constant Property [plotw_lockin](#)
- Constant Property [plotw_pickup](#)
- Constant Property [plotw_quadr](#)
- Constant Property [plotw_aux](#)
- Constant Property [plot_label_Pickup_Y](#)
Plot graph text.
- Constant Property [plot_legend_pickup](#)
- Constant Property [plot_legend_magnet](#)
- Constant Property [fctl_pickup](#)
Manual field control.

- Constant Property [fctl_quadr](#)
- Constant Property [fctl_customADC](#)
- Constant Property [fctl_customADC_unit](#)
- Constant Property [fctl_fielddtxt_dndn](#)
- Constant Property [fctl_field_dndn](#)
- Constant Property [fctl_fielddtxt_dn](#)
- Constant Property [fctl_field_dn](#)
- Constant Property [fctl_fielddtxt_up](#)
- Constant Property [fctl_field_up](#)
- Constant Property [fctl_fielddtxt_upup](#)
- Constant Property [fctl_field_upup](#)
- Constant Property [file_dataformat_header](#)

Save File.

5.16.1 Member Data Documentation

5.16.1.1 Constant Property `cfgmain_show_ask_lockin`

whether option 'ask if lockin sensitivity is set up correctly' is shown

5.16.1.2 Constant Property `cfgmain_show_ask_saveconfig_before_meas`

whether option 'ask if config shall be saved before measurement starts' is shown

5.16.1.3 Constant Property `cfgmain_signal1ch`

[Config](#) window.

5.16.1.4 Constant Property `cfgmain_signal2ch`

5.16.1.5 Constant Property `fctl_customADC`

5.16.1.6 Constant Property `fctl_customADC_unit`

5.16.1.7 Constant Property `fctl_field_dn`

5.16.1.8 Constant Property `fctl_field_dndn`

5.16.1.9 Constant Property `fctl_field_up`

5.16.1.10 Constant Property `fctl_field_upup`

5.16.1.11 Constant Property `fctl_fielddtxt_dn`

5.16.1.12 Constant Property `fctl_fielddtxt_dndn`

5.16.1.13 Constant Property `fctl_fielddtxt_up`

5.16.1.14 Constant Property `fctl_fielddtxt_upup`

5.16.1.15 Constant Property `fctl_pickup`

Manual field control.

5.16.1.16 Constant Property fctl_quadr

5.16.1.17 Constant Property file_dataformat_header

Save File.

5.16.1.18 Constant Property plot_label_Pickup_Y

Plot graph text.

5.16.1.19 Constant Property plot_legend_magnet

5.16.1.20 Constant Property plot_legend_pickup

5.16.1.21 Constant Property plotw_aux

5.16.1.22 Constant Property plotw_lockin

5.16.1.23 Constant Property plotw_pickup

5.16.1.24 Constant Property plotw_quadr

5.16.1.25 Constant Property plotw_showquadr


Plot window.

The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[GuiCfg.cpp](#)

5.17 hgsetget Class Reference

Inheritance diagram for hgsetget:




The documentation for this class was generated from the following file:


- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[MagnetizationCurve.cpp](#)

5.18 HWController Class Reference

Inheritance diagram for HWController:



Collaboration diagram for HWController:



Public Member Functions

- function [HWController](#) ()
- function [_delete](#) (in obj)
- function [beforeAllLoops](#) (in obj, in cycles)
- function [afterAllLoops](#) (in obj, in restart_measure)
- function [beforeLoop](#) (in obj)
- function [afterLoop](#) (in obj)
- function [getValue](#) (in obj)
- function [setValue](#) (in obj, in val)

Public Attributes

- Property [filename_suffix](#)

Protected Member Functions

- function [timedFunction](#) (in obj)
- function [abortMeasurement](#) (in obj)
- function [cyclicTimer](#) (in obj, in event)
- function [setupTimer](#) (in obj, in time_interval_s)
- function [stopTimer](#) (in obj)

Protected Attributes

- Property [controller_value](#)
internal storage
- Property [controller_setpoint](#)
- Property [CyclicTimer](#)
Timer for timedFunction.

5.18.1 Constructor & Destructor Documentation

5.18.1.1 function [HWController](#) ()

5.18.2 Member Function Documentation

5.18.2.1 function [_delete](#) (in *obj*)

5.18.2.2 function [abortMeasurement](#) (in *obj*) [protected]

5.18.2.3 function [afterAllLoops](#) (in *obj*, in *restart_measure*)

5.18.2.4 function [afterLoop](#) (in *obj*)

5.18.2.5 function [beforeAllLoops](#) (in *obj*, in *cycles*)

5.18.2.6 function [beforeLoop](#) (in *obj*)

5.18.2.7 function [cyclicTimer](#) (in *obj*, in *event*) [protected]

5.18.2.8 function [getValue](#) (in *obj*)

5.18.2.9 function [setupTimer](#) (in *obj*, in *time_interval_s*) [protected]

5.18.2.10 function [setValue](#) (in *obj*, in *val*)

5.18.2.11 function [stopTimer](#) (in *obj*) [protected]

5.18.2.12 function [timedFunction](#) (in *obj*) [protected]

5.18.3 Member Data Documentation

5.18.3.1 Property [controller_setpoint](#) [protected]

5.18.3.2 Property [controller_value](#) [protected]

internal storage

5.18.3.3 Property CyclicTimer [protected]

Timer for timedFunction.

5.18.3.4 Property filename_suffix


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/HWController.cpp

5.19 InterlockGUIElements Class Reference

You can add GUI elements (eg. at window initialization) and enable/disable them all together enywhre else.

Collaboration diagram for InterlockGUIElements:



Public Member Functions

- function [addControlEI](#) (in obj, in handle)
Add a GUI element.
- function [addPanel](#) (in obj, in handle)
Add all elements inside a frame/panel.
- function [disableItems](#) (in obj)
Disable all GUI elements.
- function [enableItems](#) (in obj)
Enable all GUI elements.
- function [garbageCollect](#) (in obj)
Remove Objects that do not exist anymore, but are still in list.
- function [_delete](#) (in obj)

Public Attributes

- Property [Items](#)
- Property [states](#)
- Property [items_enabled](#)

5.19.1 Detailed Description

You can add GUI elements (eg. at window initialization) and enable/disable them all together anywhere else.

5.19.2 Member Function Documentation

5.19.2.1 function `_delete` (in *obj*)

5.19.2.2 function `addControlEI` (in *obj*, in *handle*)

Add a GUI element.

Parameters

<i>handle</i>	Handle of the GUI element
---------------	---------------------------

5.19.2.3 function `addPanel` (in *obj*, in *handle*)

Add all elements inside a frame/panel.

Parameters

<i>handle</i>	Handle of the GUI element
---------------	---------------------------

5.19.2.4 function `disableItems` (in *obj*)

Disable all GUI elements.

5.19.2.5 function `enableItems` (in *obj*)

Enable all GUI elements.

5.19.2.6 function `garbageCollect` (in *obj*)

Remove Objects that do not exist anymore, but are still in list.

5.19.3 Member Data Documentation

5.19.3.1 Property `Items`

5.19.3.2 Property `items_enabled`

5.19.3.3 Property `states`


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[InterlockGUIElements.cpp](#)

5.20 LevelBar Class Reference

Extends an Axes object to act as a level bar.

Collaboration diagram for LevelBar:



Public Member Functions

- function [LevelBar](#) (in parent_axes, in min_value, in max_value, in color)
Class constructor.
- function [_delete](#) (in obj)
- function [setValue](#) (in obj, in value)
set current value of level bar
- function [setTextDisplay](#) (in obj, in TextBox, in unit_str)
set textboxes for displaying the current value
- function [setMinMaxDisplay](#) (in obj, in MinTextBox, in MaxTextBox, in unit_str)
set textboxes for displaying min/max boudaries

Public Attributes

- [PhyEditBox TextDisplay](#)
handle of Text display %
- [PhyEditBox MinDisplay](#)
handle of minimum boundary display %
- [PhyEditBox MaxDisplay](#)
handle of maximum boundary display %

Protected Attributes

- Property [laxes](#)
Axes object.
- Property [lbar](#)
the bar plot
- Property [minval](#)
minimum value of boundary
- Property [maxval](#)
maximum value of boundary

5.20.1 Detailed Description

Extends an Axes object to act as a level bar.

The Axes and Text objects can be placed using Mathworks's GUIDE(TM)

usage: `LevelBar(parent_axes, min_value, max_value, color)` `parent_axes` Handle of Axes object `min_value` Minimum value of bar `max_value` Maximum Value of bar `color` Optional; Color of the [LevelBar](#); if not specified: `ColorOrder` index 1 of Axes object ist used

```
lvlBar = PkgAdvGUI.LevelBar(handles.axVOut, range(1), range(2));
```

optionally set textboxes or lables for displaying the current value parameters are the text box handle and the unit specifier (g, V, A) `'cfg.gui_bar_magnet_vout.setTextDisplay (handles.txtVOut, 'V');'` and for minimum and maximum values `'cfg.gui_bar_magnet_vout.setMinMaxDisplay (handles.txtMin,handles.txtMax, 'V');'`

The value can be changed with: `lvlBar.setValue (value) ;`

5.20.2 Constructor & Destructor Documentation

5.20.2.1 function LevelBar (in parent_axes, in min_value, in max_value, in color)

Class constructor.

Parameters

<i>parent_axes</i>	Handle of Axes object
<i>min_value</i>	Minimum value of bar
<i>max_value</i>	Maximum Value of bar
<i>color</i>	Optional; Color of the LevelBar ; if not specified: <code>ColorOrder</code> index 1 of Axes object ist used

Returns

instance of the classDocumentationExample class.

5.20.3 Member Function Documentation

5.20.3.1 function _delete (in obj)

5.20.3.2 function setMinMaxDisplay (in obj, in MinTextBox, in MaxTextBox, in unit_str)

set textboxes for displaying min/max boudaries

Parameters

<i>MinTextBox</i>	handle to minimum text box
<i>MaxTextBox</i>	handle to maximum text box
<i>unit_str</i>	Unit of value or 0 to deactivate SI representation

5.20.3.3 function setTextDisplay (in *obj*, in *TextBox*, in *unit_str*)

set textboxes for displaying the current value

Parameters

<i>TextBox</i>	handle to text box
<i>unit_str</i>	Unit of value or 0 to deactivate SI representation

5.20.3.4 function setValue (in *obj*, in *value*)

set current value of level bar

Parameters

<i>new</i>	value of level bar
------------	--------------------

5.20.4 Member Data Documentation

5.20.4.1 Property *laxes* [protected]

Axes object.

5.20.4.2 Property *lbar* [protected]

the bar plot

5.20.4.3 Property *MaxDisplay*

handle of maximum boundary display %

5.20.4.4 Property *maxval* [protected]

maximum value of boundary

5.20.4.5 Property *MinDisplay*

handle of minimum boundary display %

5.20.4.6 Property *minval* [protected]

minimum value of boundary

5.20.4.7 PhyEditBox TextDisplay

handle of Text display %


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/LevelBar.cpp

5.21 MagField Class Reference

[MagField](#) class; magnetic field calculation, calibration.

Collaboration diagram for MagField:



Public Member Functions

- function [MagField](#) (in [calibration_file_nr](#))
Class constructor.
- function [_delete](#) (in obj)
- function [save_calib](#) (in obj)
Save calibration matrix and polynom as file.
- function [FieldToVolt](#) (in obj, in field)
Calculates required power supply control voltage in order to reach a desired field value.
- function [VoltToField](#) (in obj, in volt)
Calculates expected field value from power supply control voltage.
- function [HallprobeToField](#) (in obj, in volt)
Converts hall amplifier voltage into magn.
- function [polyFitMatrix](#) (in obj)
Calculates calibration polynom from calibration matrix.

Static Public Member Functions

- static function [calibrationFiles](#) (in num_files)
Returns List of calibration files plus description plus date of calibration.
- static function [getDescription](#) (in calib_nr)
Get description of calibration file.
- static function [getDate](#) (in calib_nr)
Get calibration date of calibration file.
- static function [setDescription](#) (in calib_nr, in [description](#))
Alters description of certain calibration file.

Public Attributes

- Constant Property [num_calibration_files](#)
Maximum number of calibration Files.
- Constant Property [not_calibrated](#)
Description text for non existing calibration.
- Constant Property [poly_order](#)
order of yoke calibration polynom
- Constant Property [max_yoke_offset](#)
Oe, offsets between field setpoint and actual field.
- Property [field_matrix](#)
:field, :volt (calibration matrix)
- Property [field_poly](#)
calibration polynom
- Property [calibration_file_nr](#)
number of used calibration file
- Property [description](#)
description (name) of calibration file
- Property [yoke_offset](#)
offset between field setpoint and actual field (current remanence or calibr error)

5.21.1 Detailed Description

[MagField](#) class; magnetic field calculation, calibration.

reads hall probe for determining current magnetic field and controls power supply

5.21.2 Constructor & Destructor Documentation

5.21.2.1 function [MagField](#) (in [calibration_file_nr](#))

Class constructor.

Parameters

calibration_file_nr	Number of calibration file (1-9) to load
-------------------------------------	--

5.21.3 Member Function Documentation

5.21.3.1 `function _delete (in obj)`

5.21.3.2 `static function calibrationFiles (in num_files) [static]`

Returns List of calibration files plus description plus date of calibration.

Parameters

<i>num_files</i>	Optional; also obsolete
------------------	-------------------------

Return values

<i>str</i>	List of calibration files plus description
------------	--

5.21.3.3 function FieldToVolt (in *obj*, in *field*)

Calculates required power supply control voltage in order to reach a desired field value.

Parameters

<i>field</i>	Desired magnetic field
--------------	------------------------

Return values

<i>volt</i>	Remote control voltage for magnet power supply
-------------	--

5.21.3.4 static function getDate (in *calib_nr*) [static]

Get calibration date of calibration file.

Parameters

<i>calib_nr</i>	Calibration file number
-----------------	-------------------------

Return values

<i>str</i>	Date string
------------	-------------

5.21.3.5 static function getDescription (in *calib_nr*) [static]

Get description of calibration file.

Parameters

<i>calib_nr</i>	Calibration file number
-----------------	-------------------------

Return values

<i>str</i>	Description (Name)
------------	--------------------

5.21.3.6 function HallprobeToField (in *obj*, in *volt*)

Converts hall amplifier voltage into magn.

field value using calibration factor and offset

Parameters

<i>volt</i>	Output voltage of hall amplifier
-------------	----------------------------------

Return values

<i>field</i>	Magnetic field
--------------	----------------

5.21.3.7 function polyFitMatrix (in *obj*)

Calculates calibration polynom from calibration matrix.

5.21.3.8 function save_calib (in *obj*)

Save calibration matrix and polynom as file.

5.21.3.9 static function setDescription (in *calib_nr*, in *description*) [static]

Alters description of certain calibration file.

Parameters

<i>calib_nr</i>	Calibration file number
<i>description</i>	Description(name)

Return values

<i>str</i>	List of calibration files plus description
------------	--

5.21.3.10 function VoltToField (in *obj*, in *volt*)

Calculates expected field value from power supply control voltage.

Parameters

<i>volt</i>	Desired magnetic field
-------------	------------------------

Return values

<i>field</i>	Remote control voltage for magnet power supply
--------------	--

5.21.4 Member Data Documentation

5.21.4.1 Property calibration_file_nr

number of used calibration file

5.21.4.2 Property description

description (name) of calibration file

5.21.4.3 Property field_matrix

:field, :volt (calibration matrix)

5.21.4.4 Property field_poly

calibration polynom

5.21.4.5 Constant Property max_yoke_offset

Oe, offsets between field setpoint and actual field.

5.21.4.6 Constant Property not_calibrated

Description text for non existing calibration.

5.21.4.7 Constant Property num_calibration_files

Maximum number of calibration Files.

5.21.4.8 Constant Property poly_order

order of yoke calibration polynom

5.21.4.9 Property yoke_offset

offset between field setpoint and actual field (current remanence or calibr error)


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[MagField.cpp](#)


5.22 MagnetizationCurve Class Reference

[MagnetizationCurve](#) class (represents data of entire magnetization curve)

Inheritance diagram for MagnetizationCurve:



Collaboration diagram for MagnetizationCurve:



Public Member Functions

- function [MagnetizationCurve](#) (in data)
Constructor.
- function set [datapoints](#) (in obj, in val)
set data points.
- function get [datapoints](#) (in obj)
Returns entire data matrix.
- function [addDataPoint](#) (in obj, in val)
Add a data point; see Measurement.magnetizationLoop() near Line 143 for an example.
- function [addElementStruct](#) (in ElStruct)
obsolete. structure to array
- function [averageLoop](#) (in obj, in loops)
Fill [MagnetizationCurve](#) with an average loop.

Public Attributes

- Constant Property [dataformat_header](#)
how .data is written in to m3dat files# moved to [GuiCfg](#); VSM/MOKE edition
- Constant Property [dataformat_descriptor](#)
format string for writing to m3dat files
- Constant Property [dataformat_descriptor_read](#)
format string for reading from m3dat files
- Property [datapoints](#)
combined value: ElStruct.signal_stdev ElStruct.customval ElStruct.customval_adc]; format = :field_setpoint :field_↔ meas :signal :signalQ :temp :field_stdev :signal_stdev :customval :customval_adc(aux_adc)

Private Attributes

- Property [field_setpoint](#)
array of setpoints
- Property [field](#)
array of measured field values
- Property [signal](#)
array of signal from detector
- Property [signal_quadr](#)
array of signal 2
- Property [temperature](#)
array of temperature values (from tempctl or manually entered)
- Property [customval](#)
array of values (from customctl)
- Property [customval_adc](#)
array of auxilliary ADC input
- Property [field_stdev](#)
array of standard deviation in field measurement
- Property [signal_stdev](#)
array of standard deviation in signal measurement

5.22.1 Detailed Description

[MagnetizationCurve](#) class (represents data of entire magnetization curve)
represents data of entire magnetization curve

5.22.2 Constructor & Destructor Documentation

5.22.2.1 function [MagnetizationCurve](#) (in *data*)

Constructor.

[MagnetizationCurve](#)([*data*])

Parameters

<i>data</i>	Optional. Data from another MagnetizationCurve
-------------	--

5.22.3 Member Function Documentation

5.22.3.1 function [addDataPoint](#) (in *obj*, in *val*)

Add a data point; see [Measurement.magnetizationLoop\(\)](#) near Line 143 for an example.

Parameters

<i>val</i>	structure of entire data matrix/array
------------	---------------------------------------

5.22.3.2 function [addElementStruct](#) (in *ElStruct*)

obsolete. structure to array

5.22.3.3 function averageLoop (in *obj*, in *loops*)

Fill [MagnetizationCurve](#) with an average loop.

Parameters

<i>loops</i>	array of MagnetizationCurve to be averaged
--------------	--

5.22.3.4 function set datapoints (in *obj*, in *val*)

set data points.

To delete all: obj.datapoints = 0; set property example

5.22.3.5 function get datapoints (in *obj*)

Returns entire data matrix.

Return values

<i>val</i>	matrix of all measured values
------------	-------------------------------

5.22.4 Member Data Documentation

5.22.4.1 Property customval [private]

array of values (from customctl)

5.22.4.2 Property customval_adc [private]

array of auxilliary ADC input

5.22.4.3 Constant Property dataformat_descriptor

format string for writing to m3dat files

5.22.4.4 Constant Property dataformat_descriptor_read

format string for reading from m3dat files

5.22.4.5 Constant Property dataformat_header

how .data is written in to m3dat files# moved to [GuiCfg](#); VSM/MOKE edition

5.22.4.6 Property datapoints

combined value: ElStruct.signal_stddev ElStruct.customval ElStruct.customval_adc]; format = :field_setpoint :field↔
_meas :signal :signalQ :temp :field_stddev :signal_stddev :customval :customval_adc(aux_adc)

5.22.4.7 Property field [private]

array of measured field values

5.22.4.8 Property field_setpoint [private]

array of setpoints

5.22.4.9 Property `field_stdev` [private]

array of standard deviation in field measurement

5.22.4.10 Property `signal` [private]

array of signal from detector

5.22.4.11 Property `signal_quadr` [private]

array of signal 2

5.22.4.12 Property `signal_stdev` [private]

array of standard deviation in signal measurement

5.22.4.13 Property `temperature` [private]

array of temperature values (from tempctl or manually entered)

The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[MagnetizationCurve.cpp](#)

5.23 Measurement Class Reference

Main functions for [Measurement](#) and Calibration; Contains main parts of measurement logic.

Static Public Member Functions

- static function [updateMagnetVout](#) (in obj, in event, in AnalIO)
Timer function that updates current measurement values in GUI display.
- static function [degaussYoke](#) (in AnalIO)
Demagnetize Yoke.
- static function [viewFile](#) ()
Load a measurement data file.
- static function [ipause](#) (in seconds)
interruptible pause; interrupts if runstate != running; exploits java sleep function (seems to be more accurate)

5.23.1 Detailed Description

Main functions for [Measurement](#) and Calibration; Contains main parts of measurement logic.

5.23.2 Member Function Documentation

5.23.2.1 static function `degaussYoke` (in *AnalIO*) [static]

Demagnetize Yoke.

Produces a decaying alternating magn. field

Parameters

<i>AnaIO</i>	Running AnalogIO object
--------------	---

5.23.2.2 static function `ipause (in seconds) [static]`

interruptible pause; interrupts if runstate != running; exploits java sleep function (seems to be more accurate)

Parameters

<i>seconds</i>	time in seconds
----------------	-----------------

5.23.2.3 static function `updateMagnetVout (in obj, in event, in AnaIO) [static]`

Timer function that updates current measurement values in GUI display.

Period about 0.3 sec

Parameters

<i>AnaIO</i>	Running AnalogIO object
--------------	---

5.23.2.4 static function `viewFile () [static]`

Load a measument data file.


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[Measurement.cpp](#)

5.24 MeasurementFile Class Reference

[MeasurementFile](#) class (saves/loads measured data); See m for an example.

Collaboration diagram for MeasurementFile:



Public Member Functions

- function [MeasurementFile](#) (in [filename](#))
Create data file handling object. Usage:
- function [_delete](#) (in obj)
- function [setFilenameSuffix](#) (in obj, in suffix)
Changes `filename.m3dat` to `filename_suffix.m3dat`;
- function [saveMagnetizationCurve](#) (in obj, in mag_curve, in simple_dat)
Saves the measured data.
- function [loadMagnetizationCurve](#) (in obj)
Loads data file.
- function [readConfigFromHeader](#) (in obj, in Header, in CfgObj, in cfgproperty, in fileproperty)
Extract parameter from Header structure and overwrites parameter in [Config](#).
- function [getHeaderValue](#) (in obj, in Header, in property)
Extract parameter from Header structure and return value.

Public Attributes

- Property [filename](#)
file path to measurement file

Private Attributes

- Property [fileId_csv](#)
data separated by tabstop 'TAB'; no header
- Property [filename_csv](#)
- Property [fileId_vsmdat](#)
header + data
- Property [filename_vsmdat](#)

5.24.1 Detailed Description

[MeasurementFile](#) class (saves/loads measured data); See m for an example.

saves/loads measured data (csv files)

5.24.2 Constructor & Destructor Documentation

5.24.2.1 function [MeasurementFile](#) (in [filename](#))

Create data file handling object. Usage:

see `fileReadSample.m` for an example usage. `MF = MeasurementFile(filename)` `filename` = path to data file or `filename = 'new'` File dialog asking for new filename `filename = 'exist'` File dialog asking for existing file

Parameters

filename	Path to data file or 'new' or 'exist'
--------------------------	---------------------------------------

5.24.3 Member Function Documentation

5.24.3.1 function `_delete` (in *obj*)

5.24.3.2 function `getHeaderValue` (in *obj*, in *Header*, in *property*)

Extract parameter from Header structure and return value.

Parameters

<i>Header</i>	Header structure returned by loadMagnetizationCurve()
<i>property</i>	name of parameter in data file

Return values

<i>strvalue</i>	value of parameter as string
<i>numvalue</i>	value of parameter as integer

5.24.3.3 function loadMagnetizationCurve (in *obj*)

Loads data file.

Return values

<i>Loops</i>	array of MagnetizationCurve contains measured data
<i>Header</i>	File header as structure

5.24.3.4 function readConfigFromHeader (in *obj*, in *Header*, in *CfgObj*, in *cfgproperty*, in *fileproperty*)

Extract parameter from Header structure and overwrites parameter in [Config](#).

Parameters

<i>Header</i>	Header structure returned by loadMagnetizationCurve()
<i>CfgObj</i>	Config object
<i>cfgproperty</i>	parameter name in Config object
<i>fileproperty</i>	name of parameter in data file

5.24.3.5 function saveMagnetizationCurve (in *obj*, in *mag_curve*, in *simple_dat*)

Saves the measured data.

Parameters

<i>mag_curve</i>	MagnetizationCurve object with measured data
<i>simple_dat</i>	1: data also saved in simple tab seperated .dat file

5.24.3.6 function setFilenameSuffix (in *obj*, in *suffix*)

Changes `filename.m3dat` to `filename_suffix.m3dat`;

Parameters

<i>suffix</i>	String to add after filename like <code>'_25°C'</code>
---------------	--

5.24.4 Member Data Documentation

5.24.4.1 Property `fileld_csv` [`private`]

data separated by tabstop 'TAB', no header

5.24.4.2 Property fileld_vsmdat [private]

header + data

5.24.4.3 Property filename

file path to measurement file

5.24.4.4 Property filename_csv [private]


5.24.4.5 Property filename_vsmdat [private]

The documentation for this class was generated from the following file:


- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[MeasurementFile.cpp](#)

5.25 NetIO Class Reference

Inheritance diagram for NetIO:



Collaboration diagram for NetIO:



Public Member Functions

- function [NetIO](#) ()
- function [_delete](#) (in obj)
- virtual function [processCommands](#) (in obj, in commands)

Public Attributes

- Constant Property [updateInterval](#)
seconds

Protected Member Functions

- function [startConnection](#) (in obj)
- function [stopConnection](#) (in obj)
- function [processReadbuffer](#) (in obj)
- function [sendData](#) (in obj, in txt)
- function [waitForAnswerTo](#) (in obj, in txt)
- function [timedReceive](#) (in obj, in event)

Protected Attributes

- Property [Connection](#)
tcpclient object
- Property [ip_addr](#)
- Property [port](#)
- Property [read_buffer](#)
- Property [ReceiveTimer](#)
- Property [await_respond](#)
wait for receiving answer

5.25.1 Constructor & Destructor Documentation

5.25.1.1 function NetIO ()

5.25.2 Member Function Documentation

5.25.2.1 function _delete (in *obj*)

5.25.2.2 virtual function processCommands (in *obj*, in *commands*) [virtual]

Reimplemented in [TempController](#).

5.25.2.3 function processReadbuffer (in *obj*) [protected]

5.25.2.4 function sendData (in *obj*, in *txt*) [protected]

5.25.2.5 function startConnection (in *obj*) [protected]

5.25.2.6 function stopConnection (in *obj*) [protected]

5.25.2.7 function timedReceive (in *obj*, in *event*) [protected]

5.25.2.8 function waitForAnswerTo (in *obj*, in *txt*) [protected]

5.25.3 Member Data Documentation

5.25.3.1 Property await_respond [protected]

wait for receiving answer

5.25.3.2 Property Connection [protected]

tcpclient object

5.25.3.3 Property ip_addr [protected]

5.25.3.4 Property port [protected]

5.25.3.5 Property read_buffer [protected]

5.25.3.6 Property ReceiveTimer [protected]

5.25.3.7 Constant Property updateInterval


seconds

The documentation for this class was generated from the following file:


- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[NetIO.cpp](#)

5.26 NoController Class Reference

Inheritance diagram for NoController:



Collaboration diagram for NoController:



Public Member Functions

- function [NoController](#) ()
- function [_delete](#) (in obj)
- function [beforeAllLoops](#) (in obj, in cycles)
- function [afterAllLoops](#) (in obj, in restart_measure)
- function [getValue](#) (in obj)
- function [beforeLoop](#) (in obj)
- function [afterLoop](#) (in obj)
- function [setValue](#) (in obj, in val)

Public Attributes

- Property [filename_suffix](#)

Protected Member Functions

- function [timedFunction](#) (in obj)
- function [abortMeasurement](#) (in obj)
- function [cyclicTimer](#) (in obj, in event)
- function [setupTimer](#) (in obj, in time_interval_s)
- function [stopTimer](#) (in obj)

Protected Attributes

- Property [controller_value](#)
internal storage
- Property [controller_setpoint](#)
- Property [CyclicTimer](#)
Timer for timedFunction.

5.26.1 Constructor & Destructor Documentation

5.26.1.1 function [NoController](#) ()

5.26.2 Member Function Documentation

5.26.2.1 function [_delete](#) (in *obj*)

5.26.2.2 function [abortMeasurement](#) (in *obj*) [protected],[inherited]

5.26.2.3 function [afterAllLoops](#) (in *obj*, in *restart_measure*)

5.26.2.4 function [afterLoop](#) (in *obj*) [inherited]

5.26.2.5 function [beforeAllLoops](#) (in *obj*, in *cycles*)

5.26.2.6 function [beforeLoop](#) (in *obj*) [inherited]

5.26.2.7 function [cyclicTimer](#) (in *obj*, in *event*) [protected],[inherited]

5.26.2.8 function [getValue](#) (in *obj*)

5.26.2.9 function [setupTimer](#) (in *obj*, in *time_interval_s*) [protected],[inherited]

5.26.2.10 function [setValue](#) (in *obj*, in *val*) [inherited]

5.26.2.11 function [stopTimer](#) (in *obj*) [protected],[inherited]

5.26.2.12 function [timedFunction](#) (in *obj*) [protected],[inherited]

5.26.3 Member Data Documentation

5.26.3.1 Property [controller_setpoint](#) [protected],[inherited]

5.26.3.2 Property controller_value [protected],[inherited]

internal storage

5.26.3.3 Property CyclicTimer [protected],[inherited]

Timer for timedFunction.

5.26.3.4 Property filename_suffix [inherited]


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[NoController.cpp](#)

5.27 PhyEditBox Class Reference

Extends a text box in order to show SI values with units (kg, mA, ...)

Collaboration diagram for PhyEditBox:



Public Member Functions

- function [PhyEditBox](#) (in [hObject](#), in [unit_str](#), in [default_val](#), in [min_val](#), in [max_val](#), in [cfg_obj](#), in [cfg_varname](#))
Class constructor.
- function [_delete](#) (in obj)
- function [internal_Callback](#) (in obj, in [hObject](#), in eventdata, in handles)
Internal callback function for the text change event; Regular callback function will be called after this method.
- function [getValue](#) (in obj)
- function [setValue](#) (in obj, in val)

Static Public Member Functions

- static function [checkVal](#) (in value_text, in use_bounds, in [min_val](#), in [max_val](#), in [unit_str](#))

Public Attributes

- Property [min_val](#)
minimum value allowed
- Property [max_val](#)
maximum value allowed
- Property [default_val](#)
default value (in case of invalid input)
- Property [unit_str](#)
Unit specifier 'V', 'Oe', '', 0 (deactivate SI)

Protected Member Functions

- function [checkValue](#) (in obj)

Protected Attributes

- Property [hObject](#)
handle of text box
- Property [ori_Callback](#)
original text changed callback of text box
- Property [value](#)
current value
- Property [cfg_obj](#)
An object or structure that should get automatically updated on value changes of the [PhyEditBox](#).
- Property [cfg_varname](#)
The property of [cfg_obj](#) which will be updated.

5.27.1 Detailed Description

Extends a text box in order to show SI values with units (kg, mA, ...)

The Axes and Text objects can be placed using Mathworks's GUIDE(TM)

usage: `PhyEditBox(hObject, unit_str, default_val, min_val, max_val, [cfg_obj], [cfg_varname])` `hObject` Handle of text box object `unit_str` Unit string (V, A, Oe, ...); 0 deactivates SI representation `default_val` Default value; ignored whe `cfg_obj` is used `min_val` Minimum value allowed `max_val` Maximum value allowed `cfg_obj` Optional; An object or structure that should get automatically updated on value changes of the [PhyEditBox](#) `cfg_varname` Optional; The property of `cfg_obj` which will be updated

```
'edBox = PkgAdvGUI.PhyEditBox(handles.field, 'Oe', 0, 0, 30000, cfg, 'magnet_max_field');
```

The value can be changed with: `lvlBar.SetValue(value);`

5.27.2 Constructor & Destructor Documentation

5.27.2.1 function [PhyEditBox](#) (in *hObject*, in *unit_str*, in *default_val*, in *min_val*, in *max_val*, in *cfg_obj*, in *cfg_varname*)

Class constructor.

Parameters

<i>hObject</i>	Handle of text box object
<i>unit_str</i>	Unit string (V, A, Oe, ...); 0 deactivates SI representation
<i>default_val</i>	Default value; ignored whe <i>cfg_obj</i> is used
<i>min_val</i>	Minimum value allowed
<i>max_val</i>	Maximum value allowed
<i>cfg_obj</i>	Optional; An object or structure that should get automatically updated on value changes of the PhyEditBox
<i>cfg_varname</i>	Optional; The property of <i>cfg_obj</i> which will be updated

Returns

instance of the class

5.27.3 Member Function Documentation

5.27.3.1 `function _delete (in obj)`

5.27.3.2 `static function checkVal (in value_text, in use_bounds, in min_val, in max_val, in unit_str)` `[static]`

5.27.3.3 `function checkValue (in obj)` `[protected]`

5.27.3.4 `function getValue (in obj)`

5.27.3.5 `function internal_Callback (in obj, in hObject, in eventdata, in handles)`

Internal callback function for the text change event; Regular callback function will be called after this method.

5.27.3.6 `function setValue (in obj, in val)`

5.27.4 Member Data Documentation

5.27.4.1 `Property cfg_obj` `[protected]`

An object or structure that should get automatically updated on value changes of the [PhyEditBox](#).

5.27.4.2 `Property cfg_varname` `[protected]`

The property of *cfg_obj* which will be updated.

5.27.4.3 `Property default_val`

default value (in case of invalid input)

5.27.4.4 `Property hObject` `[protected]`

hande of text box

5.27.4.5 `Property max_val`

maximum value allowed

5.27.4.6 Property min_val

minimum value allowed

5.27.4.7 Property ori_Callback [protected]

original text changed callback of text box

5.27.4.8 Property unit_str

Unit specifier 'V', 'Oe', ", 0 (deactivate SI)

5.27.4.9 Property value [protected]

current value


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/PhyEditBox.cpp


5.28 RangeTable Class Reference

Extends a Uitable with [PhyEditBox](#) capabilities.

Inheritance diagram for RangeTable:



Collaboration diagram for RangeTable:



Public Member Functions

- function [RangeTable](#) ()
- function [createSweepTable](#) (in obj)
generate data for sweep table
- function [firstSweep](#) (in obj)
Begin sequence.
- function [nextSweep](#) (in obj)
Walk through sequence.
- function [stopSweep](#) (in obj)
Reset sequence.
- function [setTable](#) (in obj, in [table](#))
Set Uitable object which could become a sweep table.
- function [setData](#) (in obj, in [data](#))
Insert data from float array into table.
- function [getData](#) (in obj)
Set Uitable object which could become a sweep table.
- function [internal_Callback](#) (in obj, in hObject, in eventdata, in handles)
Internal callback function for the text change event; Regular callback function will be called after this method.
- function [internal_KeyPress_Callback](#) (in obj, in hObject, in eventdata, in handles)
Internal callback function for the key pressed event; Regular callback function will be called after this method.
- function [checkFields](#) (in obj, in correct)
Set Uitable object which could become a sweep table.

- function [getFieldRange](#) (in obj)
Determine highest and lowest field value of sequence.
- function [checkFrameSize](#) (in obj, in init)
align outer frame of table to the real size of the table

Public Attributes

- Constant Property [dataformat_header](#)
data block headline of m3dat files
- Constant Property [dataformat_descriptor](#)
how .data is written to m3dat files
- Constant Property [dataformat_descriptor_read](#)
how to read .data from m3dat files
- Constant Property [use_unofficial_java](#)
restore scrollbar position after editing table, etc ...;
- Property [data](#)
this might cause problems on some systems
- Property [from](#)
start of selected field sweep
- Property [to](#)
end of selected field sweep
- Property [step](#)
step size of selected field sweep
- Property [fieldarray](#)
contains the discrete sequence of field value from..to
- Property [cycles](#)
number of loops to measure
- Property [use_rangetable](#)
only one sweep (back and forth)
- Property [scrollbar_rows](#)
show scrollbar if nr of rows above this value
- Property [autoextend](#)
add additional rows when last row gets edited
- Property [scrollbar_size](#)
size of scrollbar for adjusting the table's frame
- Property [ignore_bounds](#)
ignore allowed min/max value
- Property [lastedit_row](#)
row of last text change
- Property [back_to_zero_field](#)
- Property [frame_autosize](#)
adjust frame size automatically
- Property [got_resized](#)
1: use scrollbar, do not autosize height
- Property [active_table](#)
range and unit checking after each input (use checkValues)
- Property [minval](#)
Config.instance.magnet_min_field.
- Property [maxval](#)
Config.instance.magnet_max_field.
- Property [unit_str](#)

Protected Member Functions

- function [checkValue](#) (in obj, in value_text, in out_of_bounds, in correct)
Check if value ist valid and between min/max.
- function [fillTable](#) (in obj)
fill numeric values of .data into Table
- function [getJavaHandle](#) (in obj)
access java gui handle of uitable
- function [generateFieldArray](#) (in obj)
generate discrete field values for a field sweepusing (.from .to .step) and store them in .fieldarray

Protected Attributes

- Table [table](#)
UITable object.
- Property [sweep_idx](#)
- Property [ori_Callback](#)
original callback function of table
- Property [ori_keypress_Callback](#)
original callback function of table
- Property [table_initialized](#)
needed for restoring scrollbar position
- Property [tableJavaHandle](#)

5.28.1 Detailed Description

Extends a UITable with [PhyEditBox](#) capabilities.

5.28.2 Constructor & Destructor Documentation

5.28.2.1 function RangeTable ()

5.28.3 Member Function Documentation

5.28.3.1 function checkFields (in obj, in correct) [inherited]

Set UITable object which could become a sweep table.

Parameters

<i>correct</i>	Constrain values that are out of bounds
----------------	---

Return values

<i>out_of_bounds</i>	= 1: some values are out of bounds
----------------------	------------------------------------

5.28.3.2 function checkFrameSize (in obj, in init) [inherited]

align outer frame of table to the real size of the table

Parameters

<i>init</i>	1 = window ist opened (initial start up of sweep table)
-------------	---

5.28.3.3 function `checkValue (in obj, in value_text, in out_of_bounds, in correct)` [protected],[inherited]

Check if value ist valid and between min/max.

Parameters

<i>value_text</i>	New Text as string
<i>out_of_bounds</i>	out of bounds condition of last check
<i>correct</i>	Constrain value to min/max

Return values

<i>value</i>	Numeric value
<i>value_text</i>	Formatted string of value
<i>out_of_bounds</i>	set to 1 if values exceed min/max

5.28.3.4 function `createSweepTable (in obj)`

generate data for sweep table

Return values

<i>swdata</i>	numeric array for SweepTable.setData
---------------	--

5.28.3.5 function `fillTable (in obj)` [protected],[inherited]

fill numeric values of .data into Table

5.28.3.6 function `firstSweep (in obj)` [inherited]

Begin sequence.

Set first sweep of table as current sweep. Sweep values can be accessed via .from, .to, .step. See Measurement.↵ magnetizationLoop() Line 97

Return values

<i>sequence_end</i>	= 1 if last line of table is reached
---------------------	--------------------------------------

5.28.3.7 function `generateFieldArray (in obj)` [protected],[inherited]

generate discrete field values for a field sweep using (.from .to .step) and store them in .fieldarray

5.28.3.8 function `getData (in obj)` [inherited]

Set UITable object which could become a sweep table.

Return values

<i>data</i>	array of float
-------------	----------------

5.28.3.9 function getFieldRange (in *obj*) [inherited]

Determine highest and lowest field value of sequence.

Return values

<i>field_range</i>	2 dim array whith lowest and highest field value
--------------------	--

5.28.3.10 function getJavaHandle (in *obj*) [protected],[inherited]

access java gui handle of uitable

Return values

<i>jhandle</i>	Java hande
----------------	------------

5.28.3.11 function internal_Callback (in *obj*, in *hObject*, in *eventdata*, in *handles*) [inherited]

Internal callback function for the text change event; Regular callback function will be called after this method.

5.28.3.12 function internal_KeyPress_Callback (in *obj*, in *hObject*, in *eventdata*, in *handles*) [inherited]

Internal callback function for the key pressed event; Regular callback function will be called after this method.

5.28.3.13 function nextSweep (in *obj*) [inherited]

Walk through sequence.

Set next sweep of table as current sweep. Sweep values can be accessed via .from, .to, .step. See Measurement.magnetizationLoop()

Return values

<i>sequence_end</i>	= 1 if last line of table is reached
---------------------	--------------------------------------

5.28.3.14 function setData (in *obj*, in *data*) [inherited]

Insert data from float array into table.

Parameters

<i>data</i>	array of numbers
-------------	------------------

5.28.3.15 function setTable (in *obj*, in *table*) [inherited]

Set Uitable object which could become a sweep table.

Parameters

<i>table</i>	Existing UITable
--------------	------------------

5.28.3.16 function stopSweep (in *obj*) [inherited]

Reset sequence.

5.28.4 Member Data Documentation

5.28.4.1 Property active_table [inherited]

range and unit checking after each input (use checkValues)

5.28.4.2 Property autoextend [inherited]

add additional rows when last row gets edited

5.28.4.3 Property back_to_zero_field [inherited]

5.28.4.4 Property cycles [inherited]

number of loops to measure

5.28.4.5 Property data [inherited]

this might cause problems on some systems

:from :to :step :active

5.28.4.6 Constant Property dataformat_descriptor [inherited]

how .data is written to m3dat files

5.28.4.7 Constant Property dataformat_descriptor_read [inherited]

how to read .data from m3dat files

5.28.4.8 Constant Property dataformat_header [inherited]

data block headline of m3dat files

5.28.4.9 Property fieldarray [inherited]

contains the discrete sequence of field value from..to

5.28.4.10 Property frame_autosize [inherited]

adjust frame size automatically

5.28.4.11 Property from [inherited]

start of selected field sweep

5.28.4.12 Property got_resized [inherited]

1: use scrollbar, do not autosize height

5.28.4.13 Property ignore_bounds [inherited]

ignore allowed min/max value

5.28.4.14 Property lastedit_row [inherited]

row of last text change

5.28.4.15 Property maxval [inherited]

Config.instance.magnet_max_field.

5.28.4.16 Property minval [inherited]

Config.instance.magnet_min_field.

5.28.4.17 Property ori_Callback [protected],[inherited]

original callback function of table

5.28.4.18 Property ori_keypress_Callback [protected],[inherited]

original callback function of table

5.28.4.19 Property scrollbar_rows [inherited]

show scrollbar if nr of rows above this value

5.28.4.20 Property scrollbar_size [inherited]

size of scrollbar for adjusting the table's frame

5.28.4.21 Property step [inherited]

step size of selected field sweep

5.28.4.22 Property sweep_idx [protected],[inherited]**5.28.4.23 Table table** [protected],[inherited]

UITable object.

5.28.4.24 Property `table_initialized` `[protected],[inherited]`

needed for restoring scrollbar position

5.28.4.25 Property `tableJavaHandle` `[protected],[inherited]`

5.28.4.26 Property `to` `[inherited]`

end of selected field sweep

5.28.4.27 Property `unit_str` `[inherited]`

5.28.4.28 Property `use_rangetable` `[inherited]`

only one sweep (back and forth)

5.28.4.29 Constant Property `use_unofficial_java` `[inherited]`


restore scrollbar position after editing table, etc ..;

The documentation for this class was generated from the following file:


- `C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/RangeTable.cpp`

5.29 Singleton Class Reference

Inheritance diagram for Singleton:



Collaboration diagram for Singleton:



Public Member Functions

- function [getSingletonData](#) (in *obj*)
- function [setSingletonData](#) (in *obj*, in [singletonData](#))

Static Public Member Functions

- virtual static [instance](#) ()

Private Attributes

- Property [singletonData](#)

5.29.1 Member Function Documentation

5.29.1.1 function [getSingletonData](#) (in *obj*)

5.29.1.2 virtual static [instance](#) () [static],[virtual]

5.29.1.3 function [setSingletonData](#) (in *obj*, in *singletonData*)

5.29.2 Member Data Documentation

5.29.2.1 Property [singletonData](#) [private]


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[Singleton.cpp](#)

5.30 SweepPanel Class Reference

Holds and organizes all tables options and buttons for the sweep configuration inside the fMPlot window.

Collaboration diagram for SweepPanel:



Public Member Functions

- function [SweepPanel](#) (in [chkbox_advanced](#), in `chkbox_back_to_zerofield`, in `txt_cycles`, in `txt_current_cycle`, in `tab_fielranges`, in `tab_sweepsequence`, in [pan_easy](#), in [pan_advanced](#))
Class constructor.
- function [setCurrentCycle](#) (in `obj`, in `cycle`)
set current cycle number which is then shown in `txt_current_cycle`
- function [chkbox_Callback](#) (in `obj`, in `hObject`, in `eventdata`, in `handles`)
internal callback
- function [showSweepSequence](#) (in `obj`, in `true_false`)
Switch between simple and advanced field sweep setup.
- function [chkboxChanged](#) (in `obj`)
- function [chkboxZeroField_Callback](#) (in `obj`, in `hObject`, in `eventdata`, in `handles`)
internal callback
- function [chkboxZeroFieldChanged](#) (in `obj`)

Public Attributes

- Constant Property [str_current_cycle](#)
- Property [chkbox_advanced](#)
option (checkbox) for using a table of sweeps instead of simple min and max field of a loop
- Property [chkbox_zerofield](#)
checkbox whether magn. field should be set to zero after measurement
- Property [txtCycles](#)
text box with number of cycles to measure
- Property [txtCurrentCycle](#)
text box with "Current cycle: ..."
- [SweepTable SweepTable](#)
table with explicit field sweep sequence (advanced table)
- [RangeTable RangeTable](#)
table with field sweep ranges
- Panel [pan_easy](#)
hide and show advanced sweep feature
- Panel [pan_advanced](#)
hide and show advanced sweep feature
- Property [ori_chkbox_callback](#)
original checkbox callback
- Property [ori_chkboxZeroField_Callback](#)
original checkbox callback

5.30.1 Detailed Description

Holds and organizes all tables options and buttons for the sweep configuration inside the fMPlot window.

5.30.2 Constructor & Destructor Documentation

5.30.2.1 **function SweepPanel** (in [chkbox_advanced](#), in [chkbox_back_to_zerofield](#), in [txt_cycles](#), in [txt_current_cycle](#), in [tab_fielranges](#), in [tab_sweepsequence](#), in [pan_easy](#), in [pan_advanced](#))

Class constructor.

Parameters

chkbox_↔ advanced	checkbox field range table/advanced sweep setup
chkbox_back_↔ to_zerofield	checkbox sweep back to 0 field after measurement
txt_cycles	textbox number of sweep sequence cycles
txt_current_cycle	textbox/label for indicating current cycle number
tab_fielranges	field RangeTable
tab_↔ sweepsequence	field SweepTable
pan_easy	frame with elements for simple setup
pan_advanced	frame with elements for advanced sweep setup

Returns

instance of the [SweepPanel](#) class.

5.30.3 Member Function Documentation

5.30.3.1 function `chkbox_Callback` (in *obj*, in *hObject*, in *eventdata*, in *handles*)

internal callback

5.30.3.2 function `chkboxChanged` (in *obj*)

5.30.3.3 function `chkboxZeroField_Callback` (in *obj*, in *hObject*, in *eventdata*, in *handles*)

internal callback

5.30.3.4 function `chkboxZeroFieldChanged` (in *obj*)

5.30.3.5 function `setCurrentCycle` (in *obj*, in *cycle*)

set current cycle number which is then shown in `txt_current_cycle`

5.30.3.6 function `showSweepSequence` (in *obj*, in *true_false*)

Switch between simple and advanced field sweep setup.

Parameters

<i>true_false</i>	1 = advanced setup / 0 = simple setup
-------------------	---------------------------------------

5.30.4 Member Data Documentation

5.30.4.1 Property `chkbox_advanced`

option (checkbox) for using a table of sweeps instead of simple min and max field of a loop

5.30.4.2 Property `chkbox_zerofield`

checkbox whether magn. field should be set to zero after measurement

5.30.4.3 Property `ori_chkbox_callback`

original checkbox callback

5.30.4.4 Property `ori_chkboxZeroField_Callback`

original checkbox callback

5.30.4.5 Panel `pan_advanced`

hide and show advanced sweep feature

5.30.4.6 Panel `pan_easy`

hide and show advanced sweep feature

5.30.4.7 RangeTable RangeTable

table with field sweep ranges

5.30.4.8 Constant Property str_current_cycle

5.30.4.9 SweepTable SweepTable

table with explicit field sweep sequence (advanced table)

5.30.4.10 Property txtCurrentCycle

text box with "Current cycle: ..."

5.30.4.11 Property txtCycles

text box with number of cycles to measure


The documentation for this class was generated from the following file:

- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[SweepPanel.cpp](#)


5.31 SweepTable Class Reference

Extends a Uitable with [PhyEditBox](#) capabilities and and provides the sweep sequence.

Inheritance diagram for SweepTable:



Collaboration diagram for SweepTable:



Public Member Functions

- function [SweepTable](#) (in value_min, in value_max, in unit_string)
Class constructor.
- function [firstSweep](#) (in obj)
Begin sequence.
- function [nextSweep](#) (in obj)
Walk through sequence.
- function [stopSweep](#) (in obj)
Reset sequence.
- function [setTable](#) (in obj, in [table](#))
Set UITable object which could become a sweep table.
- function [setData](#) (in obj, in [data](#))
Insert data from float array into table.
- function [getData](#) (in obj)
Set UITable object which could become a sweep table.
- function [internal_Callback](#) (in obj, in hObject, in eventdata, in handles)
Internal callback function for the text change event; Regular callback function will be called after this method.
- function [internal_KeyPress_Callback](#) (in obj, in hObject, in eventdata, in handles)
Internal callback function for the key pressed event; Regular callback function will be called after this method.
- function [checkFields](#) (in obj, in correct)
Set UITable object which could become a sweep table.
- function [getFieldRange](#) (in obj)
Determine highest and lowest field value of sequence.
- function [checkFrameSize](#) (in obj, in init)
align outer frame of table to the real size of the table

Public Attributes

- Constant Property [dataformat_header](#)
data block headline of m3dat files
- Constant Property [dataformat_descriptor](#)
how .data is written to m3dat files
- Constant Property [dataformat_descriptor_read](#)
how to read .data from m3dat files
- Constant Property [use_unofficial_java](#)
restore scrollbar position after editing table, etc ..;
- Property [data](#)
this might cause problems on some systems
- Property [from](#)
start of selected field sweep
- Property [to](#)
end of selected field sweep
- Property [step](#)
step size of selected field sweep
- Property [fieldarray](#)
contains the discrete sequence of field value from..to
- Property [cycles](#)
number of loops to measure
- Property [use_rangetable](#)
only one sweep (back and forth)
- Property [scrollbar_rows](#)
show scrollbar if nr of rows above this value
- Property [autoextend](#)
add additional rows when last row gets edited
- Property [scrollbar_size](#)
size of scrollbar for adjusting the table's frame
- Property [ignore_bounds](#)
ignore allowed min/max value
- Property [lastedit_row](#)
row of last text change
- Property [back_to_zero_field](#)
- Property [frame_autosize](#)
adjust frame size automatically
- Property [got_resized](#)
1: use scrollbar, do not autosize height
- Property [active_table](#)
range and unit checking after each input (use checkValues)
- Property [minval](#)
Config.instance.magnet_min_field.
- Property [maxval](#)
Config.instance.magnet_max_field.
- Property [unit_str](#)

Protected Member Functions

- function [checkValue](#) (in obj, in value_text, in out_of_bounds, in correct)
Check if value ist valid and between min/max.
- function [fillTable](#) (in obj)
fill numeric values of .data into Table
- function [getJavaHandle](#) (in obj)
access java gui handle of uitable
- function [generateFieldArray](#) (in obj)
generate discrete field values for a field sweepusing (.from .to .step) and store them in .fieldarray

Protected Attributes

- Table [table](#)
UITable object.
- Property [sweep_idx](#)
- Property [ori_Callback](#)
original callback function of table
- Property [ori_keypress_Callback](#)
original callback function of table
- Property [table_initialized](#)
needed for restoring scrollbar position
- Property [tableJavaHandle](#)

5.31.1 Detailed Description

Extends a UITable with [PhyEditBox](#) capabilities and and provides the sweep sequence.

5.31.2 Constructor & Destructor Documentation

5.31.2.1 function SweepTable (in value_min, in value_max, in unit_string)

Class constructor.

Parameters

<i>value_min</i>	Minimum value allowed
<i>value_max</i>	Maximum value allowed
<i>unit_str</i>	Unit string ('V', 'A', 'Oe', ...); 0 deactivates SI representation

Returns

instance of the class

5.31.3 Member Function Documentation

5.31.3.1 function checkFields (in obj, in correct)

Set UITable object which could become a sweep table.

Parameters

<i>correct</i>	Constrain values that are out of bounds
----------------	---

Return values

<i>out_of_bounds</i>	= 1: some values are out of bounds
----------------------	------------------------------------

5.31.3.2 function checkFrameSize (in *obj*, in *init*)

align outer frame of table to the real size of the table

Parameters

<i>init</i>	1 = window ist opened (initial start up of sweep table)
-------------	---

5.31.3.3 function checkValue (in *obj*, in *value_text*, in *out_of_bounds*, in *correct*) [protected]

Check if value ist valid and between min/max.

Parameters

<i>value_text</i>	New Text as string
<i>out_of_bounds</i>	out of bounds condition of last check
<i>correct</i>	Constrain value to min/max

Return values

<i>value</i>	Numeric value
<i>value_text</i>	Formatted string of value
<i>out_of_bounds</i>	set to 1 if values exceed min/max

5.31.3.4 function fillTable (in *obj*) [protected]

fill numeric values of .data into Table

5.31.3.5 function firstSweep (in *obj*)

Begin sequence.

Set first sweep of table as current sweep. Sweep values can be accessed via .from, .to, .step. See Measurement.↔ magnetizationLoop() Line 97

Return values

<i>sequence_end</i>	= 1 if last line of table is reached
---------------------	--------------------------------------

5.31.3.6 function generateFieldArray (in *obj*) [protected]

generate discrete field values for a field sweep using (.from .to .step) and store them in .fieldarray

5.31.3.7 function getData (in *obj*)

Set UITable object which could become a sweep table.

Return values

<i>data</i>	array of float
-------------	----------------

5.31.3.8 function getFieldRange (in *obj*)

Determine highest and lowest field value of sequence.

Return values

<i>field_range</i>	2 dim array whith lowest and highest field value
--------------------	--

5.31.3.9 function getJavaHandle (in *obj*) [protected]

access java gui handle of uitable

Return values

<i>jhandle</i>	Java hande
----------------	------------

5.31.3.10 function internal_Callback (in *obj*, in *hObject*, in *eventdata*, in *handles*)

Internal callback function for the text change event; Regular callback function will be called after this method.

5.31.3.11 function internal_KeyPress_Callback (in *obj*, in *hObject*, in *eventdata*, in *handles*)

Internal callback function for the key pressed event; Regular callback function will be called after this method.

5.31.3.12 function nextSweep (in *obj*)

Walk through sequence.

Set next sweep of table as current sweep. Sweep values can be accessed via .from, .to, .step. See Measurement.↵ magnetizationLoop()

Return values

<i>sequence_end</i>	= 1 if last line of table is reached
---------------------	--------------------------------------

5.31.3.13 function setData (in *obj*, in *data*)

Insert data from float array into table.

Parameters

<i>data</i>	array of numbers
-------------	------------------

5.31.3.14 function setTable (in *obj*, in *table*)

Set UITable object which sould become a sweep table.

Parameters

<i>table</i>	Existing UITable
--------------	------------------

5.31.3.15 function stopSweep (in *obj*)

Reset sequence.

5.31.4 Member Data Documentation

5.31.4.1 Property active_table

range and unit checking after each input (use checkValues)

5.31.4.2 Property autoextend

add additional rows when last row gets edited

5.31.4.3 Property back_to_zero_field

5.31.4.4 Property cycles

number of loops to measure

5.31.4.5 Property data

this might cause problems on some systems

:from :to :step :active

5.31.4.6 Constant Property dataformat_descriptor

how .data is written to m3dat files

5.31.4.7 Constant Property dataformat_descriptor_read

how to read .data from m3dat files

5.31.4.8 Constant Property dataformat_header

data block headline of m3dat files

5.31.4.9 Property fieldarray

contains the discrete sequence of field value from..to

5.31.4.10 Property frame_autosize

adjust frame size automatically

5.31.4.11 Property from

start of selected field sweep

5.31.4.12 Property got_resized

1: use scrollbar, do not autosize height

5.31.4.13 Property ignore_bounds

ignore allowed min/max value

5.31.4.14 Property lastedit_row

row of last text change

5.31.4.15 Property maxval

Config.instance.magnet_max_field.

5.31.4.16 Property minval

Config.instance.magnet_min_field.

5.31.4.17 Property ori_Callback [protected]

original callback function of table

5.31.4.18 Property ori_keypress_Callback [protected]

original callback function of table

5.31.4.19 Property scrollbar_rows

show scrollbar if nr of rows above this value

5.31.4.20 Property scrollbar_size

size of scrollbar for adjusting the table's frame

5.31.4.21 Property step

step size of selected field sweep

5.31.4.22 Property sweep_idx [protected]

5.31.4.23 Table table [protected]

UITable object.

5.31.4.24 Property table_initialized [protected]

needed for restoring scrollbar position

5.31.4.25 Property tableJavaHandle [protected]

5.31.4.26 Property to

end of selected field sweep

5.31.4.27 Property unit_str

5.31.4.28 Property use_rangetable

only one sweep (back and forth)

5.31.4.29 Constant Property use_unofficial_java


restore scrollbar position after editing table, etc ...;

The documentation for this class was generated from the following file:


- C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/[SweepTable.cpp](#)

5.32 TempController Class Reference

Inheritance diagram for TempController:



Collaboration diagram for TempController:



Public Member Functions

- function [TempController](#) ()
- function [_delete](#) (in obj)
- function [beforeAllLoops](#) (in obj, in cycles)
- function [afterAllLoops](#) (in obj, in restart_measure)
- function [beforeLoop](#) (in obj)
- function [afterLoop](#) (in obj)
- function [getValue](#) (in obj)
- function [setValue](#) (in obj, in val)

Public Attributes

- Property [meas_continue](#)
- Property [filename_suffix](#)
- Constant Property [updateInterval](#)
seconds

Protected Member Functions

- function [processCommands](#) (in obj, in commands)
- function [timedFunction](#) (in obj)
- function [abortMeasurement](#) (in obj)
- function [cyclicTimer](#) (in obj, in event)
- function [setupTimer](#) (in obj, in time_interval_s)
- function [stopTimer](#) (in obj)

- function [startConnection](#) (in obj)
- function [stopConnection](#) (in obj)
- function [processReadbuffer](#) (in obj)
- function [sendData](#) (in obj, in txt)
- function [waitForAnswerTo](#) (in obj, in txt)
- function [timedReceive](#) (in obj, in event)

Protected Attributes

- Property [controller_value](#)
internal storage
- Property [controller_setpoint](#)
- Property [CyclicTimer](#)
Timer for timedFunction.
- Property [Connection](#)
tcpclient object
- Property [ip_addr](#)
- Property [port](#)
- Property [read_buffer](#)
- Property [ReceiveTimer](#)
- Property [await_respond](#)
wait for receiving answer

5.32.1 Constructor & Destructor Documentation

5.32.1.1 function [TempController](#) ()

5.32.2 Member Function Documentation

5.32.2.1 function [_delete](#) (in *obj*)

5.32.2.2 function [abortMeasurement](#) (in *obj*) [protected],[inherited]

5.32.2.3 function [afterAllLoops](#) (in *obj*, in *restart_measure*)

5.32.2.4 function [afterLoop](#) (in *obj*) [inherited]

5.32.2.5 function [beforeAllLoops](#) (in *obj*, in *cycles*)

5.32.2.6 function [beforeLoop](#) (in *obj*) [inherited]

5.32.2.7 function [cyclicTimer](#) (in *obj*, in *event*) [protected],[inherited]

5.32.2.8 function [getValue](#) (in *obj*) [inherited]

5.32.2.9 function [processCommands](#) (in *obj*, in *commands*) [protected],[virtual]

Reimplemented from [NetIO](#).

5.32.2.10 function processReadbuffer (in *obj*) [protected],[inherited]

5.32.2.11 function sendData (in *obj*, in *txt*) [protected],[inherited]

5.32.2.12 function setupTimer (in *obj*, in *time_interval_s*) [protected],[inherited]

5.32.2.13 function setValue (in *obj*, in *val*) [inherited]

5.32.2.14 function startConnection (in *obj*) [protected],[inherited]

5.32.2.15 function stopConnection (in *obj*) [protected],[inherited]

5.32.2.16 function stopTimer (in *obj*) [protected],[inherited]

5.32.2.17 function timedFunction (in *obj*) [protected],[inherited]

5.32.2.18 function timedReceive (in *obj*, in *event*) [protected],[inherited]

5.32.2.19 function waitForAnswerTo (in *obj*, in *txt*) [protected],[inherited]

5.32.3 Member Data Documentation

5.32.3.1 Property await_respond [protected],[inherited]

wait for receiving answer

5.32.3.2 Property Connection [protected],[inherited]

tcpclient object

5.32.3.3 Property controller_setpoint [protected],[inherited]

5.32.3.4 Property controller_value [protected],[inherited]

internal storage

5.32.3.5 Property CyclicTimer [protected],[inherited]

Timer for timedFunction.

5.32.3.6 Property filename_suffix [inherited]

5.32.3.7 Property ip_addr [protected],[inherited]

5.32.3.8 Property meas_continue

5.32.3.9 Property port [protected],[inherited]

5.32.3.10 Property read_buffer [protected],[inherited]

5.32.3.11 Property ReceiveTimer [protected],[inherited]

5.32.3.12 Constant Property `updateInterval` `[inherited]`

seconds

The documentation for this class was generated from the following file:

- `C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/TempController.cpp`

Chapter 6

File Documentation

6.1 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIn_digilent.cpp File Reference

Classes

- class [AnalogIn_digilent](#)
[AnalogIO](#) class.

6.2 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIn_mcc.cpp File Reference

Classes

- class [AnalogIn_mcc](#)
Handles analog input for cards using DAQ-toolbox leagacy interface.

6.3 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogIO.cpp File Reference

Classes

- class [AnalogIO](#)
Anaolog In/Out implementation using DAQ-Toolbox legacy interface (tested with [Measurement](#) Computing and National Instruments Hardware)

6.4 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/AnalogOut_mcc.cpp File Reference

Classes

- class [AnalogOut_mcc](#)
[AnalogIO](#) class.

6.5 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/ArduinoIO.cpp File Reference

Classes

- class [ArduinoIO](#)

6.6 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/calibrateMagnet.cpp File Reference

Functions

- function [calibrateMagnet](#) (in description)
Calibrates yoke.

6.6.1 Function Documentation

6.6.1.1 function [calibrateMagnet](#) (in *description*)

Calibrates yoke.

Parameters

<i>description</i>	New name of calibration
--------------------	-------------------------

6.7 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgBoolean.cpp File Reference

Classes

- class [CfgBoolean](#)
Extends a check box in order to automatically update the corresponding [Config](#) parameter.

6.8 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgRange.cpp File Reference

Classes

- class [CfgRange](#)
Extends a text box in order to accept numeric ranges (2 dim arrays like "-10 10") and automatically update the corresponding [Config](#) parameter.

6.9 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStr.cpp File Reference

Classes

- class [CfgStr](#)

Extends a text box in order to accept string values and automatically update the corresponding [Config](#) parameter.

6.10 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CfgStrOrNum.cpp File Reference

Classes

- class [CfgStrOrNum](#)

Extends a text box in order to accept numeric an string values and automatically update the corresponding [Config](#) parameter.

6.11 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Config.cpp File Reference

Classes

- class [Config](#)

Configuration class.

6.12 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/CustomController.cpp File Reference

Classes

- class [CustomController](#)

6.13 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/daM3System.cpp File Reference

6.14 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/DataFilter.cpp File Reference

Classes

- class [DataFilter](#)

[DataFilter](#) class (filters noise of acquired data)

6.15 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/DataPlot.cpp File Reference

Classes

- class [DataPlot](#)

[DataPlot](#) class (graphical representation of data)

6.16 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/ErrorHandler.cpp File Reference

Classes

- class [ErrorHandler](#)

Logs Exceptions and stores them in error_log.mat Most IO and other Exceptions are saved in error_log.mat and hence can be reviewed using static functions of the [ErrorHandler](#) class.

6.17 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/export2karamoke.cpp File Reference

6.18 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Fake_AnalogInAndOut.cpp File Reference

Classes

- class [Fake_AnalogInAndOut](#)
[AnalogIO](#) class.

6.19 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/figCfgMain.cpp File Reference

Functions

- function [figCfgMain](#) (in varargin)
Main configuration window.
- function [disableGUI](#) (in handles)
- function [enableGUI](#) (in handles)
- function [updateTimePerDataPoint](#) (in handles)
Update text field for acquisition time per data point.
- function [onoff](#) (in bool)
Boolean to on/off.
- function [figCfgMain_OpeningFcn](#) (in hObject, in eventdata, in handles, in varargin)
Window opening function.
- function [figCfgMain_OutputFcn](#) (in hObject, in eventdata, in handles)
- function [cmdCalibrateMagnet_Callback](#) (in hObject, in eventdata, in handles)
Run magnet Calibration; Executes on button press.
- function [lstMFieldCalib_Callback](#) (in hObject, in eventdata, in handles)
- function [lstMFieldCalib_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtMFieldDescr_Callback](#) (in hObject, in eventdata, in handles)
- function [txtMFieldDescr_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtHallFactor_Callback](#) (in hObject, in eventdata, in handles)
- function [txtHallFactor_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdSaveDescr_Callback](#) (in hObject, in eventdata, in handles)
Rename magnet calibration; Executes on button press.
- function [cmdViewCalib_Callback](#) (in hObject, in eventdata, in handles)
View magnet calibration data; Executes on button press.

- function [cmdManualField_Callback](#) (in hObject, in eventdata, in handles)

Show manual field control window; Executes on button press.

- function [txtmcal_minfield_Callback](#) (in hObject, in eventdata, in handles)
- function [txtmcal_minfield_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtmcal_maxfield_Callback](#) (in hObject, in eventdata, in handles)
- function [txtmcal_maxfield_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtmcal_Vmin_Callback](#) (in hObject, in eventdata, in handles)
- function [txtmcal_Vmin_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtmcal_Vmax_Callback](#) (in hObject, in eventdata, in handles)
- function [txtmcal_Vmax_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtmcal_Vstep_Callback](#) (in hObject, in eventdata, in handles)
- function [txtmcal_Vstep_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [lstFilter_Callback](#) (in hObject, in eventdata, in handles)
- function [lstFilter_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtPauseLargeStep_Callback](#) (in hObject, in eventdata, in handles)
- function [txtPauseLargeStep_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtPauseBetween_Callback](#) (in hObject, in eventdata, in handles)
- function [txtPauseBetween_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCsamples_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCsamples_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCduration_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCduration_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCsamplescont_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCsamplescont_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtLockin_Callback](#) (in hObject, in eventdata, in handles)
- function [txtLockin_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdTestADC_Callback](#) (in hObject, in eventdata, in handles)

Run hardware test; Executes on button press.

- function [cmdSaveSettings_Callback](#) (in hObject, in eventdata, in handles)
- function [figCfgMain_CloseRequestFcn](#) (in hObject, in eventdata, in handles)

Executes when window is getting closed.

- function [cmdMeasureWindow_Callback](#) (in hObject, in eventdata, in handles)

Show measurement setup window; Executes on button press.

- function [txtADCPickupCh_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCPickupCh_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCPickupQCh_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCPickupQCh_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCHallCh_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCHallCh_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCFeedbackCh_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCFeedbackCh_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [chkADCcuseFeedback_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCPickupVRange_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCPickupVRange_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCHallVRange_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCHallVRange_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCCardID_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCCardID_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtDACCARDID_Callback](#) (in hObject, in eventdata, in handles)
- function [txtDACCARDID_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtDACPowerSupplCh_Callback](#) (in hObject, in eventdata, in handles)
- function [txtDACPowerSupplCh_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtDACsamplerate_Callback](#) (in hObject, in eventdata, in handles)
- function [txtDACsamplerate_CreateFcn](#) (in hObject, in eventdata, in handles)

- function [txtDACramp_Callback](#) (in hObject, in eventdata, in handles)
- function [txtDACramp_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [chkDynamicYokeOffset_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCDriver_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCDriver_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtDACDriver_Callback](#) (in hObject, in eventdata, in handles)
- function [txtDACDriver_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtFileExt_Callback](#) (in hObject, in eventdata, in handles)
- function [txtFileExt_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [chkLockinSensitivity_Callback](#) (in hObject, in eventdata, in handles)
- function [cmdViewDataFile_Callback](#) (in hObject, in eventdata, in handles)
- function [txtTimePerDatapoint_Callback](#) (in hObject, in eventdata, in handles)
- function [txtTimePerDatapoint_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdCfgAdvanced_Callback](#) (in hObject, in eventdata, in handles)
Show advanced configuration window; Executes on button press.
- function [cmdCfgHardware_Callback](#) (in hObject, in eventdata, in handles)
Show hardware configuration window; Executes on button press.
- function [mnuNormal_Callback](#) (in hObject, in eventdata, in handles)
Show standard configuration window; Executes when menu entry is clicked.
- function [mnuAdvanced_Callback](#) (in hObject, in eventdata, in handles)
Show advanced configuration window; Executes when menu entry is clicked.
- function [mnuHardware_Callback](#) (in hObject, in eventdata, in handles)
Show hardware(all) configuration window; Executes when menu entry is clicked.
- function [mnuConfigView_Callback](#) (in hObject, in eventdata, in handles)
- function [cmdViewFieldoffset_Callback](#) (in hObject, in eventdata, in handles)
- function [txtLargeStep_Callback](#) (in hObject, in eventdata, in handles)
- function [txtLargeStep_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [lstDACDriver_Callback](#) (in hObject, in eventdata, in handles)
- function [lstDACDriver_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [lstDACCardID_Callback](#) (in hObject, in eventdata, in handles)
- function [lstDACCardID_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [lstADCDriver_Callback](#) (in hObject, in eventdata, in handles)
- function [lstADCDriver_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [lstADCCardID_Callback](#) (in hObject, in eventdata, in handles)
- function [lstADCCardID_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCCustomValCh_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCCustomValCh_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtADCCustomValCalib_Callback](#) (in hObject, in eventdata, in handles)
- function [txtADCCustomValCalib_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [chkCustomValADC_Callback](#) (in hObject, in eventdata, in handles)
- function [chkTempCtl_Callback](#) (in hObject, in eventdata, in handles)
- function [chkCustomCtl_Callback](#) (in hObject, in eventdata, in handles)
- function [mnuLoadSave_Callback](#) (in hObject, in eventdata, in handles)
- function [mnuSave_Callback](#) (in hObject, in eventdata, in handles)
Save current configuration; Executes when menu entry is clicked.
- function [mnuSaveAs_Callback](#) (in hObject, in eventdata, in handles)
Show current configuration in different file (not the standard VSM_config.mat); Executes when menu entry is clicked.
- function [mnuLoadFrom_Callback](#) (in hObject, in eventdata, in handles)
Load configuration from different file and restart application; Executes when menu entry is clicked.
- function [chkSaveconfigMeas_Callback](#) (in hObject, in eventdata, in handles)

6.19.1 Function Documentation

6.19.1.1 function chkADCuseFeedback_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.2 function chkCustomCtl_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.3 function chkCustomValADC_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.4 function chkDynamicYokeOffset_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.5 function chkLockinSensitivity_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.6 function chkSaveconfigMeas_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.7 function chkTempCtl_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.8 function cmdCalibrateMagnet_Callback (in *hObject*, in *eventdata*, in *handles*)

Run magnet Calibration; Executes on button press.

6.19.1.9 function cmdCfgAdvanced_Callback (in *hObject*, in *eventdata*, in *handles*)

Show advanced configuration window; Executes on button press.

6.19.1.10 function cmdCfgHardware_Callback (in *hObject*, in *eventdata*, in *handles*)

Show hardware configuration window; Executes on button press.

6.19.1.11 function cmdManualField_Callback (in *hObject*, in *eventdata*, in *handles*)

Show manual field control window; Executes on button press.

6.19.1.12 function cmdMeasureWindow_Callback (in *hObject*, in *eventdata*, in *handles*)

Show measurement setup window; Executes on button press.

6.19.1.13 function cmdSaveDescr_Callback (in *hObject*, in *eventdata*, in *handles*)

Rename magnet calibration; Executes on button press.

6.19.1.14 function cmdSaveSettings_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.15 function cmdTestADC_Callback (in *hObject*, in *eventdata*, in *handles*)

Run hardware test; Executes on button press.

6.19.1.16 function cmdViewCalib_Callback (in *hObject*, in *eventdata*, in *handles*)

View magnet calibration data; Executes on button press.

6.19.1.17 function cmdViewDataFile_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.18 function cmdViewFieldoffset_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.19 function disableGUI (in *handles*)

6.19.1.20 function enableGUI (in *handles*)

6.19.1.21 function figCfgMain (in *varargin*)

Main configuration window.

Parameters

<i>varargin</i>	not used
-----------------	----------

Return values

<i>varargout</i>	not used
------------------	----------

6.19.1.22 function figCfgMain_CloseRequestFcn (in *hObject*, in *eventdata*, in *handles*)

Executes when window is getting closed.

6.19.1.23 function figCfgMain_OpeningFcn (in *hObject*, in *eventdata*, in *handles*, in *varargin*)

Window opening function.

6.19.1.24 function figCfgMain_OutputFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.25 function IstADCCardID_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.26 function IstADCCardID_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.27 function IstADCDriver_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.28 function IstADCDriver_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.29 function IstDACCardID_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.30 function IstDACCardID_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.31 function IstDACDriver_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.32 function IstDACDriver_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.33 function IstFilter_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.34 function IstFilter_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.35 function IstMFieldCalib_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.36 function IstMFieldCalib_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.37 function mnuAdvanced_Callback (in *hObject*, in *eventdata*, in *handles*)

Show advanced configuration window; Executes when menu entry is clicked.

6.19.1.38 function mnuConfigView_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.39 function mnuHardware_Callback (in *hObject*, in *eventdata*, in *handles*)

Show hardware(all) configuration window; Executes when menu entry is clicked.

6.19.1.40 function mnuLoadFrom_Callback (in *hObject*, in *eventdata*, in *handles*)

Load configuration from different file and restart application; Executes when menu entry is clicked.

6.19.1.41 function mnuLoadSave_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.42 function mnuNormal_Callback (in *hObject*, in *eventdata*, in *handles*)

Show standard configuration window; Executes when menu entry is clicked.

6.19.1.43 function mnuSave_Callback (in *hObject*, in *eventdata*, in *handles*)

Save current configuration; Executes when menu entry is clicked.

6.19.1.44 function mnuSaveAs_Callback (in *hObject*, in *eventdata*, in *handles*)

Show current configuration in different file (not the standard VSM_config.mat); Executes when menu entry is clicked.

6.19.1.45 function onoff (in *bool*)

Boolean to on/off.

Parameters

<i>bool</i>	boolean
-------------	---------

Return values

<i>val</i>	on or off
------------	-----------

6.19.1.46 function txtADCCardID_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.47 function txtADCCardID_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.48 function txtADCCustomValCalib_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.49 function txtADCCustomValCalib_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.50 function txtADCCustomValCh_Callback (in *hObject*, in *eventdata*, in *handles*)

6.19.1.51 function txtADCCustomValCh_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

- 6.19.1.52 function txtADCDriver_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.53 function txtADCDriver_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.54 function txtADCduration_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.55 function txtADCduration_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.56 function txtADCFeedbackCh_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.57 function txtADCFeedbackCh_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.58 function txtADCHallCh_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.59 function txtADCHallCh_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.60 function txtADCHallVRange_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.61 function txtADCHallVRange_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.62 function txtADCPickupCh_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.63 function txtADCPickupCh_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.64 function txtADCPickupQCh_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.65 function txtADCPickupQCh_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.66 function txtADCPickupVRange_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.67 function txtADCPickupVRange_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.68 function txtADCsamples_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.69 function txtADCsamples_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.70 function txtADCsamplescont_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.71 function txtADCsamplescont_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.72 function txtDACCardID_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.73 function txtDACCardID_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.74 function txtDACDriver_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.75 function txtDACDriver_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.76 function txtDACPowerSupplCh_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.77 function txtDACPowerSupplCh_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.78 function txtDACramp_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.79 function txtDACramp_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

- 6.19.1.80 function txtDACsamplerate_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.81 function txtDACsamplerate_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.82 function txtFileExt_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.83 function txtFileExt_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.84 function txtHallFactor_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.85 function txtHallFactor_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.86 function txtLargeStep_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.87 function txtLargeStep_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.88 function txtLockin_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.89 function txtLockin_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.90 function txtmcal_maxfield_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.91 function txtmcal_maxfield_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.92 function txtmcal_minfield_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.93 function txtmcal_minfield_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.94 function txtmcal_Vmax_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.95 function txtmcal_Vmax_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.96 function txtmcal_Vmin_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.97 function txtmcal_Vmin_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.98 function txtmcal_Vstep_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.99 function txtmcal_Vstep_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.100 function txtMFieldDescr_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.101 function txtMFieldDescr_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.102 function txtPauseBetween_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.103 function txtPauseBetween_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.104 function txtPauseLargeStep_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.105 function txtPauseLargeStep_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.106 function txtTimePerDatapoint_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.19.1.107 function txtTimePerDatapoint_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.19.1.108 `function updateTimePerDataPoint (in handles)`

Update text field for acquisition time per data point.

Parameters

<i>handles</i>	the window handles structure
----------------	------------------------------

6.20 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/figFieldControl.cpp File Reference

Functions

- function [figFieldControl](#) (in varargin)
Manual field control and hall probe calibration window.
- function [figFieldControl_OpeningFcn](#) (in hObject, in eventdata, in handles, in varargin)
Window opening function.
- function [updateFields](#) (in obj, in event, in handles)
Timer function; Updates all measured values shown in the window.
- function [disableButtons](#) (in handles)
Disable output control buttons.
- function [enableButtons](#) (in handles)
Enable output control buttons.
- function [driveMagnet](#) (in vout, in handles)
Linear approaches DAC output to setpoint (vout)
- function [driveMagnet_relative](#) (in delta_vout, in handles)
Linear approaches DAC output to setpoint (relative) (current+delta_vout)
- function [setField](#) (in field, in handles)
Reach magnetic field.
- function [setField_relative](#) (in delta_field, in handles)
Reach magnetic field relative to previous field.
- function [figFieldControl_OutputFcn](#) (in hObject, in eventdata, in handles)
- function [txtHallFct_Callback](#) (in hObject, in eventdata, in handles)
- function [txtHallFct_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtField_Callback](#) (in hObject, in eventdata, in handles)
- function [txtField_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtVOut_Callback](#) (in hObject, in eventdata, in handles)
- function [txtVOut_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdZero_Callback](#) (in hObject, in eventdata, in handles)
Set DAC voltage to zero; Executes on button press.
- function [cmdUp_Callback](#) (in hObject, in eventdata, in handles)
Increase DAC voltage (small step); Executes on button press.
- function [cmdUpUp_Callback](#) (in hObject, in eventdata, in handles)
Increase DAC voltage (large step); Executes on button press.
- function [cmdDn_Callback](#) (in hObject, in eventdata, in handles)
Decrease DAC voltage (small step); Executes on button press.
- function [cmdDnDn_Callback](#) (in hObject, in eventdata, in handles)
Decrease DAC voltage (large step); Executes on button press.
- function [txtHallVolt_Callback](#) (in hObject, in eventdata, in handles)
- function [txtHallVolt_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdAdjustHall_Callback](#) (in hObject, in eventdata, in handles)
Runs hall probe calibration (linear calibration factor); Executes on button press.
- function [cmdAdjustVOut_Callback](#) (in hObject, in eventdata, in handles)
- function [cmdSetVOut_Callback](#) (in hObject, in eventdata, in handles)

Set DAC output voltage, which is connected to power supply remote; Executes on button press.

- function [figFieldCtl_CloseRequestFcn](#) (in hObject, in eventdata, in handles)

Executes when user attempts to close window; Executes on button press.

- function [txtVOutNow_Callback](#) (in hObject, in eventdata, in handles)
- function [txtVOutNow_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdInterrupt_Callback](#) (in hObject, in eventdata, in handles)

Immediately stops ramping up powersupply (kind of emergency stop); Executes on button press.

- function [txtPickup_Callback](#) (in hObject, in eventdata, in handles)
- function [txtPickup_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtPickupQ_Callback](#) (in hObject, in eventdata, in handles)
- function [txtPickupQ_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdFZero_Callback](#) (in hObject, in eventdata, in handles)

Set field to zero; Executes on button press.

- function [cmdFUp_Callback](#) (in hObject, in eventdata, in handles)

Raises magnetic field (small step); Executes on button press.

- function [cmdFUpUp_Callback](#) (in hObject, in eventdata, in handles)

Raises magnetic field (large step); Executes on button press.

- function [cmdFDn_Callback](#) (in hObject, in eventdata, in handles)

Decreases magnetic field (small step); Executes on button press.

- function [cmdFDnDn_Callback](#) (in hObject, in eventdata, in handles)

Decreases magnetic field (large step); Executes on button press.

- function [txtFieldSetpoint_Callback](#) (in hObject, in eventdata, in handles)
- function [txtFieldSetpoint_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdSetField_Callback](#) (in hObject, in eventdata, in handles)

Set field setpoint; Executes on button press.

- function [txtHallOff_Callback](#) (in hObject, in eventdata, in handles)
- function [txtHallOff_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdAdjustHallOffset_Callback](#) (in hObject, in eventdata, in handles)

Determine offset voltage of hall probe amplifier; Executes on button press.

- function [chkWriteProtect_Callback](#) (in hObject, in eventdata, in handles)
- function [cmdDegauss_Callback](#) (in hObject, in eventdata, in handles)

Demagnetize yoke (degauss); Executes on button press.

- function [txtCustomADC_Callback](#) (in hObject, in eventdata, in handles)
- function [txtCustomADC_CreateFcn](#) (in hObject, in eventdata, in handles)

6.20.1 Function Documentation

6.20.1.1 function [chkWriteProtect_Callback](#) (in hObject, in eventdata, in handles)

6.20.1.2 function [cmdAdjustHall_Callback](#) (in hObject, in eventdata, in handles)

Runs hall probe calibration (linear calibration factor); Executes on button press.

6.20.1.3 function [cmdAdjustHallOffset_Callback](#) (in hObject, in eventdata, in handles)

Determine offset voltage of hall probe amplifier; Executes on button press.

6.20.1.4 function [cmdAdjustVOut_Callback](#) (in hObject, in eventdata, in handles)

6.20.1.5 function [cmdDegauss_Callback](#) (in hObject, in eventdata, in handles)

Demagnetize yoke (degauss); Executes on button press.

6.20.1.6 function cmdDn_Callback (in *hObject*, in *eventdata*, in *handles*)

Decrease DAC voltage (small step); Executes on button press.

6.20.1.7 function cmdDnDn_Callback (in *hObject*, in *eventdata*, in *handles*)

Decrease DAC voltage (large step); Executes on button press.

6.20.1.8 function cmdFDn_Callback (in *hObject*, in *eventdata*, in *handles*)

Decreases magnetic field (small step); Executes on button press.

6.20.1.9 function cmdFDnDn_Callback (in *hObject*, in *eventdata*, in *handles*)

Decreases magnetic field (large step); Executes on button press.

6.20.1.10 function cmdFUp_Callback (in *hObject*, in *eventdata*, in *handles*)

Raises magnetic field (small step); Executes on button press.

6.20.1.11 function cmdFUpUp_Callback (in *hObject*, in *eventdata*, in *handles*)

Raises magnetic field (large step); Executes on button press.

6.20.1.12 function cmdFZero_Callback (in *hObject*, in *eventdata*, in *handles*)

Set field to zero; Executes on button press.

6.20.1.13 function cmdInterrupt_Callback (in *hObject*, in *eventdata*, in *handles*)

Immediately stops ramping up powersupply (kind of emergency stop); Executes on button press.

6.20.1.14 function cmdSetField_Callback (in *hObject*, in *eventdata*, in *handles*)

Set field setpoint; Executes on button press.

6.20.1.15 function cmdSetVOut_Callback (in *hObject*, in *eventdata*, in *handles*)

Set DAC output voltage, which is connected to power supply remote; Executes on button press.

6.20.1.16 function cmdUp_Callback (in *hObject*, in *eventdata*, in *handles*)

Increase DAC voltage (small step); Executes on button press.

6.20.1.17 function cmdUpUp_Callback (in *hObject*, in *eventdata*, in *handles*)

Increase DAC voltage (large step); Executes on button press.

6.20.1.18 function cmdZero_Callback (in *hObject*, in *eventdata*, in *handles*)

Set DAC voltage to zero; Executes on button press.

6.20.1.19 function disableButtons (in *handles*)

Disable output control buttons.

Parameters

<i>handles</i>	the window handles structure
----------------	------------------------------

6.20.1.20 function driveMagnet (in *vout*, in *handles*)

Linear approaches DAC output to setpoint (vout)

Parameters

<i>vout</i>	DAC output voltage setpoint
<i>handles</i>	the window handles structure

6.20.1.21 function driveMagnet_relative (in *delta_vout*, in *handles*)

Linear approaches DAC output to setpoint (relative) (current+delta_vout)

Parameters

<i>delta_vout</i>	change in DAC output voltage
<i>handles</i>	the window handles structure

6.20.1.22 function enableButtons (in *handles*)

Enable output control buttons.

Parameters

<i>handles</i>	the window handles structure
----------------	------------------------------

6.20.1.23 function figFieldControl (in *varargin*)

Manual field control and hall probe calibration window.

Parameters

<i>varargin</i>	not used
-----------------	----------

Return values

<i>varargout</i>	not used
------------------	----------

6.20.1.24 function figFieldControl_OpeningFcn (in *hObject*, in *eventdata*, in *handles*, in *varargin*)

Window opening function.

6.20.1.25 function figFieldControl_OutputFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.26 function figFieldCtl_CloseRequestFcn (in *hObject*, in *eventdata*, in *handles*)

Executes when user attempts to close window; Executes on button press.

6.20.1.27 function setField (in *field*, in *handles*)

Reach magnetic field.

Parameters

<i>field</i>	magnetic field setpoint
<i>handles</i>	the window handles structure

6.20.1.28 function setField_relative (in *delta_field*, in *handles*)

Reach magnetic field relative to previous field.

Parameters

<i>delta_field</i>	change in magnetic field setpoint
<i>handles</i>	the window handles structure

6.20.1.29 function txtCustomADC_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.30 function txtCustomADC_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.31 function txtField_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.32 function txtField_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.33 function txtFieldSetpoint_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.34 function txtFieldSetpoint_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.35 function txtHallFct_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.36 function txtHallFct_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.37 function txtHallOff_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.38 function txtHallOff_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.39 function txtHallVolt_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.40 function txtHallVolt_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.41 function txtPickup_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.42 function txtPickup_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.20.1.43 function txtPickupQ_Callback (in *hObject*, in *eventdata*, in *handles*)

6.20.1.44 function txtPickupQ_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

- 6.20.1.45 function txtVOut_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.20.1.46 function txtVOut_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.20.1.47 function txtVOutNow_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.20.1.48 function txtVOutNow_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.20.1.49 function updateFields (in *obj*, in *event*, in *handles*)

Timer function; Updates all measured values shown in the window.

6.21 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/fileReadSample.cpp File Reference

6.22 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/findjobj.cpp File Reference

Functions

- function [findjobj](#) (in container, in varargin)
- function [warnInvisible](#) (in varargin)
- function [paramSupplied](#) (in paramsList, in paramName)
- function [getCurrentFigure](#) ()
- function [getRootPanel](#) (in hFig)
- function [traverseContainer](#) (in jcontainer, in level, in parent)
- function [getXY](#) (in jcontainer)
- function [getNumMenuComponents](#) (in jcontainer)
- function [removeDuplicateNode](#) (in thisIdx)
- function [processArgs](#) (in varargin)
- function [processPrintArgs](#) (in varargin)
- function [processPositionArgs](#) (in varargin)
- function [processSizeArgs](#) (in varargin)
- function [processClassArgs](#) (in varargin)
- function [processPropertyArgs](#) (in varargin)
- function [processDepthArgs](#) (in varargin)
- function [charizeData](#) (in data)
- function [setProp](#) (in list, in name, in value, in category)
- function [getTreeData](#) (in data)
- function [getCbsData](#) (in obj, in stripStdCbsFlag)
- function [getRalativeDivlocation](#) (in jDiv)
- function [setTreeNodeIcon](#) (in treenode, in container)
- function [presentObjectTree](#) ()
- function [resizeImg](#) (in varargin)
- function [disableDbstopError](#) ()
- function [restoreDbstopError](#) (in identifiers)
- function [expandNode](#) (in progressBar, in tree, in tree_h, in parentNode, in parentRow)
- function [createJButton](#) (in nameStr, in handler, in toolTipText)
- function [flashComponent](#) (in jComps, in delaySecs, in numTimes)
- function [nodeSelected](#) (in src, in evd, in tree)
- function [iff](#) (in test, in trueVal, in falseVal)
- function [getPropsHtml](#) (in nodeHandle, in dataFields)

- function [updateNodeTooltip](#) (in nodeHandle, in uiObject)
- function [nodeExpanded](#) (in src, in evd, in tree)
- function [setIconSize](#) (in iconImage)
- function [getChildrenNodes](#) (in tree, in parentNode, in isRootHGNode)
- function [getNodeName](#) (in hndl, in charsLimit)
- function [stripStdCbs](#) (in evNames)
- function [cbHideStdCbs_Callback](#) (in src, in evd, in callbacksTable, in varargin)
- function [btWebsite_Callback](#) (in src, in evd, in varargin)
- function [btRefresh_Callback](#) (in src, in evd, in varargin)
- function [btExport_Callback](#) (in src, in evd, in varargin)
- function [btFocus_Callback](#) (in src, in evd, in varargin)
- function [btInspect_Callback](#) (in src, in evd, in varargin)
- function [btCheckFex_Callback](#) (in src, in evd, in varargin)
- function [checkVersion](#) ()
- function [getTopSelectedObject](#) (in jTree, in root)
- function [tbCallbacksChanged](#) (in src, in evd, in object, in table)
- function [revertCbTableModification](#) (in table, in modifiedRowIdx, in modifiedColIdx, in cbName, in object, in errMsg)
- function [getLabelsJavaPos](#) (in container)
- function [traverseHGContainer](#) (in hcontainer, in level, in parent)
- function [dispError](#) ()
- function [ischar](#) (in data)
- function [setTreeContextMenu](#) (in obj, in node, in tree_h)
- function [treeMousePressedCallback](#) (in hTree, in eventData, in tree_h)
- function [menuRemoveItem](#) (in hObj, in eventData, in jmenu, in item)
- function [getNodeTitleStr](#) (in obj, in node)
- function [treeMouseMovedCallback](#) (in hTree, in eventData)
- function [requestFocus](#) (in hTree, in eventData, in obj)

6.22.1 Function Documentation

6.22.1.1 function [btCheckFex_Callback](#) (in *src*, in *evd*, in *varargin*)

6.22.1.2 function [btExport_Callback](#) (in *src*, in *evd*, in *varargin*)

6.22.1.3 function [btFocus_Callback](#) (in *src*, in *evd*, in *varargin*)

6.22.1.4 function [btInspect_Callback](#) (in *src*, in *evd*, in *varargin*)

6.22.1.5 function [btRefresh_Callback](#) (in *src*, in *evd*, in *varargin*)

6.22.1.6 function [btWebsite_Callback](#) (in *src*, in *evd*, in *varargin*)

6.22.1.7 function [cbHideStdCbs_Callback](#) (in *src*, in *evd*, in *callbacksTable*, in *varargin*)

6.22.1.8 function [charizeData](#) (in *data*)

6.22.1.9 function [checkVersion](#) ()

6.22.1.10 function [createJButton](#) (in *nameStr*, in *handler*, in *toolTipText*)

6.22.1.11 function [disableDbstopError](#) ()

6.22.1.12 function [dispError](#) ()

- 6.22.1.13 function `expandNode` (in *progressBar*, in *tree*, in *tree_h*, in *parentNode*, in *parentRow*)
- 6.22.1.14 function `findjobj` (in *container*, in *varargin*)
- 6.22.1.15 function `flashComponent` (in *jComps*, in *delaySecs*, in *numTimes*)
- 6.22.1.16 function `getCbsData` (in *obj*, in *stripStdCbsFlag*)
- 6.22.1.17 function `getChildrenNodes` (in *tree*, in *parentNode*, in *isRootHGNode*)
- 6.22.1.18 function `getCurrentFigure` ()
- 6.22.1.19 function `getLabelsJavaPos` (in *container*)
- 6.22.1.20 function `getNodeName` (in *hndl*, in *charsLimit*)
- 6.22.1.21 function `getNodeTitleStr` (in *obj*, in *node*)
- 6.22.1.22 function `getNumMenuComponents` (in *jcontainer*)
- 6.22.1.23 function `getPropsHtml` (in *nodeHandle*, in *dataFields*)
- 6.22.1.24 function `getRalativeDivlocation` (in *jDiv*)
- 6.22.1.25 function `getRootPanel` (in *hFig*)
- 6.22.1.26 function `getTopSelectedObject` (in *jTree*, in *root*)
- 6.22.1.27 function `getTreeData` (in *data*)
- 6.22.1.28 function `getXY` (in *jcontainer*)
- 6.22.1.29 function `iff` (in *test*, in *trueVal*, in *falseVal*)
- 6.22.1.30 function `ischar` (in *data*)
- 6.22.1.31 function `menuRemoveItem` (in *hObj*, in *eventData*, in *jmenu*, in *item*)
- 6.22.1.32 function `nodeExpanded` (in *src*, in *evd*, in *tree*)
- 6.22.1.33 function `nodeSelected` (in *src*, in *evd*, in *tree*)
- 6.22.1.34 function `paramSupplied` (in *paramsList*, in *paramName*)
- 6.22.1.35 function `presentObjectTree` ()
- 6.22.1.36 function `processArgs` (in *varargin*)
- 6.22.1.37 function `processClassArgs` (in *varargin*)
- 6.22.1.38 function `processDepthArgs` (in *varargin*)
- 6.22.1.39 function `processPositionArgs` (in *varargin*)
- 6.22.1.40 function `processPrintArgs` (in *varargin*)

- 6.22.1.41 function processPropertyParams (in varargin)
- 6.22.1.42 function processSizeArgs (in varargin)
- 6.22.1.43 function removeDuplicateNode (in thisIdx)
- 6.22.1.44 function requestFocus (in hTree, in eventData, in obj)
- 6.22.1.45 function resizeImg (in varargin)
- 6.22.1.46 function restoreDbstopError (in identifiers)
- 6.22.1.47 function revertCbTableModification (in table, in modifiedRowIdx, in modifiedColIdx, in cbName, in object, in errMsg)
- 6.22.1.48 function setIconSize (in iconName)
- 6.22.1.49 function setProp (in list, in name, in value, in category)
- 6.22.1.50 function setTreeContextMenu (in obj, in node, in tree_h)
- 6.22.1.51 function setTreeNodeIcon (in treeNode, in container)
- 6.22.1.52 function stripStdCbs (in evNames)
- 6.22.1.53 function tbCallbacksChanged (in src, in evd, in object, in table)
- 6.22.1.54 function traverseContainer (in jcontainer, in level, in parent)
- 6.22.1.55 function traverseHGContainer (in hcontainer, in level, in parent)
- 6.22.1.56 function treeMouseMovedCallback (in hTree, in eventData)
- 6.22.1.57 function treeMousePressedCallback (in hTree, in eventData, in tree_h)
- 6.22.1.58 function updateNodeTooltip (in nodeHandle, in uiObject)
- 6.22.1.59 function warnInvisible (in varargin)

6.23 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/fMPlot.cpp File Reference

Functions

- function [fMPlot](#) (in varargin)
Measurement, Calibration and Plot window usage:
- function [fMPlot_OpeningFcn](#) (in hObject, in eventdata, in handles, in varargin)
Window opening function.
- function [resizeOriginal](#) (in hObject, in wsize)
- function [resizeGUI](#) (in winObject, in handles, in oldsize, in newsize)
- function [setUnitPixels](#) (in handles)
- function [setUnitNormalized](#) (in handles)
- function [fMPlot_OutputFcn](#) (in hObject, in eventdata, in handles)
- function [figMPlot_CloseRequestFcn](#) (in hObject, in eventdata, in handles)
Executes when user attempts to close figMPlot.

- function [cmdStopMeasurement_Callback](#) (in hObject, in eventdata, in handles)
Immediately abort measurement; Executes on button press.
- function [chkAdvanced_Callback](#) (in hObject, in eventdata, in handles)
- function [chkZeroField_Callback](#) (in hObject, in eventdata, in handles)
- function [txtCycles_Callback](#) (in hObject, in eventdata, in handles)
- function [txtCycles_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtEasyTo_Callback](#) (in hObject, in eventdata, in handles)
- function [txtEasyTo_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtEasyStep_Callback](#) (in hObject, in eventdata, in handles)
- function [txtEasyStep_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [grpAdvancedSweeps_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtVOut_Callback](#) (in hObject, in eventdata, in handles)
- function [txtVOut_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtFieldOut_Callback](#) (in hObject, in eventdata, in handles)
- function [txtFieldOut_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtPickup_Callback](#) (in hObject, in eventdata, in handles)
- function [txtPickup_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtPickupQ_Callback](#) (in hObject, in eventdata, in handles)
- function [txtPickupQ_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [cmdMeasure_Callback](#) (in hObject, in eventdata, in handles)
Start measuement (advanced sweep table); Executes on button press.
- function [cmdCreateSequenceAndMeasure_Callback](#) (in hObject, in eventdata, in handles)
Start simplified measurement; Create sweep table from field range table and start; Executes on button press.
- function [figMPlot_SizeChangedFcn](#) (in hObject, in eventdata, in handles)
- function [figMPlot_WindowButtonDownFcn](#) (in hObject, in eventdata, in handles)
- function [figMPlot_ButtonDownFcn](#) (in hObject, in eventdata, in handles)
- function [txtLockin_Callback](#) (in hObject, in eventdata, in handles)
- function [txtLockin_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtFieldSetpoint_Callback](#) (in hObject, in eventdata, in handles)
- function [txtFieldSetpoint_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [chkShowPoints_Callback](#) (in hObject, in eventdata, in handles)
- function [txtComment_Callback](#) (in hObject, in eventdata, in handles)
Update measurement comment; Executes on button press.
- function [txtComment_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [lstLoopPlots_Callback](#) (in hObject, in eventdata, in handles)
- function [lstLoopPlots_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [chkShowAverage_Callback](#) (in hObject, in eventdata, in handles)
- function [chkShowQuadrature_Callback](#) (in hObject, in eventdata, in handles)
- function [togglebutton1_Callback](#) (in hObject, in eventdata, in handles)
- function [chkLoopPause_Callback](#) (in hObject, in eventdata, in handles)
- function [figMPlot_KeyPressFcn](#) (in hObject, in eventdata, in handles)
- function [txtTempManual_Callback](#) (in hObject, in eventdata, in handles)
- function [txtTempManual_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [txtCustomADC_Callback](#) (in hObject, in eventdata, in handles)
- function [txtCustomADC_CreateFcn](#) (in hObject, in eventdata, in handles)
- function [chkTemperatureCtl_Callback](#) (in hObject, in eventdata, in handles)
- function [chkCustomCtl_Callback](#) (in hObject, in eventdata, in handles)

6.23.1 Function Documentation

6.23.1.1 function `chkAdvanced_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.2 function `chkCustomCtl_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.3 function `chkLoopPause_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.4 function `chkShowAverage_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.5 function `chkShowPoints_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.6 function `chkShowQuadrature_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.7 function `chkTemperatureCtl_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.8 function `chkZeroField_Callback` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.9 function `cmdCreateSequenceAndMeasure_Callback` (in *hObject*, in *eventdata*, in *handles*)

Start simplified measurement; Create sweep table from field range table and start; Executes on button press.

6.23.1.10 function `cmdMeasure_Callback` (in *hObject*, in *eventdata*, in *handles*)

Start measuement (advanced sweep table); Executes on button press.

6.23.1.11 function `cmdStopMeasurement_Callback` (in *hObject*, in *eventdata*, in *handles*)

Immediately abort measurement; Executes on button press.

6.23.1.12 function `figMPlot_ButtonDownFcn` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.13 function `figMPlot_CloseRequestFcn` (in *hObject*, in *eventdata*, in *handles*)

Executes when user attempts to close figMPlot.

6.23.1.14 function `figMPlot_KeyPressFcn` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.15 function `figMPlot_SizeChangedFcn` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.16 function `figMPlot_WindowButtonDownFcn` (in *hObject*, in *eventdata*, in *handles*)

6.23.1.17 function `fMPlot` (in *varargin*)

[Measurement](#), Calibration and Plot window usage:

Usage:

```
plotwindow = fMPlot(window_appearance, additional_windowtext);

@b additional_windowtext: additional text displayed in window title
@b window_appearance:
    window_appearance = Config.mplotview_magnetcalibration : window
    layout for magnet calibration
    window_appearance = Config.mplotview_measure : window layout for a
    normal measurement
    window_appearance = Config.mplotview_viewcalibration : window layout
    for viewing a voltage-field plot of a magnet calibration file
```

```

window_appearance = Config.mplotview_viewdatafile      : window layout for
viewing previously measured magnetization data

```

Parameters

<i>varargin</i>	Integer defined by Config.mplotview_*; specifying use case related window appearance
-----------------	--

Return values

<i>varargout</i>	not used
------------------	----------

6.23.1.18 function fMPlot_OpeningFcn (in *hObject*, in *eventdata*, in *handles*, in *varargin*)

Window opening function.

6.23.1.19 function fMPlot_OutputFcn (in *hObject*, in *eventdata*, in *handles*)

6.23.1.20 function grpAdvancedSweeps_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.23.1.21 function IstLoopPlots_Callback (in *hObject*, in *eventdata*, in *handles*)

6.23.1.22 function IstLoopPlots_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.23.1.23 function resizeGUI (in *winObject*, in *handles*, in *oldsize*, in *newsize*)

6.23.1.24 function resizeOriginal (in *hObject*, in *wsizes*)

6.23.1.25 function setUnitNormalized (in *handles*)

6.23.1.26 function setUnitPixels (in *handles*)

6.23.1.27 function togglebutton1_Callback (in *hObject*, in *eventdata*, in *handles*)

6.23.1.28 function txtComment_Callback (in *hObject*, in *eventdata*, in *handles*)

Update measurement comment; Executes on button press.

6.23.1.29 function txtComment_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.23.1.30 function txtCustomADC_Callback (in *hObject*, in *eventdata*, in *handles*)

6.23.1.31 function txtCustomADC_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.23.1.32 function txtCycles_Callback (in *hObject*, in *eventdata*, in *handles*)

6.23.1.33 function txtCycles_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.23.1.34 function txtEasyStep_Callback (in *hObject*, in *eventdata*, in *handles*)

6.23.1.35 function txtEasyStep_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.23.1.36 function txtEasyTo_Callback (in *hObject*, in *eventdata*, in *handles*)

6.23.1.37 function txtEasyTo_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

- 6.23.1.38 function txtFieldOut_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.39 function txtFieldOut_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.40 function txtFieldSetpoint_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.41 function txtFieldSetpoint_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.42 function txtLockin_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.43 function txtLockin_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.44 function txtPickup_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.45 function txtPickup_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.46 function txtPickupQ_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.47 function txtPickupQ_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.48 function txtTempManual_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.49 function txtTempManual_CreateFcn (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.50 function txtVOut_Callback (in *hObject*, in *eventdata*, in *handles*)
- 6.23.1.51 function txtVOut_CreateFcn (in *hObject*, in *eventdata*, in *handles*)

6.24 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/GuiCfg.cpp File Reference

Classes

- class [GuiCfg](#)

6.25 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/HWController.cpp File Reference

Classes

- class [HWController](#)

6.26 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/InterlockGUIElements.cpp File Reference

Classes

- class [InterlockGUIElements](#)

You can add GUI elements (eg. at window initialization) and enable/disable them all together enywhre else.

6.27 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/LevelBar.cpp File Reference

Classes

- class [LevelBar](#)
Extends an Axes object to act as a level bar.

6.28 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/M3System.cpp File Reference

6.29 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MagField.cpp File Reference

Classes

- class [MagField](#)
[MagField](#) class; magnetic field calculation, calibration.

6.30 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Magnetization↔Curve.cpp File Reference

Classes

- class [MagnetizationCurve](#)
[MagnetizationCurve](#) class (represents data of entire magnetization curve)

6.31 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/magnetization↔Loop.cpp File Reference

Functions

- function [magnetizationLoop](#) ()
Main measurement function; proceeds a complete measurement.

6.31.1 Function Documentation

6.31.1.1 function magnetizationLoop ()

Main measurement function; proceeds a complete measurement.

6.32 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Measurement.cpp File Reference

Classes

- class [Measurement](#)

Main functions for [Measurement](#) and Calibration; Contains main parts of measurement logic.

6.33 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/MeasurementFile.cpp File Reference

Classes

- class [MeasurementFile](#)

*[MeasurementFile](#) class (saves/loads measured data); See *m* for an example.*

6.34 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/NetIO.cpp File Reference

Classes

- class [NetIO](#)

6.35 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/NoController.cpp File Reference

Classes

- class [NoController](#)

6.36 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/PhyEditBox.cpp File Reference

Classes

- class [PhyEditBox](#)

Extends a text box in order to show SI values with units (kg, mA, ...)

6.37 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/prepareDocu.cpp File Reference

Functions

- function [prepareDocu](#) ()
- function [prepareMFiles](#) ()
- function [postprocessCPP](#) ()
- function [recursdir](#) (in baseDir, in searchExpression)

6.37.1 Function Documentation

6.37.1.1 function `postprocessCPP` ()

6.37.1.2 function `prepareDocu` ()

6.37.1.3 function `prepareMFiles` ()

6.37.1.4 function `recursdir` (in *baseDir*, in *searchExpression*)

6.38 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/RangeTable.cpp File Reference

Classes

- class [RangeTable](#)
Extends a `UITable` with [PhyEditBox](#) capabilities.

6.39 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/show_stddev.cpp File Reference

6.40 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/Singleton.cpp File Reference

Classes

- class [Singleton](#)

6.41 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/sinum.cpp File Reference

Functions

- function [sinum](#) (in *str*, in *uni*)

6.41.1 Function Documentation

6.41.1.1 function `sinum` (in *str*, in *uni*)

6.42 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/sipre.cpp File Reference

Functions

- function [sipre](#) (in *val*, in *sgf*, in *pfx*, in *trz*)

6.42.1 Function Documentation

6.42.1.1 function sipre (in *val*, in *sgf*, in *px*, in *tr*)

6.43 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/SweepPanel.cpp File Reference

Classes

- class [SweepPanel](#)

Holds and organizes all tables options and buttons for the sweep configuration inside the fMPlot window.

6.44 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/SweepTable.cpp File Reference

Classes

- class [SweepTable](#)

Extends a UITable with [PhyEditBox](#) capabilities and provides the sweep sequence.

6.45 C:/Users/VSM/Documents/MATLAB/VSM-Prog/doxygen/cpp_vs/TempController.cpp File Reference

Classes

- class [TempController](#)